



KERN & Sohn GmbH

Ziegelei 1

D-72336 Balingen

E-Mail: info@kern-sohn.com

Tel: +49-[0]7433- 9933-0

Fax: +49-[0]7433-9933-149

Internet: www.kern-sohn.com

Operating instruction Analytical balance

KERN ALJ_N/ALS_N

Version 2.2

10/2007

GB



ALJ_N/ALS_N-BA-e-0722



KERN ALJ_N/ALS_N

Version 2.2 10/2007

Operating Manual

Electronic Analytical Balance

Table of Contents

1	Technical data	4
2	Declaration of conformity	7
3	Basic Information (General)	9
3.1	Proper use	9
3.2	Improper Use	9
3.3	Warranty	9
3.4	Monitoring of Test Resources	9
4	Basic Safety Precautions	10
4.1	Pay attention to the instructions in the Operation Manual	10
4.2	Personnel training	10
5	Transport and storage	10
5.1	Testing upon acceptance	10
5.2	Packaging	10
6	Unpacking, Setup and Commissioning	10
6.1	Installation Site, Location of Use	10
6.2	Unpacking	11
6.2.1	Placing	11
6.2.2	Scope of delivery	11
6.3	Mains connection	12
6.4	Connection of peripheral devices	12
6.5	Initial Commissioning	12
6.5.1	Stability display	12
6.5.2	Balance zero display	12
7	Adjustment	13
7.1.1	Automatic adjustment with internal weight (ALJ only)	13
7.1.2	Adjustment functions "P1 CAL"	15
7.2	Verification	22
8	Operating elements	24
8.1	Backlit display	24
8.2	Keyboard overview	24
9	User menu	25
9.1	Navigation in the menu	27

10	Weighing functions P6 func	28
10.1	Taring	28
10.2	Standard weighing unit	29
10.3	Weighing units switch-over	29
10.4	Underfloor weighing	30
10.5	Selecting the operating mode	31
10.5.1	Settings for menu item P6.1 Ffun for „ALL“	32
11	Main function “P3 rEAd”	33
11.1	Filter settings	33
11.2	Rest position display	34
11.3	Auto Zero	35
11.4	Changing readability – decimal place	36
12	Further useful functions “P7 othEr”	37
12.1	Display background illumination	37
12.2	Acoustic signal for key operation	38
12.3	Printing balance parameters	39
13	GLP/ISO log “P2 GLP”	40
14	Operating modes “P6 Func”	43
14.1	Parts counting	44
14.2	Weighing with tolerance range	47
14.3	Percent determination	49
14.3.1	Determining the reference weight by weighing (function F4 PrcA)	49
14.3.2	Determining the reference weight by entering numeric value (function F4 Prcb)	50
14.4	Density determination – functions “d_Co“ and “d_Li“	51
15	Data output RS 232C “P4 Print”	52
15.1	Technical data	52
15.2	Pin allocation of the balance output plug (front view)	52
15.3	Interface cable	53
15.4	Parameter for RS 232C interface	54
15.4.1	Navigation in the menu	54
15.4.2	Menu overview	54
15.5	Communication protocol / remote control commands	56
15.5.1	Manual output	56
15.5.2	PC controlled output	57
15.5.3	Output of date/time	57
16	Error messages	58
17	Service, maintenance, disposal	58
17.1	Cleaning	58
17.2	Service, maintenance	58
17.3	Disposal	58
18	Instant help	59

1 Technical data

KERN	ALJ 160-4NM	ALJ 220-4NM
<i>Weighing range (max)</i>	160 g	220 g
<i>Readability (d)</i>	0,1 mg	0,1 mg
<i>Minimum load</i>	10 mg	10 mg
<i>Verification value</i>	1 mg	1 mg
<i>Verification class</i>	I	I
<i>Reproducibility</i>	0,2 mg	0,2 mg
<i>Linearity</i>	± 0,2 mg	± 0,2 mg
<i>Stabilization time (typical)</i>	4 sec	4 sec
<i>Minimum unit weight at piece counting</i>	> 0,5 mg	> 0,5 mg
<i>Warm-up time</i>	8 hours	
<i>Adjustment weight</i>	internal	
<i>Reference quantities at piece counting</i>	10, 20, 50, freely selectable	
<i>Weighing Units</i>	mg, g, ct	
<i>Electric Supply</i>	230V/50 Hz (Euro)11V AC	
<i>Operating temperature</i>	+ 18° C + 30° C	
<i>Humidity of air</i>	max. 80 % (not condensing)	
<i>Underfloor weighing</i>	Clevis type eyelet, standard	
<i>Housing (B x D x H) mm</i>	206 x 335 x 335	
<i>Dimensions glass windscreen</i>	168 x 160 x 225 (weighing space)	
<i>Weighing plate (stainless steel) mm</i>	85 mm	
<i>Weight kg (net)</i>	6.3 kg	
<i>Interface</i>	RS 232C	

KERN	ALJ 220-5 DNM	ALJ 310-4N
Weighing range (max)	60 g/220 g	310 g
Readability (d)	0,01 mg/0,1 mg	0,1 mg
Minimum load	1 mg	-
Verification value	1 mg	-
Verification class	I	-
Reproducibility	0,02 mg/0,1 mg	0,1 mg
Linearity	± 0,1 mg/0,2 mg	± 0,3 mg
Minimum unit weight at piece countin	> 0,1 mg	> 0,5 mg
Stabilization time (typical)	13 sec/5 sec	4 sec
Warm-up time	8 hours	
Adjustment weight	internal	
Reference quantities at piece counting	10, 20, 50, freely selectable	
Weighing Units	mg, g, ct	
Electric Supply	230V/50 Hz (Euro)11V AC	
Operating temperature	+ 18° C + 30° C	
Humidity of air	max. 80 % (not condensing)	
Underfloor weighing	Clevis type eyelet, standard	
Housing (B x D x H) mm	206 x 335 x 335	
Dimensions glass windscreen	168 x 160 x 225 (weighing space)	
Weighing plate (stainless steel)	70 mm	85 mm
Weight kg (net)	6.3 kg	
Interface	RS 232C	

KERN	ALS 120-4N	ALS 220-4N
<i>Weighing range (max)</i>	120 g	220 g
<i>Readability (d)</i>	0.1 mg	
<i>Reproducibility</i>	0.2 mg	
<i>Linearity</i>	± 0.2 mg	
<i>Warm-up time</i>	8 hours	
<i>Recommended adjustment weight, not added (class)</i>	100 g (E2)	200 g (E2)
<i>Minimum unit weight at piece counting</i>	> 0.5 mg	
<i>Reference quantities at piece counting</i>	10, 20, 50, freely selectable	
<i>Weighing Units</i>	mg, g, ct, oz, ozt, GN, mom, dwt, teal	
<i>Stabilization time (typical)</i>	4 sec.	
<i>Electric Supply</i>	230V/50 Hz (Euro) 11V AC	
<i>Operating temperature</i>	+ 18° C + 30° C	
<i>Humidity of air</i>	max. 80 % (not condensing)	
<i>Underfloor weighing</i>	Clevis type eyelet, standard	
<i>Housing (B x D x H) mm</i>	206 x 335 x 335	
<i>Dimensions glass windscreen</i>	168 x 160 x 225 (weighing space)	
<i>Weighing plate (stainless steel) mm</i>	85 mm	
<i>Weight kg (net)</i>	6.3 kg	
<i>Interface</i>	RS 232C	

2 Declaration of conformity



KERN & Sohn GmbH
 D-72322 Balingen-Frommern
 Postfach 4052
 E-Mail: info@kern-sohn.de

Tel: 0049-[0]7433- 9933-0
 Fax: 0049-[0]7433-9933-149
 Internet: www.kern-sohn.de

Konformitätserklärung

EC-Konformitätserklärung
EC- Déclaration de conformité
EC-Dichiarazione di conformità
EC- Declaração de conformidade
EC-Deklaracja zgodności

EC-Declaration of -Conformity
EC-Declaración de Conformidad
EC-Conformiteitverklaring
EC- Prohlášení o shode
ЕС-Заявление о соответствии

D	Konformitäts- erklärung	Wir erklären hiermit, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Normen übereinstimmt.
GB	Declaration of conformity	We hereby declare that the product to which this declaration refers conforms with the following standards.
CZ	Prohlášení o shode	Tímto prohlašujeme, že výrobek, kterého se toto prohlášení týká, je v souladu s níže uvedenými normami.
E	Declaración de conformidad	Manifetamos en la presente que el producto al que se refiere esta declaración está de acuerdo con las normas siguientes
F	Déclaration de conformité	Nous déclarons avec cela responsabilité que le produit, auquel se rapporte la présente déclaration, est conforme aux normes citées ci-après.
I	Dichiarazione di conformità	Dichiariamo con ciò che il prodotto al quale la presente dichiarazione si riferisce è conforme alle norme di seguito citate.
NL	Conformiteit- verklaring	Wij verklaren hiermede dat het product, waarop deze verklaring betrekking heeft, met de hierna vermelde normen overeenstemt.
P	Declaração de conformidade	Declaramos por meio da presente que o produto no qual se refere esta declaração, corresponde às normas seguintes.
PL	Deklaracja zgodności	Niniejszym oświadczamy, że produkt, którego niniejsze oświadczenie dotyczy, jest zgodny z poniższymi normami.
RUS	Заявление о соответствии	Мы заявляем, что продукт, к которому относится данная декларация, соответствует перечисленным ниже нормам.

Electronic Balance: KERN ALS...N, ALJ...NM, ALJ...N

Mark applied	EU Directive	Standards
CE	89/336EEC EMC	EN 61000-4-2 :1999 EN 61000-4-3 :1996 EN 61000-4-4 : 1999 EN 61000-4-5 : 1998 EN 61000-4-6 : 1999 EN 61000-4-11 : 1997 EN 55022 :2000

Date: 27.02.2007

Signature: _____

**Gottl. KERN & Sohn GmbH
 Management**

Gottl. KERN & Sohn GmbH, Ziegelei 1, D-72336 Balingen, Tel. +49-[0]7433/9933-0, Fax +49-[0]7433/9933-149

Konformitätserklärung

EC-Konformitätserklärung
EC- Déclaration de conformité
EC-Dichiarazione di conformità
EC- Declaração de conformidade
EC-Deklaracja zgodności

EC-Declaration of -Conformity
EC-Declaración de Conformidad
EC-Conformiteitverklaring
EC- Prohlášení o shode
ЕС-Заявление о соответствии

D	Konformitäts- erklärung	Wir erklären hiermit, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Normen übereinstimmt. Diese Erklärung gilt nur in Verbindung mit der Konformitätsbescheinigung einer benannten Stelle.
GB	Declaration of conformity	We hereby declare that the product to which this declaration refers conforms with the following standards. This declaration is only valid with the certificate of conformity by a notified body.
CZ	Prohlášení o shode	Tímto prohlašujeme, že výrobek, kterého se toto prohlášení týká, je v souladu s níže uvedenými normami. Toto prohlášení platí pouze ve spojitosti s deklarací o souladu uvedeného pracoviště se směrnicemi EU.
E	Declaración de conformidad	Manifestamos en la presente que el producto al que se refiere esta declaración está de acuerdo con las normas siguientes. Esta declaración solo será válida acompañada del certificado de conformidad de una institución renombrada.
F	Déclaration de conformité	Nous déclarons avec cela responsabilité que le produit, auquel se rapporte la présente déclaration, est conforme aux normes citées ci-après. Cette déclaration est valide seulement avec un certificat de conformité d'un organisme notifié.
I	Dichiarazione di conformità	Dichiariamo con ciò che il prodotto al quale la presente dichiarazione si riferisce è conforme alle norme di seguito citate. Questa dichiarazione sarà valida solo se accompagnata dal certificato di conformità della parte nominale.
NL	Conformiteit- verklaring	Wij verklaren hiermede dat het product, waarop deze verklaring betrekking heeft, met de hierna vermelde normen overeenstemt. Deze verklaring geldt uitsluitend in verbinding met het certificaat van overeenstemming vanwege een daarmee belaste instantie.
P	Declaração de conformidade	Declaramos por meio da presente que o produto no qual se refere esta declaração, corresponde às normas seguintes. Esta declaração vale só em combinação com um certificado de conformidade duma instituição nomeada.
PL	Deklaracja zgodności	Niniejszym oświadczamy, że produkt, którego niniejsze oświadczenie dotyczy, jest zgodny z poniższymi normami. Niniejsze oświadczenie obowiązuje wyłącznie w połączeniu z oświadczeniem o zgodności danego miejsca.
RUS	Заявление о соответствии	Мы заявляем, что продукт, к которому относится данная декларация, соответствует перечисленным ниже нормам. Эта декларация действует совместно с удостоверением соответствия названной лаборатории.

Electronic Balance: KERN ALS...N, ALJ...NM

EU Directive	Standards	EC-type-approval certificate no.	Issued by	Modell
90/384/EEC	EN 45501	TCM 128/07 - 4519	CMI	ALJ..NM

Date: 27.02.2007

Signature:



Gottl. KERN & Sohn GmbH
Management

Gottl. KERN & Sohn GmbH, Ziegelei 1, D-72336 Balingen, Tel. +49-[0]7433/9933-0, Fax +49-[0]7433/9933-149

3 Basic Information (General)

3.1 Proper use

The balance you purchased is intended to determine the weighing value of material to be weighed. It is intended to be used as a “non-automatic” balance, i.e. the material to be weighed is manually and carefully placed in the centre of the weighing plate. As soon as a stable weighing value is reached the weighing value can be read.

3.2 Improper Use

Do not use balance for dynamic weighings. In the event that small quantities are removed or added to the material to be weighed, incorrect weighing results can be displayed due to the “stability compensation” in the balance. (Example: Slowly draining fluids from a container on the balance.)

Do not leave permanent load on the weighing plate. This may damage the measuring system.

--

Impacts and overloading exceeding the stated maximum load (max) of the balance, minus a possibly existing tare load, must be strictly avoided. Balance may be damaged by this.

Never operate balance in explosive environment. The serial version is not explosion protected.

The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.

The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.

3.3 Warranty

Warranty claims shall be voided in case

- Our conditions in the operation manual are ignored
- The appliance is used outside the described uses
- The appliance is modified or opened
- Mechanical damage or damage by media, liquids, natural wear and tear
- The appliance is improperly set up or incorrectly electrically connected
- The measuring system is overloaded

3.4 Monitoring of Test Resources

In the framework of quality assurance the measuring-related properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (www.kern-sohn.com) with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited DKD calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

4 Basic Safety Precautions

4.1 Pay attention to the instructions in the Operation Manual

Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.

4.2 Personnel training

The appliance may only be operated and maintained by trained personnel.

5 Transport and storage

5.1 Testing upon acceptance

When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

5.2 Packaging

Keep all parts of the original packaging in case you need to return the appliance.

Only use original packaging for returning.

Before sending, disconnect all connected cables and loose/movable parts.

Attach possibly existing transport safeguards. Secure all parts, e.g. weighing plate, power unit etc., to prevent slipping and damage.

6 Unpacking, Setup and Commissioning

6.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

Therefore, observe the following for the installation site:

- Place the balance on a firm, level surface;
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight;
- Protect the balance against direct draughts due to open windows and doors;
- Avoid jarring during weighing;
- Protect the balance against high humidity, vapours and dust;
- Do not expose the device to extreme dampness for longer periods of time. Non-permitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for ca. 2 hours at room temperature.
- Avoid static charge of goods to be weighed or weighing container.


Major display deviations (incorrect weighing results) may be experienced should electromagnetic fields (e.g. due to mobile phones or radio equipment), static electricity accumulations or instable power supply occur. Change location or remove source of interference.

6.2 Unpacking

Carefully remove the balance from the packaging, remove plastic cover and setup balance at the intended workstation.

6.2.1 Placing

Design layout of weighing balance:



Level balance with foot screws until the air bubble of the water balance is in the prescribed circle.

6.2.2 Scope of delivery

Serial accessories:

- Balance (see chapter 6.2.1)
- Mains power supply
- Operating Manual

6.3 Mains connection

Power is supplied via the external mains adapter. The stated voltage value must be the same as the local voltage.

Only use original KERN mains adapters. Using other makes requires consent by KERN.

6.4 Connection of peripheral devices

Before connecting or disconnecting of additional devices (printer, PC) to the data interface, always disconnect the balance from the power supply.

With your balance, only use accessories and peripheral devices by KERN, as they are ideally tuned to your balance.

6.5 Initial Commissioning



In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap. 1).

During this warming up time the balance must be connected to the power supply (mains, accumulator or battery).


The accuracy of the balance depends on the local acceleration of gravity.

Strictly observe hints in chapter Adjustment.

6.5.1 Stability display

The appearance of the stability symbol [] on the display indicates that the balance is in a stable state. If the status is instable the [] display disappears.

6.5.2 Balance zero display

If an exact zero reading is not displayed on the balance in spite of the weighing dish being empty, press the TARE key and the balance will start resetting to zero []

.

7 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each balance must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the balance has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the balance periodically in weighing operation.

7.1.1 Automatic adjustment with internal weight (ALJ only)

With the internal adjustment weight, the weighing accuracy can be checked and re-adjusted at any time.


The automatic adjustment function is always enabled. You can start adjustment at any time by pressing the **F**-key 3 times.

Automatic adjustment is started each time


- after the weighing balance was disconnected from the mains
- when temperature changes take place
- at the end of a time interval

Temperature / time controlled adjustment is taking place:

5 minutes before the start of automatic adjustment, this will be announced by a "°C", (change of temperature) or a "▶" (after a certain time interval ends) symbol on the display.



The user must complete his/her weighing process within this time. After 5 minutes [**CAL 30**] appears on the display.




This starts a "count down" of 30 seconds [**CAL 30**] → [**CAL 0**]. During these 30 seconds it is possible to cancel the adjustment by pressing the **TARE** key. This makes the balance return to weighing mode in order to e. g. complete an unfinished measurement.

After a further 5 minutes the automatic adjustment will restart and [**CAL 30**] appears.

Selection options of internal adjustment:


1. Temperature controlled adjustment

A change in temperature of 0.8 °C triggers the weighing balance to send a message announcing internal adjustment:



2. Time-controlled adjustment

After a set time interval (1 – 12 h) in the menu (**P1.5 CALt**) has passed the weighing balance will send a message about internal adjustment:



3. By menu invocation "P1.1 iCAL", see chpt. 6.5.4

4. Press the F-key three times

7.1.2 Adjustment functions “P1 CAL“

Menu “P1 CAL“ [adjustment]

P1.1	iCAL		[internal adjustment] ALJ only
P1.2	ECAL		[external adjustment] *
P1.3	tCAL		[adjustment test]
P1.4	ACAL	both/nonE/tenno/tinnE	[automatic adjustment] ALJ only *
P1.5	CALt	1 h ÷ 12h	[time setting for automatic adjustment] ALJ only *
P1.6	CALr	YES/no	[printout adjustment log]

* = locked for verifiable appliances [Er 9 lock].

1. Function iCAL: Manual start of internal adjustment (ALJ only)

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again. Ensure that there are no objects on the weighing plate.	P1.1 iCAL
⇒ Confirm by pressing the F -key; adjustment takes place automatically. After successful adjustment the balance automatically returns to weighing mode. Information: Adjustment errors are indicated by an acoustic signal and an error message will appear on the display. Press the TARE -key and repeat the adjustment process.	CAL ↓ 0.0000 g
⇒ You can cancel adjusting by pressing the TARE -key.	


2. Function ECAL: External adjustment

Carry out adjustment with the help of the recommended adjustment weight (see chapter 1 “Technical Specifications”).



Procedure when adjusting:

Observe stable environmental conditions. A warming up time (see chapter 1) is required for stabilization.


Menu activation:



Models ALJ



Models ALS




Modelle ALS:

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again	P1.1 ECAL
⇒ Press F -key. Ensure that there are no objects on the weighing plate.	unload
⇒ Press the PRINT -key and wait until the weight value of the required adjustment weight pops up.	CAL ↓ load ↓ 200,0000 g
⇒ Place the required calibration weight carefully in the centre of the weighing plate and press the PRINT key. CAL appears on the display, adjustment is started. The display shows UNLOAD , adjustment is completed.	CAL ↓ unload
⇒ Remove adjusting weight, the balance will automatically return to weighing mode. Information: Adjustment errors are indicated by an acoustic signal and an error message will appear on the display. Press the TARE -key and repeat the adjustment process.	- - - - - ↓ 0.0000 g

3. Function tCAL: Adjustment test

Here, deviation from the last adjustment is determined. This is only a check, i.e. no values are changed.

Menu activation:



Models ALJ:

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again.	P1.1 iCAL
⇒ Press the UNIT -key repeatedly until "P1.3 tCAL" appears.	P1.3 tCAL
⇒ Press the F -key; automatic adjustment test is carried out. The result is displayed	CAL ↓ diff ↓ d 0,0042 g
⇒ Press down and hold the PRINT -key until the balance returns to the menu.	P1.3 tCAL
⇒ Press the TARE -key repeatedly until " SAVE ?" appears. Store executed changes by pressing the PRINT -key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	SAVE? ↓ 0.0000 g

Models ALS:

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again.	P1.1 ECAL
⇒ Press the UNIT -key repeatedly until "P1.2 tCAL" appears.	P1.2 tCAL
⇒ Press F key	unload
⇒ Press the PRINT -key and wait until the weight value of the required adjustment weight pops up.	CAL ↓ load ↓ 200,000 g
⇒ Place the required calibration weight carefully in the centre of the weighing plate and press the PRINT key. Adjustment test is carried out. The result is displayed	CAL ↓ diff ↓ d 0,0042 g
⇒ Press down and hold the PRINT -key until the balance returns to the menu.	P1.2 tCAL
⇒ Press the TARE -key repeatedly until " SAVE "? appears. Store executed changes by pressing the PRINT -key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	SAVE? ↓ 0.0000 g

4. Function ACAL: Automatic internal adjustment (only ALJ)

This menu item is used to set criteria for the start of automatic adjustment .

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again.	P1.1 ECAL
⇒ Press the UNIT -key repeatedly until the “ACAL“ appears.	P1.4 ACAL
⇒ Press the F -key; current setting is flashing	nonE
⇒ Use the UNIT key to select one of the settings below: nonE Automatic adjustment disabled tEmp Start automatic Adjustment after change in temperature timE Start automatic Adjustment dependent on setting of function P1.5 CALt both Start automatic Adjustment dependent on temperature and time	timE
⇒ Confirm setting by pressing the PRINT key. Weighing balance returns to menu.	P1.4 ACAL
⇒ Press the TARE -key repeatedly until “SAVE“? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

5. Function CAL t:
Time setting for start of automatic internal adjustment (ALJ only)


Set the time interval after which automatic adjustment is to be started in this menu item as follows:

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again.	P1.1 ECAL
⇒ Press the UNIT -key repeatedly until the “ P1.5 CALt “ function appears	P1.5 CALt
⇒ Press the F -key; current setting is flashing	1 h
⇒ Use the UNIT -key to select a time interval of 1 – 12 h	12 h
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P1.5 CALt
⇒ Press the TARE -key repeatedly until “ SAVE “? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

6. Function CALr: “print adjustment log“

Go to this menu item to enable the function used to produce printouts for adjustment data.

Menu activation:



Enable/disable function:

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.0000 g
⇒ Press F key	P1 CAL
⇒ Press the F -key again.	P1.1 ECAL
⇒ Press the UNIT -key repeatedly until the “CALr“ function appears	P1.6 CALr Models ALJ
	P1.3 CALr Models ALS
⇒ Press the F -key; current setting is flashing	no
⇒ Use the UNIT key to select one of the settings below: no Data output disabled yes Data output enabled	YES
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P1.6 CALr Models ALJ
	P1.3 CALr Models ALS
⇒ Press the TARE -key repeatedly until “SAVE“? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

Data output adjustment log:

The content of the data output is determined in menu **P2 GLP**. All parameters set to "YES" will be issued.

Example:

P2 GLP

P2.1	uSr		
P2.2	PrJ		
P2.3	Ptin		YES
P2.4	PdAt		YES
P2.5	PuS		YES
P2.6	PPrJ		YES
P2.7	PId		YES
P2.8	PFrn		YES

```
*** Automatic calibration report ***
Date   : 09/02/2007
Time   : 11:21:39
User Id : 12345678
Project Id: 87654321
Balance Id: 114493

Calibr. : Automatic
Difference: - 0.0002 g

Name .....
```

7.2 Verification

General introduction:

According to EU directive 90/384/EEC balances must be officially verified if they are used as follows (legally controlled area):

- For commercial transactions if the price of goods is determined by weighing.
- For the production of medicines in pharmacies as well as for analyses in the medical and pharmaceutical laboratory.
- For official purposes.
- For manufacturing final packages.

In cases of doubt, please contact your local trade in standard.

Verification instructions

An EU type approval exists for balances described in their technical data as verifiable. If a balance is used where obligation to verify exists as described above, it must officially verified and re-verified in regular intervals.


Re-verification of a balance is carried out according to the respective national regulations. The validity for verification of balances in Germany is e.g. 2 years.

The legal regulation of the country where the balance is used must be observed!

After verification the balance is sealed at the indicated positions.

Verification of the balance is invalid without the "seal".

Position of the “official seals“:



Balances with obligation to verify must be taken out of operation if:

- **The weighing result of the balance is outside the error limit.** Therefore, in regular intervals load balance with known test weight (ca. 1/3 of the max. load) and compare with displayed value.
- **The reverification deadline has been exceeded.**







8 Operating elements

8.1 Backlit display

Very contrastful display which can also be read in the darkness.



8.2 Keyboard overview

Key	Function
	<ul style="list-style-type: none">▪ Turn on/off
	<ul style="list-style-type: none">▪ Function key
	<ul style="list-style-type: none">• Switch-over key into operating modes
	<ul style="list-style-type: none">• Switching-over to other weighing units
	<ul style="list-style-type: none">• Printout of the weighed value on an external appliance (printer or PC)• Confirm/save settings
	<ul style="list-style-type: none">• Taring• Leave menu• Set weight display at zero

9 User menu

The user menu has seven main menus (P1 – P7), arranged in the following sub-menus:

P1 CAL [adjustment]

P1.1	iCAL		[internal adjustment] ALJ only
P1.2	ECAL		[external adjustment] *
P1.3	tCAL		[adjustment test]
P1.4	ACAL	both/nonE/tenno/tinnE	[automatic adjustment] ALJ only *
P1.5	CALt	1 h ÷ 12h	[time setting for automatic adjustment] ALJ only *
P1.6	CALr	YES/no	[printout adjustment log]

* = locked for verifiable appliances

P2 GLP [good lab practice]

P2.1	USr	–	[user]
P2.2	PrJ	–	[project]
P2.3	Ptin	YES/no	[printout time]
P2.4	PdAt	YES/no	[printout date]
P2.5	PUSr	YES/no	[printout user]
P2.6	PPrJ	YES/no	[printout project]
P2.7	PId	YES/no	[printout serial number weighing balance]
P2.8	PFr	YES/no	[frame printout]

P3 rEAd [main functions]

P3.1	AuE	Stand/Slouu/FASt	[filter settings]
P3.2	ConF	FASt_rEL/Fast/rEL	[rest control display]
P3.3	Auto	On/OFF	[Auto zero]
P3.4	Ldi9	ALuuAYS/neper/uu_StAb	[delete last decimal place]

P4 Print [parameter for serial interface RS 232]

P4.1	bAud	2400/4800/9600/19200	[Baud rate]
P4.2	CntA	YES/no	[continuous output standard weighing unit]
P4.3	Cntb	YES/no	[continuous output currently set weighing unit]
P4.4	rEPL	YES/no	[manual (press key) or automatic output]
P4.5	PStb	YES/no	[output stable/instable weighing values]
P4.6	Lo	000.0000	[input minimal weight for automatic output]







P5 unit		[weighing unit]
P5.1	StUn g/mg/ct/oz/ozt/dwt/t/mom/G	[standard weighing unit, e.g. "g"]
P5.2	mg	YES/no [mg - milligram]
P5.3	Ct	YES/no [ct – carat]
P5.4	oZ	YES/no [oz – ounce]
P5.5	oZt	YES/no [ozt – Troy ounce]
P5.6	dwt	YES/no [dwt – pennyweight]
P5.7	t	YES/no [t – Tael]
P5.8	nno	YES/no [mom - Momme]
P5.9	Gr	YES/no [gr – Grain]

P6 Func		[operating modes]
P6.1	FFun ALL/PcS/HiLo/PrcA/Prcb/d_Co/d_Li	[selection enabled operating modes]
P6.2	PcS	YES/no [piece counting]
P6.3	HiLo	YES/no [weighing with tolerance control]
P6.4	PrcA	YES/no [percentage weighing, reference weight determination "weighing"]
P6.5	Prcb	YES/no [percentage weighing, reference weight determination "numeric"]
P6.6	d_Co	YES/no [determination density "solids"]
P6.7	d_Li	YES/no [determination density "liquids"]

P7 othEr		[additional useful functions]
P7.1	bL	On/Aut/OFF Display background illumination
P7.2	bBEEP	On/OFF [touch tone]
P7.3	PrnS	[printout "weighing parameters"]


9.1 Navigation in the menu


Keyboard overview in menu:


Key	Function
	<ul style="list-style-type: none">• Call main menu• Call first parameter of a function• Shift and select number to be changed to the right
	<ul style="list-style-type: none">• Parameter selection (scroll down)• Decreasing the value for selected numeral or number
	<ul style="list-style-type: none">• Parameter selection (scroll up)• Increasing the value for selected numeral or number
	<ul style="list-style-type: none">• Return to main menu• Shift and select number to be changed to the left
	<ul style="list-style-type: none">• Confirm/save settings
	<ul style="list-style-type: none">• Exit menu without saving

Storing / jumping back to weighing mode

Any changes made in the balance memory will only be saved when the storing process is complete.


To achieve this, press the  key several times until “**SAVE**“? appears.

Any changes carried out are stored by pressing the  key.

To cancel changes, press the  key.

Afterwards the balance automatically jumps back to weighing mode.

10 Weighing functions P6 func

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key. The balance will carry out a self-test	
⇒ Your balance is ready to weigh as soon as the " 0.00 " display appears.	0.00 g
⇒ Put on items to be weighed, weighed value is displayed. Wait until the stability display appears  . Read the weighing value.	19.68 g
⇒ To turn off the balance press the ON/OFF key	

10.1 Taring

The dead weight of any weighing container may be tared away by pressing a button, so that the following weighings show the net weight of the goods to be weighed.

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	
⇒ Your balance is ready to weigh as soon as the " 0.00 " display appears.	0.00 g
⇒ Put on items to be weighed, weighed value is displayed.	19.68 g
⇒ To start the taring process press the TARE key. The weight of the container is now saved internally.	0.00 g
⇒ Place the goods to be weighed into the tare container. Read the weight of the goods on the display.	53.25 g

The taring process can be repeated any number of times, e.g. when adding several components for a mixture (adding).

The limit is reached when the whole weighing range is exhausted.

After removing the taring container the total weight is displayed as negative display.

10.2 Standard weighing unit

Selected weighing unit will be retained even after disconnection from the mains.

Operation:	Display:
⇒ Turn on balance by pressing the ON/OFF key	0.00 g
⇒ Press F key	P1 CAL
⇒ Press the UNIT -key several times until " P5 Unit " appears	P5 unit
⇒ Press F key	P5.1 StUn
⇒ Press the F -key again; the currently set weighing unit is flashing	„g“ (example)
⇒ Press the UNIT -key repeatedly until the desired weighing unit appears. [g] → [mg] → [ct] → [oz] → [ozt] → [dwt] → [t] → [mom] → [G] → [g] Note: In calibrated model selection is restricted to options [g] → [mg] and → [ct] .	mg (example)
⇒ Confirm setting by pressing the PRINT key. Weighing balance returns to menu.	P5.1 StUn
⇒ Press the TARE -key repeatedly until " SAVE "? appears	„SAVE“
Store changes carried out by pressing the PRINT key. The balance returns to weighing mode; the display show the set weighing unit. The set weighing unit remains even after disconnection from the mains.	0.0 mg

10.3 Weighing units switch-over

A weighing unit may be changed by pressing the **UNIT**-key several times:

[g] → [mg] → [ct] → [oz] → [ozt] → [dwt] → [t] → [mom] → [G] → [g]

Information:

In calibrated models selection is restricted to options **[g] →** and **[mg] → [ct]**.

The different weighing models have integrated different foreign weighing units. For details please refer to chpt. 1, "Technical specifications".

10.4 Underfloor weighing

Objects unsuitable for placing on the weighing scale due to size or shape may be weighed with the help of the flush-mounted platform.

Proceed as follows:

- Switch off balance.
- Open the closing lid on the bottom of your balance.
- Suspend hook for underfloor weighing carefully and completely.
- Place your balance over an opening.
- Suspend the goods to be weighed from the hook and carry out the weighing.




Fig. 1: Setup of balance for underfloor weighing



CAUTION

- **Always ensure that all suspended objects are stable enough to hold the desired goods to be weighed safely (danger of breaking).**
- **Never suspend loads that exceed the stated maximum load (max) (danger of breaking)**


Always ensure that there are no persons, animals or objects that might be damaged underneath the load.




NOTICE


After completing the underfloor weighing the opening on the bottom of the balance must always be closed (dust protection).


10.5 Selecting the operating mode

You can enable and disable functions in menu “**P6 Func**“ which are then available to the user without having to go back to the menu each time. All activated operating modes can be called directly by pressing the  key.


Operation:	Display:	Description:
⇒ Turn on balance by pressing the ON/OFF key	0.00 g	
⇒ Press F key	P1 CAL	
⇒ Press the UNIT -key repeatedly until “P6 Func“ appears.	P6 Func	
⇒ Press the F -key	P6.1 Ffun	
⇒ Press the F -key again	ALL	All operating modes enabled in menu P6.2 to P6.7 are callable by pressing the  key
⇒ Press the UNIT -Taste.	PcS	Parts counting
⇒ Press the UNIT -Taste.	HiLo	Tolerance weighing
⇒ Press the UNIT -Taste.	PrcA	Display of percentage
⇒ Press the UNIT -Taste.	Prcb	Percentage display via manual default
⇒ Press the UNIT -Taste.	d_Co	Density determination “solids“
⇒ Press the UNIT -Taste.	d_LI	Density determination “liquids“

Confirm settings by pressing the **PRINT** key. Balance jumps back to sub-menu **6.1.Ffun**.

Is the parameter set to **ALL** in menu **6.1.Ffun**, access to the functions set to **YES** in the **P4FFunc** menu is available via the  key.
(see chpt. 10.5.1)


However, if in menu **6.1.Ffun** one of the above parameters is selected, e. g. **PcS**, the  key can merely be used to call just this particular menu item.

10.5.1 Settings for menu item P6.1 Ffun for „ALL“

This is where you select the menu items that will be from then on selectable via the  key. **no** – Function not available.
YES – Function available.

Operation:	Display:
⇒ Selecting the menu item	P6 Func
⇒ Press the F -key	P6.1 Ffun
⇒ Press the UNIT key until the operating mode “piece counting” appears	P6.2 PcS
⇒ Press the F -key; the current setting: “no” = deactivated, “yes” = activated	no
⇒ To change the setting, press the UNIT key.	YES
⇒ Setting will be imported by pressing the PRINT key.	P6.2 PcS
⇒ Press the UNIT key until the operating mode “tolerance weighing” appears. Activation as described for piece counting	P6.3 HiLo

Please repeat this sequence of operations for any other operating mode available.



For this, press the **TARE** key repeatedly until **Save ?** appears on the display. Store changes carried out by pressing the **PRINT** key.

11 General functions “P3 rEAd”

11.1 Filter settings

This menu item allows the balance to be set according to specific ambient conditions and measuring purposes.

Operation:	Display:
⇒ Menu call „P3 rEAd“: Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until “P3 rEAd“ appears.	P3 rEAd
⇒ Confirm by pressing the F -key; function “P3.1 AuE“ appears.	P3.1 AuE
⇒ Press the F -key again; the current setting is flashing.	FASt
⇒ Use the UNIT key to select one of the filter settings below: FASt = Sensitive and fast (very quiet site) StAnd = Standard Slouu = Insensitive but slow (site exposed to vibration)	StAnd
⇒ Confirm setting by pressing the PRINT key. Weighing balance returns to menu.	P3.1 AuE
⇒ Press the TARE -key repeatedly until “SAVE“? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

11.2 Rest position display

Operation:	Display:
⇒ Call menu item "P3 rEAd": Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until "P3 rEAd" appears.	P3 rEAd
⇒ Press F key	P3.1 AuE
⇒ Press the UNIT -key until "P3. 2 ConF" appears	P3. 2 ConF
⇒ Press the F -key; current setting is flashing	FASt_rEL
⇒ Use the UNIT key to select one of the filter settings below: FASt_rEL = rest position control fast FASt = rest position control fast and accurate rEL = rest position control accurate	FASt
⇒ Confirm setting by pressing the PRINT key. Weighing balance returns to menu.	P3. 2 ConF
⇒ Press the TARE -key repeatedly until "SAVE" ? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

11.3 Auto Zero

This function is used to tare small variations in weight automatically.

In the event that small quantities are removed or added to the material to be weighed, incorrect weighing results can be displayed due to the “stability compensation” in the balance. (e.g. slow flow of liquids from a container placed on the balance, evaporating processes).

When apportioning involves small variations of weight, it is advisable to switch off this function.

Operation:	Display:
⇒ Call menu item "P3 rEAd": Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until "P3 rEAd" appears.	P3 rEAd
⇒ Press F key	P3.1 AuE
⇒ Press the UNIT -key until „P3.3 Auto“ appears.	P3.3 Auto
⇒ Press the F -key; current setting is flashing	On
⇒ Use the UNIT key to select one of the filter settings below: on = Auto zero function turned on off = auto zero function turned off	OFF
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P3.3 Auto
⇒ Press the TARE -key repeatedly until "SAVE" appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

11.4 Changing readability – decimal place

Readability may be reduced by 1 digit on the weighing balances, as required. The last decimal place will be rounded and removed from the display.

Operation:	Display:
⇒ Call menu item "P3 rEAd": Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until "P3 rEAd" appears.	P3 rEAd
⇒ Press F key	P3.1 AuE
⇒ Press the UNIT -key until „P3.4 Ldi9“ appears.	P3.4 Ldi9
⇒ Press the F -key; current setting is flashing	AlwAYS
⇒ Use the UNIT key to select one of the filter settings below: AlwAYS =shows last decimal place nEuEr =does not show last decimal place uu_StAb =last decimal place only displayed at stable weighing values	nEuEr
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P3.4 Ldi9
⇒ Press the TARE -key repeatedly until "SAVE" appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

12 Further useful functions “P7 othEr“

Here, you can set the parameters that influence the operation of the balance, such as background lighting and key sounds.

12.1 Display background illumination

Operation:	Display:
⇒ Calling menu item “ P7 othEr ”: Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until “ P7 othEr “ appears.	P7 othEr
⇒ Press F key	P7.1 bl
⇒ Press the F -key; current setting is flashing	On
⇒ Use the UNIT key to select one of the filter settings below: ON = Background illumination on OFF = Background illumination off Aut = Backlighting automatically switched off 10 seconds after achieving stable weighing value	OFF
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P7.1 bl
⇒ Press the TARE -key repeatedly until “ SAVE “? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

12.2 Acoustic signal for key operation

Operation:	Display:
⇒ Calling menu item " P7 othEr ": Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until " P7 othEr " appears.	P7 othEr
⇒ Press F key	P7.1 bl
⇒ Press the UNIT -key	P7.2 bEEP
⇒ Press the F -key; current setting is flashing	On
⇒ Use the UNIT key to select one of the filter settings below: ON = Acoustic signal turned on OFF = Acoustic signal turned off	OFF
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P7.2 bEEP
⇒ Press the TARE -key repeatedly until " SAVE "? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

12.3 Printing balance parameters

Operation:	Display:
⇒ Calling menu item " P7 othEr ": Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until " P7 othEr " appears.	P7 othEr
⇒ Press F key	P7.1 bl
⇒ Press the UNIT -key until " P7.3 PrnS " appears.	P7.3 PrnS
⇒ Press the F -key; output of balance parameter is achieved via RS 232 interface.	P7.3 PrnS
⇒ Press the TARE -key repeatedly until "SAVE" ? appears.	SAVE ?
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g

13 GLP/ISO log “P2 GLP“

Quality assurance systems require printouts of weighing results as well as of correct adjustment of the balance stating date and time and balance identification. The easiest way is to have a printer connected.

The content of the data output is determined in menu “P2 GLP“. All parameters set to “YES“ will be issued.

Examples:

P2 GLP

P2.1	uSr		max. 8 digits
P2.2	PrJ		max. 8 digits
P2.3	Ptin		YES
P2.4	PdAt		YES
P2.5	PuS		YES
P2.6	PPrJ		YES
P2.7	Pid		YES
P2.8	PFrn		YES

Date	: 09/02/2007
Time	: 11:21:39
User Id	: 12345678
Project Id	: 87654321
Balance Id	: 114493
100.0216 g	

P2.8 PFrn: YES	


Date	:20.03.07
Time	:11.31.07
UserID	:Mustermann
Balance ID	:180151
19.3406 g	

P2.8 PFrn: no	
Date	:20.03.07
Time	:11.31.07
UserID	:Mustermann
Balance ID	:180151
19.3406 g	

Defination of a standard log:

Operation:	Display:
⇒ Call menu point „P2 GLP“ : Press F key	0.0000 g ↓ P1 CAL
⇒ Press the UNIT -key repeatedly until „P2 GLP“ appears.	P2 GLP
⇒ F -key. The menu item “P2.1 Usr“ for entering the user name appears.	P2.1 Usr
⇒ Press the F -key; the first digit is flashing.	—
⇒ Enter user name or number (max. 8 digits) with the help of the arrow keys	Example: Data input: FERN
	Data output: KERN

Control via arrow keys:



Overview data input / data output:

.	0	1	2	3	4	5	6	7	8	9		
-	0	1	2	3	4	5	6	7	8	9		
A	b	c	d	E	F	G	H	I	J	K	L	M
A	B	C	D	E	F	G	H	I	J	K	L	M
n	o	P	q	r	S	t	U	V	W	X	Y	Z
N	O	P	Q	R	S	T	U	V	W	X	Y	Z
ä	b	c	d	e	F	g	h	i	J	k	l	M
a	b	c	d	e	f	g	h	i	j	k	l	M
n	o	P	q	r	S	t	u	v	w	x	y	Z
n	o	P	q	r	s	t	u	v	w	x	y	Z

⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P2.1 Usr
⇒ Press the UNIT -key; the next menu item " P2.2 PrJ " used to enter project name appears.	P2.2 PrJ


⇒ Press the F -key; the first digit is flashing.	—																																
⇒ Press the arrow keys to enter project name or number (max. 8 digits).																																	
⇒ Confirm settings by pressing the PRINT key. Weighing balance returns to menu.	P2.2 PrJ																																
⇒ Press the UNIT -key; the next menu item “ P2.3 Ptin “ used to display the time appears.	P2.3 Ptin																																
⇒ Press the F -key; current setting is flashing	no																																
⇒ Use the UNIT key to select one of the settings below: no = time not displayed YES = time displayed	YES																																
⇒ Confirm setting by pressing the PRINT key. Weighing balance returns to menu.	P2.3 Ptin																																
⇒ Press the UNIT -key; the next menu item “ P2.4 PdAt “ used to display the date appears.	P2.4 PdAt																																
<p>Use the same steps to enter further GLP-Parameter.</p> <p>Menu overview</p> <table border="0"> <tr> <td>P2.1</td> <td>Usr </td> <td>—</td> <td>[user]</td> </tr> <tr> <td>P2.2</td> <td>PrJ </td> <td>—</td> <td>[project]</td> </tr> <tr> <td>P2.3</td> <td>Ptin </td> <td>YES/no</td> <td>[printout time]</td> </tr> <tr> <td>P2.4</td> <td>PdAt </td> <td>YES/no</td> <td>[printout date]</td> </tr> <tr> <td>P2.5</td> <td>PUSr </td> <td>YES/no</td> <td>[printout user]</td> </tr> <tr> <td>P2.6</td> <td>PPrJ </td> <td>YES/no</td> <td>[printout project]</td> </tr> <tr> <td>P2.7</td> <td>Pid </td> <td>YES/no</td> <td>[printout serial number weighing balance]</td> </tr> <tr> <td>P2.8</td> <td>PFr </td> <td>YES/no</td> <td>[frame printout]</td> </tr> </table>		P2.1	Usr	—	[user]	P2.2	PrJ	—	[project]	P2.3	Ptin	YES/no	[printout time]	P2.4	PdAt	YES/no	[printout date]	P2.5	PUSr	YES/no	[printout user]	P2.6	PPrJ	YES/no	[printout project]	P2.7	Pid	YES/no	[printout serial number weighing balance]	P2.8	PFr	YES/no	[frame printout]
P2.1	Usr	—	[user]																														
P2.2	PrJ	—	[project]																														
P2.3	Ptin	YES/no	[printout time]																														
P2.4	PdAt	YES/no	[printout date]																														
P2.5	PUSr	YES/no	[printout user]																														
P2.6	PPrJ	YES/no	[printout project]																														
P2.7	Pid	YES/no	[printout serial number weighing balance]																														
P2.8	PFr	YES/no	[frame printout]																														
⇒ Press the TARE -key repeatedly until “ SAVE “? appears.	SAVE ?																																
⇒ Any changes carried out are stored by pressing the PRINT key. To cancel changes, press the TARE key. Afterwards the balance automatically jumps back to weighing mode.	0.0000 g																																

14 Operating modes “P6 Func“

The following operating modes may be enabled/disabled in this menu item:

- Parts counting
- Weighing with tolerance range
- Percent determination
- Density determination of solids/liquids

Call of operating modes (see also chpt. 10.5):




14.1 Parts counting

Before the balance can count parts, it must know the average part weight (i.e. reference). Proceed by putting on a certain number of the parts to be counted. The balance determines the total weight and divides it by the number of parts (the so-called reference quantity). Counting is then carried out on the basis of the calculated average piece weight.

As a rule:


The higher the reference quantity the higher the counting exactness.

- **Call add-up function**



Press the **UNIT**-key to select the desired reference piece number.


- **Reference piece number 10, 20 or 50**



Confirm selected reference piece number by pressing the **PRINT**-key (e.g. 20)


- **Reference**

Place as many pieces to add-up as required by the set reference piece number.




Remove reference weight. The balance is now in parts counting mode counting all units on the weighing plate.

- for selection "optional reference piece number" FrEE




– Select the place to be changed by pressing the **F**-key



– Select the digit by pressing the **UNIT**-key

- Confirm the entered reference piece number by pressing the **PRINT**-key
- "LoAd" appears on the display.



- Place as many counting parts on the balance as the set reference quantity requires, confirm by pressing the **PRINT** key.


Information:

If there is no load on the weighing plate when the **PRINT**-key is pressed, "Er5 out" will appear briefly on the display before the display of the balance returns automatically to weighing mode.



- Remove reference weight. The balance is now in parts counting mode counting all units on the weighing plate.

- **Return to weighing mode**




14.2 Weighing with tolerance range

For weighing with tolerance ranges you can enter individual upper and lower limits. For tolerance controls such as dosaging, apportioning or sorting the scale will display violated upper or lower limits and show the tolerance tag.


The triangular tolerance marker (▲) in the upper part of the display shows whether the goods to be weighed are within the two tolerance limits. The tolerance marker is only in operation during operating mode tolerance weighing; it is otherwise not visible.

The tolerance marker provides the following information:




- lo** Goods to be weighed below tolerance limit
- ok** Goods to be weighed within tolerance range
- Hi** Goods to be weighed above tolerance limit

- **Call function**




- **Setting the lower tolerance limit "Lo"**



- Select the digit to be changed by pressing the **F**-key; each digit currently active will be flashing.
- Select the digit by pressing the **UNIT**-key
- Confirm the entered lower tolerance tag by pressing the **PRINT**-key.

- **Setting the upper tolerance limit "Hi"**




- Select the digit to be changed by pressing the **F**-key; each digit currently active will be flashing.
- Select the digit by pressing the **UNIT**-key
- Confirm the entered lower tolerance tag by pressing the **PRINT**-key.

The balance is now in checkweighing mode Put on goods to be weighed, tolerance control is started

Information:

If invalid values are entered such as lower tolerance limit greater than upper tolerance limit, the balance will issue the error message “Er8 outr” and return automatically to weighing mode.

- **Return to weighing mode**




14.3 Percent determination

Percent determination allows weight display in percent, in relation to a reference weight.


14.3.1 Determining the reference weight by weighing (function F4 PrcA)

- **Call function**



- **Make reference**

- Place the reference weight
- Confirm by pressing the **PRINT** key.




- This weight is adopted as reference (100%).

Remove reference weight. The balance is now in percent determining mode. Place the load on the balance; percentage value in relation to reference body is shown on display:

14.3.2 Determining the reference weight by entering numeric value (function F4 Prcb)

- **Call function**



- **Make reference**

- You will be asked to enter the reference weight by a flashing message.




- Select the digit to be change by pressing the **F**-key; each time the currently active digit will be blinking.



- Select the digit by pressing the **UNIT**-key
- Confirm the entered reference weight by pressing the **PRINT**-key.

Now you can place the test objects onto the weighing plate; the percentage to the reference body is displayed

- **Return to weighing mode**



14.4 Density determination – functions “d_Co“ and “d_Li“

For procedure of density determination please see the operating instructions "Set density"

15 Data output RS 232C “P4 Print“


15.1 Technical data

- 8-bit ASCII Code
- 8 data bits, 1 stop bit, no parity bit
- Baud rate optional from 2400 - 19200 Baud
- For operation with interface faultless operation is only ensured with the correct KERN – interface cable (max. 2m)

Various transmission modes are available:




- Manually after pressing the **PRINT** key
- Continuously, according to setting
- Automatically according to stability display
- Prompted by external device
(For remote control commands see chpt. 15.5))

15.2 Pin allocation of the balance output plug (front view)



Pin 2: Receive data
Pin 3: Transmit data
Pin 5: Signal ground

15.3 Interface cable



Balance
 2 (RxD)
 3 (TxD)
 4 (DTR)
 5 (GND)
 6 (DSR)

15.4 Parameter for RS 232C interface



15.4.1 Navigation in the menu

- ⇒ Turn on balance by pressing the **ON/OFF** key
- ⇒ Press the **F**-key; first menu item "**P1 CAL**" appears.
- ⇒ Press the **UNIT**-key repeatedly until "**P4 Print**" appears.
- ⇒ Confirm by pressing the **F**-key; the first sub-menu for setting the Baud rate "**P4.1 b Aud**" appears.
- ⇒ Press and hold the **UNIT**-key until the selection you wish to use appears.
P4.1 b Aud → **P4.2 b CntA** → **P4.3 b Cntb** → **P4.4 rEPL** → **P4.5 PStb** → **P4.6 Lo** → **P4.1 b Aud**
- ⇒ Press the **F**-key again; the current setting is flashing.
- ⇒ Press the **UNIT**-key repeatedly until the desired setting appears (see chpt. 15.4.2).
- ⇒ Confirm settings by pressing the **PRINT** key. Weighing balance returns to menu. If desired, carry out settings for further menu items as described above.
- ⇒ Press the **TARE**-key repeatedly until "**SAVE**"? appears.

Store executed changes by pressing the **PRINT** -key. To cancel changes, press the **TARE** key. Afterwards the balance automatically jumps back to weighing mode.

15.4.2 Menu overview

Function	Choice	Description of options
F-key	UNIT-key	
P4.1 b Aud	2400	2400 bps
	4800	4800 bps
	9600	9600 bps
	19200	19200 bps
P4.2 b CntA	no	Continuous output disabled in standard weighing unit, Output requires key operation
	yes	Continuous output in standard weighing unit

P4.3 b Cntb	no	Continuous output disabled in current weighing unit, Output requires key operation
	yes	Continuous output in current weighing unit
P4.4 rEPL Models in non- verifiable setting only	no	Manual issue after pressing the PRINT key.
	yes	Automatic issue of first stable weighing value Sequence of operations: 1. Taring 2. Place weight, issue of first stable weighing value 3. Renewed output only possible after weight was removed. Condition: Display +/- 50 display steps from zero 4. Place next weight.
P4.5 PStb	no	Output even for unstable weighing value
	yes	Output for stable weighing value only
P4.6 Lo	000,005	<p>Input minimum weight for automatic output:</p>  <ul style="list-style-type: none"> Select the place to be changed by pressing the F-key  <ul style="list-style-type: none"> Select the digit by pressing the UNIT-key Confirm by pressing the PRINT key <p>Weighing value issued automatically if current weighing value exceeds entered minimum value. The next weighing value will not be issued unless the weighing value has meanwhile dropped below the entered weighing value.</p>

15.5 Communication protocol / remote control commands

Instruction:	Meaning of Instruction:
Z	Set weight display at zero
T	Taring
S	Send weight value immediately
SI	Send stable weight value
SU	Send stable weighing value in current weighing value
SUI	Send weighing value immediately in current weighing unit
C1	Turn on continuous transmission in standard weighing unit
C0	Turn off continuous transmission in standard weighing unit
CU1	Turn on continuous transmission in current weighing unit
CO1	Turn off continuous transmission in current weighing unit
PC	Send all implemented instructions

Complete each instruction with **CR LF**.

15.5.1 Manual output

The user can start output manually by pressing the **PRINT**-key (for settings see chapter 15.4.2).

Data record format:

1	2	3	4 - 12	13	14 - 16	17	18
Stability indicator	Blank	Signs	Weight	Blank	Unit	CR	LF

Stability indicator Space if stable,
 ? if not stable
 ^ if overload
 v if underload

Presign: Space, if positive
 negative sign, if negative

Weight: 9 signs, right justified

Unit: 3 signs, left justified

15.5.2 PC controlled output

Response message of balance after remote instruction was sent:

XX_ Instruction:
 XX_A CR LF Instruction accepted; will be executed
 XX_I CR LF Instruction received; impossible to carry out
 XX_^ CR LF Instruction received but *time overflow error* occurred
 XX_v CR LF Instruction received, but insufficient load
 XX_E CR LF Error during execution, timeout for stable weighing value exceeded

Data record format:

1 - 3	4	5	6	7	8 - 16	17	18 - 20	21	22
Remote control command	Blank	Stability indicator	Blank	mark	Weight	Blank	Unit	CR	LF

Instruction: 1. up to 3 signs

Stability indicator: Space if stable,
 ? if not stable
 ^ if overload
 v if underload

Presign: Space, if positive
 negative sign, if negative

Weight: 9 signs, right justified

Unit: 3 signs, left justified

15.5.3 Output of date/time

Output of date and time is enabled in menu item "P2 GLP":

- PdAt – yes
- Ptin - yes

16 Error messages

Er1 Hi	Initial weight error
Er2 nuLL	Value below allowed range
Er3 FuL1	Value above allowed range
Er4 FuL2	Weighing range exceeded
Er5 rout	Value outside allowed range e.g. tare value ≤ 0 , Reference weight = 0
Er7 tout	Disconnecting time too short
Er8 outr	Input outside range e.g. for tolerance control: Input upper limit < lower limit
Er9 Lock	Function blocked
Er10 cal	Adjustment error e.g. incorrect adjustment weight}

17 Service, maintenance, disposal

17.1 Cleaning

Before cleaning, please disconnect the appliance from the operating voltage.

Please do not use aggressive cleaning agents (solvents or similar agents), but a cloth dampened with mild soap suds. Ensure that no liquid penetrates into the device and wipe with a dry soft cloth.

Loose residue sample/powder can be removed carefully with a brush or manual vacuum cleaner.

Spilled weighing goods must be removed immediately.

17.2 Service, maintenance

The appliance may only be opened by trained service technicians who are authorized by KERN.

Before opening, disconnect from power supply.

17.3 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

18 Instant help

In case of an error in the program process, briefly turn off the balance and disconnect from power supply. The weighing process must then be restarted from the beginning.

Help:

Fault

Possible cause

The displayed weight does not glow.

- *The balance is not switched on.*
- *The mains supply connection has been interrupted (mains cable not plugged in/faulty).*
- *Power supply interrupted.*
- *(Rechargeable) batteries are inserted incorrectly or empty*
- *No (rechargeable) batteries inserted.*

The displayed weight is permanently changing

- *Draught/air movement*
- *Table/floor vibrations*
- *Weighing plate has contact with other objects.*
- *Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)*

The weighing result is obviously incorrect

- *The display of the balance is not at zero*
- *Adjustment is no longer correct.*
- *Great fluctuations in temperature.*
- *Warm-up time was ignored.*
- *Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)*

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.