- 3 Betriebsanleitung  
Elektronische Zählwaage
- 36 Operating Instructions  
Counting Balances
- 70 Mode d'emploi  
Balances de comptage
- 103 Instrucciones de funcionamiento  
Balanzas cuentapiezas
- 136 Istruzioni d'uso  
Bilance contapezzi

# Inhalt

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# Contents

36	Contents
36	Intended Use
37	Warnings and Safety Precautions
38	General View of the Balance
40	Getting Started
44	Operating the Balance
44	Weighing
46	Calibration/Adjustment
47	Counting
51	Toggle between Units
53	Weigh Averaging
55	Settings
55	Setting the Parameters (Menu)
56	Parameter Settings (Overview)
60	Generating a Printout
62	Data Interface
63	Error Codes
64	Care and Maintenance
65	Instructions for Recycling
66	Overview
66	Specifications
68	Dimensions (Scale Drawings)
69	CE Marking

## Supplement:

Key label strip for the control panel showing the procedure for "counting"

# Intended Use

The models of the CGB Series are precision balances designed for the measurement of mass covering a range from 0.01 g and 64 kg.

CGB models meet the highest requirements on the accuracy and reliability of weighing results through the following features:

- Efficient filtering-out of unfavorable ambient conditions, such as vibration, drafts, etc.
- Stable and reproducible weighing results
- Excellent readability under any lighting conditions
- Rugged, durable weighing system

CGB balances save work and speed up simple routine applications through these features:

- Reference sample quantity adjustable using the selector
- Mass unit conversion by toggling between units (second weight unit)
- Extremely fast response times of approx. 1 second
- Total ease of operation
- Runs on 6 size C alkaline manganese batteries rated to 1.5 V, max. 8,100 mAh
- Serial RS-232 port for connection to a PC
- Two lines in the printout are configurable to show your company name, for example

# Warnings and Safety Precautions

The balance has been constructed in accordance with the European Directives as well as international regulations and standards for operation of electrical equipment, electromagnetic compatibility, and stipulated safety requirements. Improper use or handling, however, can result in property damage or personal injury.

Read these operating instructions thoroughly before using your balance to prevent damage to the equipment. Keep these instructions in a safe place for future reference.

Follow the instructions below to ensure safe and trouble-free operation of your balance:

- ⚠ Do not use this balance in a hazardous area/location
- ⚠ Only for use with LISTED direct plug in power supply.
- ⚠ CAUTION. Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.
- ⚠ Make sure that the voltage rating printed on the AC adapter is identical to your local line voltage.

⚠ Warning When Using Prewired RS-232 Connecting Cables:  
RS-232 cables purchased from other manufacturers often have incorrect pin assignments for use with KERN balances. Be sure to check the pin assignment against the respective chart before connecting the cable, and disconnect any lines marked "Internally Connected" (e.g., pin 6). Failure to do so may damage or even completely ruin your balance and/or peripheral device.


- The only way to turn the power off completely is to disconnect the AC adapter
- Connect only KERN accessories and options, as these are optimally designed for use with your CGB balance
- Protect the AC adapter and the balance from moisture


When cleaning your balance, make sure that no liquid enters the balance housing; use only a slightly moistened cloth to clean the balance.

Do not open the balance housing. In case you have any problems with your balance:

- Contact your local KERN office, dealer or service center

# General View of the Balance

Item		Item	
1	Weighing pan	12	 key: Sets the display to zero. Zeroing is only allowed within the $\pm 2\%$ range of the balance's maximum capacity.
2	Battery compartment (below the control unit)	13	[ON/OFF] key: Switches the display on or off. (The balance remains energized in the standby mode – depending on the menu setting.)
3	Control unit	14	Display: average piece weight
4	Foot	15	Weight display in the basic unit selected
5	[UNIT] key: Toggles between 2 weight units or to counting mode	16	Shows that a value is in the tare memory
6	[START] key: starts measurement of the average piece weight	17	Symbol for standby operation
7	Selector: selects reference sample quantity	18	Symbol for zeroing
8	[OK] key: saves reference sample quantity	19	"Busy" symbol
9	[DELETE] key: cancels operation This key is generally used for the following: – End counting – Cancel a calibration/adjustment routine in progress	20	Weight unit or piece count (pcs)
10	[PRINT] key: generates a printout This key sends displayed values via the built-in data interface to an on-line Data Printer or to a computer.	21	Manufacturer's label
11	[TARE] key: tares (subtracts weight) The weight of an empty container is subtracted so that the net weight of a sample in the container is always shown afterwards.	22	DC jack
		23	Data interface port
		24	Main display

# Getting Started

## Guarantee

The guarantee shall expire in the event that

- the directions in the operating instructions are not followed properly
- the device is used for applications beyond those described in these operating instructions
- the device is altered or opened
- the device is mechanically damaged or is damaged due to exposure to media and/or liquids
- normal wear of the device has occurred
- the device has not been adequately set up or its electrical installation is incorrect
- the measuring equipment is overloaded

## Storage and Shipping Conditions

Do not expose the balance to shock, vibration, moisture or extreme temperature.

## Unpacking the Balance

- After unpacking the balance, check it immediately for any visible damage
- If you see any sign of damage, proceed as directed in the chapter entitled "Care and Maintenance" under the section on "Safety Inspection"

Save the box and all parts of the packaging until you have successfully installed your balance. Only the original packaging provides the best protection for shipment. Before packing your balance, unplug all connected cables to prevent damage.

## Equipment Supplied

The equipment supplied includes the components listed below:

- Balance with weighing pan mounted in place
- Plug-in AC adapter
- Operating instructions


## Installation Instructions

CGB balances are designed to provide reliable weighing results under normal ambient conditions. When choosing a location to set up your balance, observe the following so that you will be able to work with added speed and accuracy:

- Set up the balance on a stable, even surface
- Avoid placing the balance next to a heater or exposing the balance to direct sunlight
- Protect the balance from direct drafts that come from open windows or doors
- Avoid exposing the balance to strong vibration during weighing
- Protect the balance from aggressive chemical vapors
- Avoid exposing the balance to extreme moisture

## Conditioning the Balance:

Moisture in the air can condense on the surfaces of a cold balance whenever it is brought into a substantially warmer place. Therefore, condition a cold balance for approx. 2 hours at room temperature, leaving it unplugged from AC power. Afterwards, keep the balance continuously connected to AC power.


### AC Power Connection/Safety Precautions

Use only

- Original AC adapters
- Insert the right-angle plug into the jack on the balance
- Then plug the Class 2 AC adapter into any electrical outlet. No additional safety precautions are required.


### Information on Radio Frequency Interference Warning!

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause interference to radio communications. It has been tested and found to comply with the limits for a Class A computing device pursuant to Subpart J of Part 15 of FCC Rules, which are designed to provide reasonable protection against such interference, when operated in a commercial environment. Operation of this equipment in a residential area is likely to cause interference, in which case the user, at his own expense, will be required to take whatever measures may be required to correct the interference.


### Installing the Batteries

- Batteries are not included with the equipment supplied.
- △ Use only commercially available size C alkaline manganese batteries rated to 1.5 V max. 8,100 mAh.
- Lay the balance on its side.
- To open the battery compartment, lift the compartment cover.
- Install six 1.5-volt size C batteries in the compartment.
- Make sure to connect the positive and negative poles correctly.
- △ Used batteries are classified as waste that requires special handling (not "household" waste). Dispose of rechargeable batteries according to your country's applicable special waste disposal regulations.
- To close the battery compartment: Press down on the cover until it clicks into place.

For automatic shutoff of a battery-operated balance:

- Select "automatic shutoff" in the menu: Please see the chapter on "Settings" to select code **B. 4. 1**
- > The balance will then shut off automatically when it has not been operated for more than approx. 2 minutes.


### Connecting Electronic Peripheral Devices

- Make absolutely sure to unplug the balance from AC power before you connect or disconnect a peripheral device (printer or PC) to or from the interface port.

For battery operation:

- > Press the [ON/OFF] key to shut off the power completely (balance does not go into standby mode).

### Warmup Time

To deliver exact results, the balance must warm up for at least 30 minutes after initial connection to AC power. Only after this time will the balance have reached the required operating temperature.

# Operating the Balance

## Weighing

### Features

- Zeroing the balance  
Depending on the prevailing ambient conditions, the display may not show a zero readout even though there is no load on the balance. If the weight shown is less than 2% of the balance's maximum weighing capacity, you can zero the balance.
- Taring the balance (saving a [container] weight in the memory by subtracting it from the display)  
Tare the balance with an empty container on the weighing pan to obtain a readout of the net weight after filling the container.
- Printing weights

### Preparation






- Turn on the balance:  
Press [ON/OFF]
- > A self-test is performed
- To change settings: please see the chapter on "Settings"
- To load the factory settings: please see the chapter on "Settings," parameter 9 - 1

### Example

Determine a weight

Settings: factory settings

Step	Press key	Display/Printout									
1. Turn on balance	[ON/OFF]										
2. If necessary, zero balance											
3. Place empty container on balance (in this example, 11.5 g)											
4. Tare balance	[TARE]										
5. Place sample in container on balance (in this example, 120.5 g)											
6. Print weight	[PRINT]	<table><tbody><tr><td>N</td><td>+</td><td>120.5 g</td></tr><tr><td>T</td><td>+</td><td>11.5 g</td></tr><tr><td>G#</td><td>+</td><td>132.0 g</td></tr></tbody></table>	N	+	120.5 g	T	+	11.5 g	G#	+	132.0 g
N	+	120.5 g									
T	+	11.5 g									
G#	+	132.0 g									

# Calibration/Adjustment

## Purpose

Technically, calibration means determining the difference between the weight readout and the true weight (mass) of a sample. Calibration does not entail making any changes within the balance. Adjustment is the correction of this difference between the measured value displayed and the true weight (mass) of a sample, or the reduction of the difference to an allowable level within the maximum permissible error limits. Because the CGB balance automatically performs both calibration and adjustment, we use the term "calibrate/ adjust" to mean both in this manual.

## Features

Calibration/adjustment can only be performed when

- there is no load on the balance,
  - the balance is set to zero, and
  - the internal signal is stable.
- If these conditions are not met, an error code is displayed.

The value of the weight on the balance must not differ from the nominal weight by more than 2%.

You can use any of the following weight units to calibrate/adjust: g, kg, lb (1.4)

You can block calibration/adjustment of the balance (1.5)

## Example

Calibrate/adjust the balance  
Settings: factory settings

Step	Press key	Display/Printout
1. Turn on balance, if power is off	[ON/OFF]	0.0 g
2. If necessary, zero the balance	→0←	0.0 g
3. Start calibration/adjustment Calibration weight is displayed without wt. unit	[TARE] (> 2 sec)	+ 500.0
4. Place the indicated calibration weight on balance (in this case, 500 g)		500.0
After calibration/adjustment, calibration weight is displayed with wt. unit		+ 500.0 g
5. Remove calibration weight		0.0 g

# Counting

## Purpose

With the counting program, you can determine the number of parts that have approximately the same average piece weight.

## Features

- The minimum capacity is one digit related to the resolution of the particular weight unit selected.
- The resolution can be selected for the average piece weight to be saved as a reference and calculation of the piece count.
- The criterion for saving the average sample weight ("Wref") and calibration of the piece count can be selected. This criterion applies to the following key functions: [→0←], [TARE], [START] and [OK].
- Reinitialize without quitting the counting application.
- The average piece weight and reference sample quantity (piece count = "nRef") are automatically output via the data interface port after initialization, if you have selected the setting for printing several lines of data with nRef/wRef in the menu.
- Press [UNIT] to toggle between the piece count and weight.

△ Piece weights that show high variations or are too low will negatively affect the counting accuracy.

## Function Keys

[START] key:

- Press [START] to begin determining the average piece weight.
- The current weight is shown as a reference for the average piece weight that is saved when [OK] is pressed.
- Place the reference sample parts on the balance when "PIECE WEIGHT": **SAMPLE** is displayed.
- Depending on the menu code previously set, the balance is either automatically tared or not tared when the "Counting" program is started.
- The "Auto Zero" function is de-activated until the average piece weight has been saved.


[OK] key:

- Initializes the counting application using the reference sample quantity selected in the menu.
- Shows the average piece weight as "PIECE WEIGHT" after you have pressed the key once.
- If you press the [OK] key again, this will reinitialize the counting application.

If you press the [OK] key >2 sec.:

- The number of the last reference sample quantity "nRef" will be displayed (Init. or Opt.).

## Selector

- Turn this selector to set the reference sample quantity.
- Each time you change the selector setting, the new reference sample quantity will be briefly displayed.

### Reference Sample Updating

Automatic reference sample updating optimizes the counting accuracy. You can activate or de-activate this function in the menu.

The abbreviation opt for "optimizing" briefly appears in the "PIECE WEIGHT" display to indicate that reference sample updating has been completed if

- the criterion for the stability parameter selected in the menu has been met
- the current piece count is less than twice the original piece count
- the current piece count is less than 1,000
- the internally calculated piece count (such as 17.24 pcs) differs by less than  $\pm 0.3$  pcs from the whole number (17 pcs in this example)

## Factory Settings of the Parameters

Program selection:

Counting without weight unit toggling  
(2. 1. 4)

Resolution:

Standard: according to the displayed accuracy  
(3. 4. 1)

Storage parameter:

At stability (3. 5. 1)

Reference sample updating:

Automatic (3. 6. 2)

Autotare when [START] is pressed:

On (3. 7. 2)

### Additional Functions










In addition to the basic functions (power off, zeroing, taring and printing), you can also access the following functions from this application:

- Cancel initialization: [DELETE] key
- Show weight: [UNIT] key
- Calibrate/adjust balance: press [TARE] > 2 sec.

### Example 1:

Count parts into a container

Menu settings: factory settings

Step	Press key	Display/Printout
1. Turn on balance, if power is off.	[ON/OFF]	
2. Place empty container on balance and start counting.	[START]	 
3. Select reference sample quantity; for example, turn selector to "20".		Brief display of: 
4. Place reference sample quantity on balance (20 parts in this case).		
5. Save reference sample quantity. The balance now calculates the average piece weight.	[OK]	 
6. If necessary, update the reference sample using the optimizing function; in this case, add 2 to 20 additional parts to the weighing pan.		 
7. Now place unknown number of parts on pan (in this example, 260).		 
8. Print piece count. If you have selected "nRef" and "wRef" in the menu, these lines will be printed:	[PRINT]	<pre> nRef +      32 wRef +  1.1438 g Qnt  +      260 pcs G    +      297.4 g                     </pre>
9. Unload balance		
10. Repeat counting starting from step 7, if desired.		

### Example 2:

Counting parts into a full container

Settings: (changes in factory settings required for this example)

Autotare when the [START] is key is pressed: Off (Code 3. 7. 1)

Step	Press key	Display/Printout
1. Turn on balance, if power is off.	[ON/OFF]	
2. Place empty container on balance and tare.	[TARE]	0.0 g NET
3. Place full container on balance and initialize counting.	[START]	+ 267.3 g NET + 267.3 SAMPLE
4. Set reference sample quantity; in this example, turn selector to "20".		Brief display of: rEF 20
5. Remove selected quantity of reference parts from container.		+ 249.1 g NET + 249.1 SAMPLE
6. Save reference piece count. The balance now calculates the number of parts and displays the piece count.	[OK]	+ 274 pcs NET + 249.1 + 0.9100
7. Print piece count. If you have selected "nRef" and "wRef" in the menu, these lines will be printed:	[PRINT]	nRef + 20 wRef + 0.9100 g Qnt + 274 pcs N + 249.1 g T + 50.0 g G# + 299.1 g
8. Unload balance		0 pcs NET
9. Repeat counting starting from step 6, if desired.		

# Toggle between Units

## Purpose

With this application program, you can toggle the display of a weight back and forth between two units.

## Features

- Mass unit conversion of a displayed weight by toggling
- Other features as for the basic weighing function

## Factory Settings of the Parameters

Toggle between weights:

Counting without toggling weight unit  
(2. 1. 4)

Weight unit 1: grams (1. 7. 2)

Weight unit 2: pounds (3. 1. 5)

## Additional Functions

In addition to the functions for:

- turning off balance power: [ON/OFF] key
- zeroing the balance: [→0←] key
- taring the balance: [TARE] key
- printing: [PRINT] key

you can also access the following functions from this application:

- toggle between weight units 1 and 2: [UNIT] key
- calibrate/adjust balance: press [TARE] for > 2 sec.

Menu Code	Unit	Conversion Factor	Printout
(1. 7. 1) (3. 1. 1)	Grams (o)	1.0000000000	o
(1. 7. 2) (3. 1. 2)	Grams (g)	1.0000000000	g
(1. 7. 3) (3. 1. 3)	Kilograms	0.0010000000	kg
(1. 7. 4) (3. 1. 4)	Carats	5.0000000000	ct
(1. 7. 5) (3. 1. 5)	Pounds	0.00220462260	lb
(1. 7. 6) (3. 1. 6)	Ounces	0.03527396200	oz
(1. 7. 7) (3. 1. 7)	Troy ounce	0.03215074700	ozt
(1. 7. 8) (3. 1. 8)	Hong Kong taels	0.02671725000	tlh
(1. 7. 9) (3. 1. 9)	Singapore taels	0.02645544638	tlS
(1. 7. 10) (3. 1. 10)	Taiwanese taels	0.02666666000	tlt
(1. 7. 11) (3. 1. 11)	Grains	15.43235835000	GN
(1. 7. 12) (3. 1. 12)	Pennyweights	0.64301493100	dwt
(1. 7. 14) (3. 1. 14)	Parts per pound	1.12876677120	/lb
(1. 7. 15) (3. 1. 15)	Chinese taels	0.02645547175	tlc
(1. 7. 16) (3. 1. 16)	Mommes	0.26670000000	mom
(1. 7. 17) (3. 1. 17)	Austrian carats	5.00000000000	K
(1. 7. 18) (3. 1. 18)	Tola	0.08573333810	tol
(1. 7. 19) (3. 1. 19)	Baht	0.06578947437	bat
(1. 7. 20) (3. 1. 20)	Mesghal	0.21700000000	MS

### Example

Toggle unit from grams [g] (1st unit) to pounds [lb] (2nd unit)

Settings (changes in the factory settings required for this example):

Menu: Counting with toggling weight unit (2. 1. 15)

Step	Press key	Display/Printout
1. Turn on balance, if power is off	[ON/OFF]	0.0 g
2. Load balance (in this example, with 314.3 g)		+ 314.3 g
3. Toggle to pounds [lb] unit	[UNIT]	+ 0.6930 lb
4. Print weight	[PRINT]	G + 0.6930 lb
5. Toggle to grams [g] unit	[UNIT]	+ 314.3 g

# Weigh Averaging

## Purpose

Use this application program to determine weights of moving samples (such as animals) or weights under unstable ambient conditions. In this application, the balance calculates the average from a defined number of individual weight measurements called “subweighing operations.”

## Features

- The measured result displayed is the arithmetic mean shown in the preselected weight unit; a triangle indicates that this is a calculated value.
- You can adjust the selector to the number of subweighing operations before starting weigh averaging using the [OK] key.
- The number of subweighing operations remaining to be performed is displayed during a measurement.
- To display the number of subweighing operations for averaging, press the [OK] key for more than 2 sec.
- Press the [UNIT] key to toggle between the calculated result and the weight readout.
- The results are automatically output via the interface port if you have selected the printout option with data ID codes in the menu.

## Function Keys


[OK] key:

- Saves the number of subweighing operations you selected for averaging and starting weigh averaging.
- If you press the [OK] key again, averaging will be reinitialized.

[DELETE] key:

- Cancels initialization.

## Selector

- Turn this selector to set the number of subweighing operations.
- Each time you change the selector setting, the new number of subweighing operations will be briefly displayed.

## Preparation


- Select the “weigh averaging” program in the menu: see chapter on “Settings” to select code **2. 1. 12**

### Example

Determine the weight of a sample under extremely unstable ambient conditions by calculating the average of 20 subweighing operations.

Settings (changes in the factory settings required for this example):

Menu: "Weigh averaging" application program, code 2.1.12

Step	Press key	Display/Printout
1. Turn on balance, if power is off.	[ON/OFF]	
2. Tare balance	[TARE]	0.0 g
3. Select number of subweighing operations (20 measurements in this example)		rEF 20 (briefly)
4. Place sample on balance (weight readout fluctuates; here, for example, by about 275 g)		+ 8888
5. Start measurement	[OK]	+ 8888 20 19 ... 1
After 20 subweighing operations		+ 275.5 g 
If the print format is set to include data ID codes, the following is printed:		RES + 275.5 g
6. Unload balance		+ 275.5 g  (stable display)
7. Delete result	[DELETE]	
8. Repeat procedure starting from step 4, if desired.		

# Settings

## Setting the Parameters (Menu)

This means you configure, or adapt, the balance to your user requirements by selecting the parameter settings from a menu.

Example: Adapt balance to place of installation with "Extreme vibration"

Select (Code **i. i. 4**)

Step	Press key	Display
1.	Turn off balance	[ON/OFF]
2.	Turn balance back on and, while all segments are displayed, press	[ON/OFF] [TARE] briefly
○	To navigate within a menu level; the last menu option is followed by the the first option	<b>i.</b> 2. ... 9. <b>i.</b>
3.	Select the 2nd menu level	[PRINT] <b>i. i.</b>
4.	Select the 3rd menu level	[PRINT] <b>i. i. 2 o</b>
5.	In menu level 3: Select the desired option repeatedly	Press [TARE] repeatedly <b>i. i. 4</b>
6.	Confirm new setting; "o" indicates the currently set option	Press [PRINT] for 2 sec. <b>i. i. 4 o</b>
○	Return to next higher menu level (from the 3rd menu level)	[PRINT] <b>i.</b>
○	Set other menu codes, if desired	[PRINT], [TARE]
7.	Save parameter settings and exit menu or	Press [TARE] for 2 sec.
○	Exit menu without saving changes	[ON/OFF]
>	Restart the application	<b>0.0 g</b>



## Parameter Settings (Overview)


o Factory setting

√ User setting

Menu	1st Menu Level	2nd Menu Level	3rd Menu Level	Factory Setting	Menu Item
Menu	1. Weighing	1.1. Adapt filter	1.1.1	o	Minimum vibration
			1.1.2		Normal vibration
			1.1.3		Strong vibration
			1.1.4		Extreme vibration
		1.3. Stability range	1.3.1	o	1/4 digit
			1.3.2		1/2 digit
			1.3.3		1 digit
			1.3.4		2 digits
			1.3.5		4 digits
		1.4. Unit for calibration weight	1.4.1	o	Grams
			1.4.2		Kilograms
			1.4.3		Pounds
		1.5. Calibration/adjustment	1.5.1	o	Accessible
			1.5.2		Blocked
		1.6. Auto zero	1.6.1	o	On
			1.6.2		Off
		1.7. Weight unit 1 (basic unit)	1.7.1	o	Grams (o)
			1.7.2		Grams (g)
	1.7.3		Kilograms		
	1.7.4		Carats		
1.7.5	Pounds				
1.7.6	Ounces				
1.7.7	Troy ounces				
1.7.8	Hong Kong taels				
1.7.9	Singapore taels				
1.7.10	Taiwanese taels				
1.7.11	Grains				
1.7.12	Pennyweights				
1.7.14	Parts per pound				
1.7.15	Chinese taels				
1.7.16	Mommes				
1.7.17	Austrian carats				
1.7.18	Tola				
1.7.19	Baht				
1.7.20	Mesghal				
	1.8. See next page				
	2. through 9.	See following pages			

# Generating a Printout

## Purpose

You can generate a printout of weights as well as other measured values and identification codes for documentation purposes. You can format the printout to meet individual requirements.

## Features

You can choose to print an individual weight or net, gross and tare weights.

Line format: You can configure a data ID code of up to 6 characters each for the values to be printed; this data ID code is printed at the beginning of a line.

You can generate printouts either automatically or by pressing the [PRINT] key; the printout can be dependent or independent of the balance's stability parameter.

You can have the following values output automatically if menu code **7. 1. 4** (printout with data ID codes) is selected:

- Average piece weight (wRef)
- Reference sample quantity (nRef)

## Factory Settings of the Parameters

Manual/automatic print mode:

In the manual mode, an individual value is printed when you press the print key; in the automatic mode, values are printed depending on the stability parameter selected for the balance:

Manual after stability (**5. 1. 2**).

Line format:

You can configure a data ID code of up to 6 characters for each of the values to be printed: printout of the net, tare and gross values; reference sample quantity, average piece weight with data ID codes (**7. 1. 4**).

- To set the parameters, please see the chapter on "Settings"

Header Lines:

Examples

Printout without Data ID Codes:

The value currently displayed is printed (weight or calculated value with unit)

```
+ 1530.0 g
+ 58.562 ozt
+ 253 pcs
```

Weight in grams  
Weight in Troy ounces  
Piece count

Printout with Data ID Codes:

The value currently displayed can be printed with a data ID code of up to 6 characters at the beginning of each line. You can use this data ID code to identify a weight as a net value (N) or as a calculated piece count (Qnt).

```
N + 153.0 g
T + 023.4 g
G + 155.3 g
G# + 163.0 g
Qnt + 253 pcs
Res + 275.5 pcs
```

Current net weight  
Value saved in tare memory  
Current gross weight  
Calculated gross weight  
Calculated piece count  
Averaged value: Result

Print Application Parameters (for Counting):

You can print the values configured for initialization of an application

```
wRef + 0.1400 g
nRef + 10
```

Average piece weight  
Reference sample quantity

Auto Print:

You can have the weight readout printed automatically. The display update interval depends on the operating status of the balance and on the balance model.

```
N + 153.0 g
Stat
Stat L
Stat H
```

Net weight  
Display blank  
Display underload  
Display overload

# Data Interface

## Purpose

Your CGB balance comes equipped with an interface port for connection to a computer or other peripheral device.

You can use an on-line computer to change, start and/or monitor the functions of the balance and the application program (such as counting).

## Features

Type of interface: serial interface port

Operating mode: full duplex

Standard: RS-232

Transmission rates:

150; 300; 600; 1,200; 2,400; 4,800; 9,600 baud

Parity: mark, space, odd, even

Character format:

1 start bit, 7-bit ASCII, parity, 1 or 2 stop bits

Handshake:

2-wire interface: via software (XON/XOFF)

4-wire interface: via hardware (CTS/DTR)

Communications mode: ASCII

Data output format: 16 or 22 characters

For details on the data interface

(such as the data output or input formats, pin assignment, etc.) contact KERN.

# Error Codes

Error codes are shown on the main display for approx. 2 seconds, after which the program automatically returns to the weighing mode.

Display	Cause	Solution
No segments appear on the display	No AC power is available	Check the AC power supply
	The AC adapter is not plugged in	Plug the AC adapter into an electrical outlet
	Automatic shutoff is configured	Turn on the balance
	Non-rechargeable/rechargeable batteries are drained	Replace or recharge batteries using external charger
H	The load exceeds the balance capacity	Unload the weighing pan
L	Something is touching the weighing pan	Move the object that is touching the weighing pan away
E 01	Data output is not compatible with the input format	Change the setting to the compatible one in the menu
E 02	Calibration/adjustment condition not met, e.g.: – balance not zeroed – balance is loaded	Calibrate/adjust only when zero is displayed Press <b>[→0←]</b> to zero balance display Unload balance
E 08	The balance was zeroed outside the zero range	Only zero the balance when in the allowable zero range of ± 2% of the maximum capacity
E 09	Taring not possible when the gross value ≤ zero	Press <b>[→0←]</b> to zero the balance
E 22	The weight is too light or there is no sample on the weighing pan	Increase the reference sample quantity
E 30	The interface port for printer output is blocked	Set code to reset menu or contact your local KERN Service Center
The weight readout changes constantly	Unstable ambient conditions (too much vibration or the balance) is exposed to a draft	Set up the balance in a different area or change the menu setting to adapt the balance
	A foreign object is caught between the weighing pan and the balance frame	Remove the foreign object
The weight readout is obviously wrong	The balance was not calibrated/adjusted; the balance was not zeroed before weighing	Calibrate/adjust the balance; zero before weighing

If any other errors occur, please contact your local KERN service center.

# Care and Maintenance

## Service

Regular servicing by a trained Kern dealer will extend the service life of your scale and ensure its continued weighing accuracy. Kern can offer you service contracts, with your choice of regular maintenance intervals.

The frequency of maintenance intervals depends on the operating conditions and your tolerance requirements.

Our trained dealers are happy to advise you on the best maintenance schedule for your needs.

## Repairs

Repair work must be performed by trained service technicians. Any attempt by untrained persons to perform repairs may lead to considerable hazards for the user.

## Cleaning

- ⚠ Disconnect the balance from the AC adapter and unplug any interface cables that are connected to the balance.
- ⚠ Make sure that no liquids enter the balance housing.
- ⚠ Do not use any aggressive cleaning agents (solvents or similar agents).
- Unplug the AC adapter from the wall outlet (mains supply).
- If you have a data cable connected to the interface, unplug it from the balance.
- Clean the balance using a piece of cloth that has been wet with a mild detergent (soap).
- After cleaning, wipe down the balance with a soft, dry cloth.

## Cleaning Stainless Steel Surfaces

Clean all stainless steel parts regularly. Remove the stainless steel weighing pan and thoroughly clean it separately. Use a damp cloth or sponge to clean the stainless steel parts on the balance. You can use any commercially available household cleaning agent that is suitable for use on stainless steel. Clean stainless steel surfaces by wiping them down. Then clean the weighing pan thoroughly, making sure to remove all residues. Use a damp cloth or sponge to wipe down the stainless steel parts on the balance again. Afterwards, allow the balance to dry. If desired, you can apply oil to the cleaned surfaces as additional protection.

- ⚠ Do not use stainless steel cleaning agents that contain soda lye (caustic), acetic acid, hydrochloric acid, sulphuric acid or citric acid. The use of scrubbing sponges made with steel wool is not permitted. Solvents are permitted only for cleaning stainless steel parts.

# Recycling

## Safety Inspection

If there is any indication that safe operation of the balance with the AC adapter is no longer warranted:

- Turn off the power and disconnect the power cord from an electrical outlet (mains supply) immediately
- > Lock the balance and AC adapter in a secure place to ensure that the equipment cannot be used during this time

Safe operation of the balance with the AC adapter is no longer ensured when:

- there is visible damage to the AC adapter or power cord
- the AC adapter no longer functions properly
- the AC adapter has been stored for a relatively long period under unfavorable conditions

In any of these cases, notify your KERN service center. Maintenance and repair work may only be performed by service technicians who are authorized by KERN and who

- have access to the required maintenance manuals
- have attended the relevant service training courses

To ensure safe shipment, your balance has been packaged to the extent necessary using environmentally friendly materials. After successful installation of the balance, you should return this packaging for recycling because it is a valuable source of secondary raw material.

For information on recycling old weighing equipment, contact your communal or municipal waste disposal center or local recycling depot.


# Overview

## Specifications

Model		CGB0.6K0.1	CGB1.2K0.2	CGB3K0.5	CGB6K1
Weighing capacity	kg	0.6	1.2	3	6
Readability	g	0.1	0.2	0.5	1
Counting resolution	g	0.01	0.02	0.05	0.1
Tare range (subtractive)	kg	0.6	1.2	3	6
Repeatability	≤± g	0.2	0.4	1.0	2.0
Linearity	≤± g	0.3	0.6	1.5	3.0
Operating temperature range		0 ... +40°C (32°F ... 104°F)			
Sensitivity drift within -10 ... +40°C	ppm/K	50			
Stabilization time (average)	s	1.5			
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels			
Display update rate (depends on the filter level selected)		0.1 – 0.4			
External calibration weight (of at least accuracy class...)	kg lb	0.5 (M1) 1 (M1)	1 (M1) 2 (M1)	2 (M1) 4 (M1)	5 (M1) 10 (M1)
Net weight, approx.	kg/lb	2.0/4.4			
AC power source/power requirements		via AC adapter 230 V or 115 V, +15% ... -20%			
Frequency		48 – 60 Hz			
Power source, constant voltage	V	10 to 20			
Power consumption (average)	W	0.75			
Hours of operation with 6 size C cells:					
– Alkaline manganese batteries, approx.:	h	114			
– Zinc carbon batteries, approx.:	h	41			
– Nickel metal hydride rechargeable batteries: fully charged, approx.:	h	30			
– Nickel cadmium rechargeable batteries: fully charged, approx.:	h	22			
Selectable weight units		Grams, kilograms, carats, pounds, ounces, Troy ounces, Hong Kong taels, Singapore taels, Taiwa- nese taels, grains, pennyweights, parts per pound, Chinese taels, mommes, Austrian carats, tola, baht, mesghal			
Built-in interface		RS-232C			
Format:		7-bit ASCII, 1 start bit, 1 or 2 stop bits			
Parity:		Mark, odd, even or space			
Transmission rates:		150 to 9,600 baud			
Handshake:		Software or hardware			

Model		CGB16K2	CGB34K5	CGB64K10
Weighing capacity	kg	16	34	64
Readability	g	2	5	10
Counting resolution	g	0.2	0.5	1
Tare range (subtractive)	kg	16	34	64
Repeatability	≤±g	4	10	20
Linearity	≤±g	6	15	30
Operating temperature range		0 ... +40°C (32°F ... 104°F)		
Sensitivity drift within -10 ... +40°C	ppm/K	50		
Stabilization time (average)	s	1.5		
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels		
Display update rate (depends on the filter level selected)		0.1 – 0.4		
External calibration weight (of at least accuracy class...)	kg lb	10 (M1) 20 (M1)	10 (M2) 20 (M2)	20 (M2) 40 (M2)
Net weight, approx.	kg	6.0		
AC power source/power requirements		Via AC adapter 230 V or 115 V, +15% ... -20%		
Frequency		48 – 60 Hz		
Power source, constant voltage	V	10 to 20		
Power consumption (average)	W	0.75		
Hours of operation with 6 size C cells:				
– Alkaline manganese batteries, approx.:	h	114		
– Zinc carbon batteries, approx.:	h	41		
– Nickel metal hydride rechargeable batteries; fully charged, approx.:	h	30		
– Nickel cadmium rechargeable batteries; fully charged, approx.:	h	22		
Selectable weight units		Grams, kilograms, carats, pounds, ounces, Troy ounces, Hong Kong taels, Singapore taels, Taiwa- nese taels, grains, pennyweights, parts per pound, Chinese taels, mommes, Austrian carats, tola, baht, mesghal		
Built-in interface		RS-232C		
Format:		7-bit ASCII, 1 start bit, 1 or 2 stop bits		
Parity:		Mark, odd, even or space		
Transmission rates:		150 to 9,600 baud		
Handshake:		Software or hardware		

# Dimensions (Scale Drawings)

Dimensions in millimeters

Model	a	b	c	d	e	f	g	h	i
CGB weighing capacity $\leq 6$ kg	82	10	24	149	329	202	252	252	199
CGB weighing capacity $\geq 16$ kg	121	12	24	238	429	302	252	402	338

Dimensions in inches

Model	a	b	c	d	e	f	g	h	i
CGB weighing capacity $\leq 6$ kg	3.2	0.4	0.9	5.9	13	8	9.9	9.9	7.8
CGB weighing capacity $\geq 16$ kg	4.8	0.5	0.9	9.4	16.9	11.9	9.9	15.8	13.3

# CE Marking

The equipment complies with the following EC Directives and European Standards:

## **Council Directive 89/336/EEC "Electromagnetic compatibility (EMC)"**

Applicable European Standards:

1. Electromagnetic Compatibility
- 1.1 Reference to 89/336/EEG: Official Journal of the European Communities, No. 2001/C105/03

EN 61326-1 Electrical equipment for measurement, control and laboratory use EMC requirements

Part 1: General requirements  
Defined immunity to interference:  
Industrial areas, continuous nonmonitored operation  
Limitation of emissions:  
Residential areas, Class B

Important Note:

The operator shall be responsible for any modifications to KERN equipment and must check and, if necessary, correct these modifications. On request, KERN will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

## **73/23/EEC "Electrical equipment designed for use within certain voltage limits"**

Applicable European Standards:

EN 60950

Safety of information technology equipment including electrical business equipment

EN 61010

Safety requirements for electrical equipment for measurement, control and laboratory use

Part 1: General requirements

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

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→

→

- Place empty container on scale
- Start counting

→ Select reference sample quantity

→ Place reference sample quantity on scale

→ Save reference sample quantity

→ Fill container

→ Select reference sample quantity

→ Place reference sample quantity on scale

→ Save reference sample quantity

→ Fill container

- Place full container on balance
- Start counting

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E-Mail: [info@kern-sohn.de](mailto:info@kern-sohn.de)

Internet: [www.kern-sohn.de](http://www.kern-sohn.de)

Copyright by Kern, Balingen, Germany.  
All rights reserved. No part of this publication  
may be reprinted or translated in any form or by any means  
without the prior written permission of Kern.

The status of the information, specifications and  
illustrations in this manual is indicated by the date  
given below. Kern reserves the right to  
make changes to the technology, features,  
specifications and design of the equipment  
without notice.

Status: April 2003, Kern, Balingen, Germany

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## **KERN & Sohn GmbH**