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Ziegelei 1  
D-72336 Balingen-Frommern  
Postfach 40 52  
72332 Balingen

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**Operating Instructions  
Display devices**

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### **CAUTION!**

Do not use KME-TM / KMN-TM in hazardous areas!

Our product range includes special devices for hazardous areas.

### **CAUTION!**

Use only weighing terminals with Protection Class IP65 if:

- the weighing terminal is used in wet areas
- wet cleaning is necessary
- the weighing terminal is used in a dusty environment

Even with Protection Class IP65, the weighing terminal must not be used in environments with corrosion risk.

- ▲ Never flood the weighing terminal or immerse it in liquid.

### **DANGER!**

Electric shock hazard!

- ▲ Always pull out the mains plug before commencing any work on the weighing terminal.

### **DANGER!**

Hazard of electric shock if the mains cable is damaged!

- ▲ Check the mains cable for damage regularly. Switch off the weighing terminal immediately if the mains cable is damaged.
- ▲ Maintain a clearance of at least 3 cm on the rear side of the weighing terminal in order to prevent the mains cable bending too much.

### **CAUTION!**

Do not at all open the weighing terminal!

The warranty is void if this stipulation is ignored. The weighing terminal may only be opened by authorized persons.

## Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of in domestic waste. This also applies to countries outside the EU as per their specific regulations.

- Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.


Thank you for your contribution to environmental protection.

If the weighing terminal has a rechargeable **battery**:

The battery contains heavy metals.

- Observe the local regulations on the disposal of environmentally hazardous materials.

## Display

- 1 6-digit weight display
- 2 Status indicators
- 3 Keypad

## Status indicators

LED	Meaning
<b>Under / OK / Over</b> alternatively <b>Count / PCS / APW</b>	Indicators for checkweighing Indicators for counting
~	Motion indicator
<b>Net</b>	The displayed weight value is a net weight value
<b>&gt;I1&lt; / &gt;I2&lt;</b>	Display of the current weighing range of the connected weighing platform.
<b>lb / kg</b>	Currently chosen weight unit
	Condition of the rechargeable battery

## Keys

Key	Operating mode	Menu	Key	Operating mode	Menu
	Switching on/off; abort	–		Function key	Back to the next higher menu item
	Zeroing	Scrolling back		Clear key	Back to the previous menu item
	Taring	Scrolling forward		Transfer key Long key-press: Calling up menu	Activating menu item Accepting selected setting

## 2 Putting into operation

### 2.1 Opening weighing terminal

#### CAUTION!

- ▲ Before opening the terminal, switch the terminal off and disconnect the power supply.

#### Opening KME-TM

- Unscrew 4 screws and lift off cover.

#### Opening KMN-TM

The cover of the KMN-TM is locked by 4 spring clips.

1. Insert the tip of a flat-blade screwdriver into one of the two slots located on the bottom of the cover, and gently push in towards the enclosure until a "pop" sound is heard.
2. Lift off cover.

### 2.2 Connecting weighing platform

1. Insert the weighing platform cable through the cable gland into the weighing terminal.
2. Connect the weighing platform cable to the 7-pin terminal strip J2 according to the following table.

Terminal	1	2	3	4	5	6	7
Assignment	+EXC	+SEN	+SIG	Shield	-SIG	-SEN	-EXC

- With 4-wire load cells make bridges with: +ECS and +SEN, -EXC and -SEN.

### 2.3 Connecting serial interface

#### KME-TM

With KME-TM, the interface connection is carried out as a 9-pin D-sub male connector.

- Plug in the 9-pin D-sub cable on the connector at the weighing terminal.

#### KMN-TM

With KMN-TM, the serial interface must be connected inside the weighing terminal.

1. Insert the interface cable through the cable gland into the weighing terminal.
2. Connect the interface cable to the 3-pin terminal strip J3 according to the following table.

Terminal	1	2	3
Assignment	TXD	RXD	GND

## 2.4 Closing the terminal

### Closing KME-TM

- Mount the cover and tighten the 4 screws.

### Closing KMN-TM

- Mount the cover and press on the bottom housing until the spring clips engage.

## 2.5 Connecting the power supply

### CAUTION!

Before connecting the weighing terminal to the mains, check whether the voltage value printed on the rating plate corresponds with the local mains voltage.

- ▲ Do not at all connect the weighing terminal if the voltage value printed on the rating plate is different to the local mains voltage.

- Plug the mains plug into the socket.  
After connection, the weighing terminal performs a self-test.  
When the zero display appears, the weighing terminal is ready for operation.
- Adjust the weighing platform in order to obtain the highest possible precision (see Service Manual).

## 2.6 Inserting/replacing battery

### NiMH rechargeable battery


1. Open the battery case on the bottom of the terminal.
2. Undo the old rechargeable battery from the velcro strip, and disconnect plug if necessary.
3. Connect the (new) rechargeable battery to the terminal, and insert the rechargeable battery in the battery case.
4. Close the battery case.

## 2.7 Notes on certified weighing systems


On certified weighing systems, the weighing platform connection at the weighing terminal must be sealed with a sealing wire or a sealing sticker.

Please contact your local dealer your local weights and measures authorities.

### KME-TM

### KMN-TM

## 3 Basic functions

### 3.1 Switching on and off

#### Switching on

→ Press  .

The display lights up and then shows the software number.

When the weight display appears, the weighing terminal is ready for operation.

#### Switching off


→ Press and hold  until **-OFF-** is displayed.

### 3.2 Zeroing

Zeroing corrects the influence of slight soiling on the load plate.

#### Setting to zero manually

1. Unload weighing platform.

2. Press  .

The zero display appears.

#### Automatic zeroing

In case of non-certified weighing platforms, the automatic zero point correction can be deactivated in the supervisor menu (F1.4.1).

As standard, the zero point of the weighing platform is automatically corrected when the weighing platform is unloaded.

### 3.3 Simple weighing

1. Place weighing sample on the weighing platform.

2. Wait until the motion indicator goes out.

3. Read weighing result.

### 3.4 Weighing with tare

#### Taring

→ Place the empty container on the weighing platform and press  .

The zero display and the **Net** indicator appear.

## Clearing the tare

→ Press  .

The **Net** indicator goes out, the gross weight appears in the display.

- If automatic clearing of the tare weight is set in the supervisor menu (F1 . 5 . 2=On), the tare weight is cleared automatically as soon as the weighing platform is unloaded.
- If tare interlock is set in the supervisor menu (F1 . 5 . 3=On), the tare weight can only be cleared when the weighing platform is unloaded.

## Automatic taring

This function must be activated in the supervisor menu (F1 . 5 . 1=On).

→ Place the empty container on the weighing platform.

The weight applied on the weighing platform is automatically saved as the tare weight.

The zero display and the **Net** indicator appear.

## Printing/transferring data

→ Press  .

The display contents are printed out or transferred to a computer.

## 3.5 Notes on battery operation

With a fully loaded new rechargeable battery, you dispose of approx. 35 operating hours.

The  indicator shows the state of the rechargeable battery.

 constantly red

Approx. 10 % of full voltage

 slowly blinking, red

Approx. 5 % of full voltage

 fast blinking, red

Less than 5 % of full voltage,  
battery has to be loaded immediately

 green

Battery is loaded.

Loading takes place as soon as the weighing terminal is  
connected to the mains.

## 3.6 Cleaning

### **DANGER!**

Electric shock hazard due to ingress of moisture!

- ▲ Before cleaning the weighing terminal, pull out the mains plug to disconnect the unit from the power supply.

## Further notes on cleaning

- Use a damp cloth.
- Do not use any acids, alkalis or strong solvents.

## KME-TM

- Do not clean the weighing terminal using a high-pressure cleaning unit or under running water.
- Follow all the relevant instructions regarding cleaning intervals and permissible cleaning agents.

## KMN-TM

- Do not clean the weighing terminal using a high-pressure cleaning unit.

## 4 Applications

Depending on the setting of F2.1 parameter in the operator menu, different applications can be activated using the **F** key.

### 4.1 Displaying weight values with a higher resolution

For this purpose F2.1=MULT must be set in the operator menu (factory setting).

→ Press **F**.

The weight value is displayed with a higher resolution for about 20 seconds.

#### Note

The weight value in higher resolution cannot be printed.

### 4.2 Switching weight unit

For this purpose, F2.1=Unit must be set in the operator menu.

→ Press **F**.

The weight value is displayed in the second weight unit.

#### Note

The displayed weight unit remains until it is switched again.

### 4.3 Check weighing

For this purpose, F2.1=OVER and F2.2.1=CHECK (factory setting) must be set in the operator menu. In the factory setting, the check weighing function is working with upper and lower tolerances of 10 d. With parameters F2.2.3 and F2.2.4, these tolerances can be customized.

#### Setting target weight


1. Press **[F]** to activate the check weighing function.
  2. Press and hold **[F]** until **tARGET** and the 3 indicators **Under**, **OK** and **Over** appear. If **F2.2.2=WEIGHT** (factory setting) is set in the operator menu, the weight display appears.
  3. Put the target weight on the weighing platform and save with **[F]**.  
The **OK** indicator lights.
- If **F2.2.2=MANUAL** is set in the operator menu, the weight display with blinking last digit appears.
3. Enter target weight using the **[TARE]**, **[→0←]** and **[F]** keys and confirm with **[PRINT]** (see page 15).
  4. Save entered weight value as target weight using the **[F]** key.

#### Check weighing

Example: Target weight = 1.000 kg

- Weight less than target weight and below lower tolerance.  
The **Under** indicator lights.

- Weight within tolerances.  
The **OK** indicator lights.

- Weight more than target weight and above upper tolerance.  
The **Over** indicator lights.



#### Switching between checkweighing and normal weighing

- Press **[F]** to switch between checkweighing and normal weighing.

## 4.4 Classifying




For this purpose, F2.1=OVER and F2.1.1=CLASS must be set in the operator menu. In the factory setting, the classifying function is working with upper and lower tolerances of 10 d. With parameters F2.2.3 and F2.2.4, these tolerances can be customized.

### Setting target weight

1. Press  to activate the classifying function.
  2. Press and hold  until **tARGET** and the 3 indicators **Under**, **OK** and **Over** appear. If **F2.2.2=WEIGHT** (factory setting) is set in the operator menu, the weight display appears.
  3. Put the target weight on the weighing platform and save with  .  
The **OK** indicator lights.
- If **F2.2.2=MANUAL** is set in the operator menu, the weight display with blinking last digit appears.
3. Enter target weight using the ,  and  keys and confirm with  (see page 15).
  4. Save entered weight value as target weight using the  key.

### Classifying

Example: Target weight = 1.000 kg

- Weight less than target weight and below lower tolerance.  
The **Under** indicator lights.
- Weight within tolerances.  
The **OK** indicator lights.
- Weight more than target weight and above upper tolerance.  
The **Over** indicator lights.

### Switching between classifying and normal weighing

- Press  to switch between classifying and normal weighing.

## 4.5 Counting

For this purpose,  $F2.1=Count$  must be set in the operator menu.

When the weighing terminal is used predominantly for counting operations, KERN & Sohn recommends to stick the delivered label (**Count, APW, PCS**) over the left 3 indicators (**Under, OK, Over**).

### Determining piece weight and counting


1. Press  to activate the counting function.  
The **Count** and **PCS** indicators light.
2. Press and hold  until **PCS ...** appears.
3. Press  repeatedly until the desired piece number (**5, 10, 20, 50**) for the determination of the piece weight appears.
4. Put the displayed number of pieces on the weighing platform and confirm with .  
The number of pieces is displayed and the **PCS** indicator lights.
5. Put additional pieces on the weighing platform.  
The current piece number is displayed.

### Switching between number of pieces and weight


Using the  key, you can switch between the following values:

- Display of the number of pieces.  
The **PCS** indicator lights.

- Display of the piece weight.  
The **APW** indicator lights.

- Display of the total weight.  
None of the 3 indicators for counting is lighting.

## 5 Operator menu

The operator menu consists of the following blocks:

F2 – F key menu settings


F3 – terminal menu settings

F4 – communication menu settings

F6 – exit menu

### 5.1 Entering the operator menu

→ In gross mode, press and hold  until **MAStEr** appears.

→ Enter pass word    and confirm with  .  
**SEtUP** appears.


→ Press  .  
**F2** appears.

### 5.2 Operating the menu

#### Keys and their function in the menu

-  Selecting next parameter.
-  Back to the previous parameter.
-  Confirming selection.
-  Back to the previous menu item.
-  Back to the next higher menu item.

#### Numeric entry

1. Press  for editing the displayed value.  
The (last) digit is blinking.
2. Increase the displayed digit using the  key.  
– or –  
Decrease the displayed digit using the  key.

3. When entering multi-digit numbers, use the **F** key to move the cursor one place to the left.
4. Change the digit as described in step 2.
5. Repeat steps 3 and 4 if necessary.
6. When all digits are entered, use the **PRINT** key to confirm the entry.

**Note**

With **C**, you can clear the entry.

### 5.3 F2 – F key menu

Factory settings are printed with **bold** characters.

#### F2.1 – Function of the F key

4 different functions can be assigned to the F key:

- MUL10** When pressing the F key, the weight value is displayed in **10 times higher resolution**
- Unit** When pressing the F key, the weight unit switches between kg and lb
- OVER** Plus/Minus weighing  
Additional settings, see F2.2
- Count** Counting  
Additional settings, see F2.3

#### F2.2 – Plus/Minus weighing

These parameters only appear if F2.1=OVER is set.

##### F2.2.1 – Operating mode

- CHECh** **Check weighing**
- CLASS** Classifying

##### F2.2.2 – Setting the target weight

- WEIGHT** By **weighing in**
- MANUAL** By numeric entry

### F2.2.3 – Upper tolerance

After selecting the parameter, the currently set upper tolerance is displayed.

1. If necessary, use the **F** key to activate editing.
2. Change tolerance using the **←0→**, **TARE** and **F** keys.

Factory setting      **10 d**  
Possible settings    0 ... full load

### F2.2.4 – Lower tolerance

After selecting the parameter, the currently set lower tolerance is displayed.

1. Use the **F** key to activate editing if necessary.
2. Change tolerance using the **←0→**, **TARE** and **F** keys.

Factory setting      **10 d**  
Possible settings    0 ... full load

### F2.3 – Reference optimization

This parameter only appears if **F2.1=Count** is set.

**OFF**      **No reference optimization**

**ON**      Reference optimization enabled. The weighing terminal automatically determines the piece weight again when the number of parts has increased.

### F2.10 – Reset F key settings

Reset all parameters **F2.x(.x)** to factory setting.

## 5.4 F3 – terminal menu

Factory settings are printed with **bold** characters.

### F3.1 – Display settings

#### F3.1.1 – Timeout

The weighing terminal switches back to the weighing mode if during the set time no action was in the menu.

Factory setting              **60 (seconds)**  
Function disabled          0  
Possible settings          10 ... 999 (seconds)

### F3.1.2 – Brightness with battery operation

`Lo`      **low brightness**

`MEd`      high brightness

To save battery power, we recommend the `Lo` setting.

### F3.2 – Auto power off

The weighing terminal is switched off if during the set time no action was on the weighing terminal or on the weighing platform.

Factory setting      **5 (minutes)**

Function disabled      0

Possible settings      0.5 ... 60 (minutes)

### F3.3 – Battery type

This menu item is only available with weighing terminals in battery operation.

`drY`      **dry battery**

`ni-MH`      NiMH rechargeable battery

`LEAd-A`      lead-acid rechargeable battery

### F3.10 – Reset terminal settings

Reset all parameters F3.x(x) to factory setting.

## 5.5 F4 – communication menu

Factory settings are printed in **bold** characters.

### F4.1 – Connections

`Print`      When **pressing**  , the current display is printed

`APrint`      Stable weight values are printed automatically  
Additional settings: F4.2.5 and F4.2.6

`SICS`      Communication via the METTLER TOLEDO **S**tandard **I**nterface **C**ommand **S**et

`Contin`      Toledo continuous Mode

### F4.2 – Format

#### F4.2.1 – Line format

`MULTi`      **Multi line**

`SINGLE`      Single line

#### **F4.2.2 – Print format**

**stAndr Standard** (current display)

OVEr over / good / under

Count Piece number

#### **F4.2.3 – Print language**

**EnG English**

CHn Chinese

#### **F4.2.4 – Add line feed**

Factory setting **3 (lines)**

Possible settings 0 ... 9 (lines)

#### **F4.2.5 – Auto print threshold**

This menu item is only available if F4.1=APrint is set.

A stable weight value which is higher than the set value is printed automatically.

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Possible settings 0 ... max. load

#### **F4.2.6 – Auto print reset threshold**

This menu item is only available if F4.1=APrint is set.

The scale must be unloaded below the set value before a new weight value can be printed automatically.

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Possible settings 0 ... max. load

#### **F4.3 – Parameters**

##### **F4.3.1 – Baudrate**

1200

2400

4800

**9600**

19200

### **F4.3.2 – Data bits / parity**

7-odd 7 bits, parity odd

7-even 7 bits, parity even

8-none **8 bits, no parity**

### **F4.3.3 – Xon/Xoff**

On Xon/Xoff enabled

**OFF** Xon/Xoff **disabled**

### **F4.3.4 – Checksum**

On Checksum enabled

**OFF** Checksum **disabled**

### **F4.10 – Reset communication settings**

Reset all parameters F4.x(x) to factory setting.

## **5.6 F6 – ending menu**

1. Press  .

**F6** appears.

2. To save changes: Press  .

**SAVE ?** appears.

Then press  again.

– or –

To reject changes: Press  .

**AbOrt** appears.

Press  .

## 6 Supervisor menu

The supervisor menu consists of the following blocks:

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F5 – terminal settings

F6 – exit menu

### 6.1 Entering Supervisor menu

- In gross mode, press and hold  until **MAStEr** appears in the display.
- Enter password     and confirm with  .  
**SEtUP** appears in the display.
- Press  .  
**F1** appears in the display. All parameters can be modified.

#### Note on certified weighing systems (OIML or NTEP)

On certified weighing systems, the parameters F1, F5.1 and F5.4 are locked.

In order to change these parameters, proceed as follows:

1. Switch off the weighing terminal and open it.
2. Press button S1 on the mainboard and switch on weighing terminal.  
**SEtUP** is displayed, and all parameters can be modified.
3. When configuration is finished, seal the weighing terminal.

### 6.2 Operating the Supervisor menu

Operating the Supervisor menu means the same as operating the Operator menu, see page 15.

### 6.3 Block F1 – Scale

Factory settings are printed in **bold** letters.

#### F1.1 – Approval

**no**      **no approval**

**OIML**    approval according to OIML

**nTEP**    approval according to NTEP

**otHEr**    for other approvals

### **F1.2.1 – Weight units**

- 1 weight unit: kg
- 2 weight unit: lb 1 lb  $\approx$  0.454 kg

### **F1.2.2 – Weighing ranges**

- 1  $r$  **for single range weighing platforms**
- 2  $r$  for double range weighing platforms

### **F1.2.3 – Capacity of the first weighing range (coarse range)**

Possible capacities and the factory settings depend on the weighing platform connected.

→ If necessary, modify the displayed value.

### **F1.2.4 – Resolution of the first weighing range (coarse range)**

Possible resolutions and the factory settings depend on the weighing platform connected.

→ If necessary, modify the displayed value.

### **F1.2.5 – Capacity of the second weighing range (fine range)**

Possible capacities and the factory settings depend on the weighing platform connected.

This parameter only appears if  $F1.2.2=2r$  is set.

→ If necessary, modify the displayed value.

### **F1.2.6 – Resolution of the second weighing range (fine range)**

Possible resolutions and the factory settings depend on the weighing platform connected.

This parameter only appears if  $F1.2.2=2r$  is set.

→ If necessary, modify the displayed value.

### **F1.3.1 – Geo value**

Adaptation of the weighing platform to the geographical location, see table in the annex.

Possible settings 0 ... 31



### **F1.3.2 – Linearization during adjustment**

**LinOFF** Linearization disabled

LinOn Linearization enabled

### F1.3.3 – Adjustment

The steps with grey background only appear if parameter F1.3.2=LinOn is set.

Display	Key	Description
<b>E SCL</b>		Unload weighing platform
		Confirm empty weighing platform
<b>10 CAL</b> ... <b>0 CAL</b>		The weighing terminal counts down from 10 to 0. Zero is determined
<b>Add Ld</b>		Load half maximum load
		Confirm half load
<b>000000</b>		Enter weight value for half maximum load
	  	Enter weight value
<b>003000</b>		Weight value for half maximum load entered
		Confirm weight value
<b>10 CAL</b> ... <b>0 CAL</b>		The weighing terminal counts down from 10 to 0. Half maximum load is adjusted
<b>FULL Ld</b>		Load maximum load
		Confirm maximum load
<b>000000</b>		Enter weight value of maximum load
	  	Enter weight value
<b>006000</b>		Weight value for maximum load entered
		Confirm weight value
<b>10 CAL</b> ... <b>0 CAL</b>		The weighing terminal counts down from 10 to 0. Maximum load is adjusted
<b>done</b>		Adjustment finished. This message is displayed for about 2 seconds
<b>F1.4</b>		Next block in the supervisor menu

#### **F1.4.1 – Automatic zero setting**

- OFF** Automatic zero setting disabled (not available in NTEP-Mode)
- 0.5 d** Automatic zero setting within **+/-0.5 d**
- 1 d** Automatic zero setting within +/-1.0 d (not available in OIML-Mode)
- 3 d** Automatic zero setting within +/-3 d (not available in OIML-Mode)

#### **F1.4.2 – Power up zero**

- OFF** Power up zero disabled
- 2** Power up zero within +/-2 %
- 10** Power up zero within **+/-10 %**
- 20** Power up zero within +/-20 % (not available in OIML and NTEP-Mode)

#### **F1.4.3 – Pushbutton zero**

- 0** Pushbutton zero disabled
- 2** Pushbutton zero with **+/-2 % zero setting range**
- 10** Pushbutton zero with +/-10 % zero setting range (not available in OIML and NTEP-Mode)
- 20** Pushbutton zero with +/-20 % zero setting range (not available in OIML and NTEP-Mode)

#### **F1.5.1 – Automatic taring**

- On** Automatic taring enabled
- OFF** Automatic taring **disabled**

#### **F1.5.2 – Auto clear tare**

- On** Clearing tare automatically enabled
- OFF** Clearing tare automatically **disabled**

#### **F1.5.3 – Tare Interlock**

- On** The weighing platform must be unloaded to zero before the tare weight can be cleared.
- OFF** Function **disabled**

#### **F1.5.4 – Auto tare threshold**

This menu item is only available if F1.5.1=On is set.

The weighing platform must be loaded to the set value before the weight value is automatically tared.

- Factory setting **10 d**
- Possible settings 0 ... maximum load

### **F1.5.5 – Auto clear tare threshold**

This menu item is only available if F1.5.1=On is set.

The weighing platform must be unloaded below the set value before a new weight value can be tared automatically.

If F1.5.2=On is set, the weighing platform must be unloaded to the set value before the tare value is cleared automatically.

Factory setting      **10 d**

Possible settings    0 ... maximum load

### **F1.6.1 – Digital filter**

The digital filter stabilize+s the weight display when the load is moving or vibrating.

Lo            Low filter

**MEd**        **Medium** filter

HIGH        High filter

### **F1.6.2 – Motion detection**

#### **0.5 d    Motion detection within +/-0.5 d**

1 d        Motion detection within +/-1 d (not available in OIML and NTEP-Mode)

3 d        Motion detection within +/-3 d (not available in OIML and NTEP-Mode)

### **F1.10 – Resetting parameters 1.x(x) to factory setting**

Reset all F1 parameters to factory setting, except for calibration values.

In OIML and NTEP mode the Geo value is not reset.

## 6.4 Block F5 – Maintenance

Factory settings are printed in **bold** letters.

### F5.1 – Calibration values

#### F5.1.1 – Show zero-counts

#### F5.1.2 – Show half load weight value

#### F5.1.3 – Show half load counts

#### F5.1.4 – Show full load weight value

#### F5.1.5 – Show full load weight counts

### F5.2– Keypad test

The terminal shows **PrESS**.

→ Press      .

→ Press  to exit keypad test.

### F5.3 – Display test

All display segments light up.

### F5.4 – Display internal resolution

### F5.5 – COM1 test

To do so, the terminal must be connected to a computer.

### F5.6 – Print setup

Print all setup parameters

### F5.10 – General reset

Reset all parameters of groups F1 to F4 to factory settings.

## 7 Error messages

Error code	Error	Remedy
Err 3	EEPROM error	→ Switch weighing terminal off and on
Err 4	Number of reference parts too small	→ Put on additional reference parts
Err 6	EEPROM read/write error	→ Contact your dealer or local representative.
Err 35	Weighing platform in motion when calibrating	→ Ensure that the weighing platform is quiet
Err 70	Keypad error	→ Contact your dealer or local representative.
	Underload	→ Unload weighing platform → Press  → If the message reappears, contact your dealer or local representative.
	Overload	→ Decrease load
	Zero setting outside zero setting range	→ Unload weighing platform
Weighing terminal switches off automatically	<ul style="list-style-type: none"> <li>• Automatic power off enabled</li> <li>• Battery voltage too low</li> </ul>	→ Switch on weighing terminal → Load battery
Weighing terminal cannot be switched on	Fuse blown	→ Contact your dealer or local representative.



## 8 Technical data

### 8.1 General technical data



Display	<ul style="list-style-type: none"> <li>• weight value: 7-digit display, 6 numbers, 30 mm high</li> <li>• Status indicator: 10 indicators</li> </ul>
Operating temperature	• -10 ... +40 °C
Storage temperature	• -20 ... +60 °C
Relative humidity	• 10 ... 85 %, non condensing
Weight (incl. packing)	<ul style="list-style-type: none"> <li>• KME-TM: approx. 1.3 kg / 2.9 lb</li> <li>• KMN-TM: approx. 2.9 kg / 6.4 lb</li> </ul>

### Dimensions

#### KME-TM

#### KMN-TM

## 8.2 Applications

Weighing functions	<ul style="list-style-type: none"> <li>• Zero setting, taring, clearing tare</li> <li>• Unit switching</li> <li>• Display in higher resolution</li> <li>• Simple check weighing/classifying</li> <li>• Counting</li> </ul>
Additional functions	<ul style="list-style-type: none"> <li>• Printing</li> <li>• Automatic printout</li> <li>• English or Chinese printout formats selectable</li> <li>• Strip printer PQ16 supported</li> <li>• Energy saving technology, battery indicator</li> <li>• Automatic power off</li> </ul>

## 8.3 Technical data A/D converters


Supply	<ul style="list-style-type: none"> <li>• +5 VDC</li> </ul>
Connectable weighing cells	<ul style="list-style-type: none"> <li>• Max. 4 350 <math>\Omega</math> weighing cells</li> </ul>
Resolution	<ul style="list-style-type: none"> <li>• Max. 30.000 d</li> <li>• Internal max. 1.000.000 d</li> </ul>
Update rate	<ul style="list-style-type: none"> <li>• 30/s</li> </ul>
Signal input	<ul style="list-style-type: none"> <li>• 0 ... 5 mV (Zero)</li> <li>• 1 ... 10 mV (Span)</li> </ul>

### Terminal assignment for load cell connection

Terminal	1	2	3	4	5	6	7
Assignment	+EXC	+SEN	+SIG	Shield	-SIG	-SEN	-EXC

## 8.4 Serial interface data

### Assignment of the serial interface connection (KME-TM)

Pin	1	2	3	4	5	6	7
Assignment	-	RXD	TXD	-	GND	-	-

For KMN-TM, refer to page 6.

## SICS interface commands

The weighing terminal supports the MT-SICS (METTLER TOLEDO **S**tandard **I**nterface **C**ommand **S**et) command set. With SICS commands, it is possible to configure, query and operate the terminal from a PC. SICS commands are divided up into various levels.

For further information about the MT-SICS command set, see MT-SICS Manual (Order No. 22 011 459) or contact the KERN & Sohn Customer Service.

	<b>Command</b>	<b>Meaning</b>
LEVEL 0	@	Reset the scale
	I0	Inquiry of all available SICS commands
	I1	Inquiry of SICS level und SICS version
	I2	Inquiry of scale data
	I3	Inquiry of scale software version
	I4	Inquiry of serial number
	S	Send stable weight value
	SI	Send weight value immediately
	SIR	Send weight value immediately and repeatedly
	Z	Zero the scale
	ZI	Zero immediately
LEVEL 1	T	Tare
	TAC	Clear tare
	TI	Tare immediately

## Toledo Continuous commands

<b>Command</b>	<b>Meaning</b>
<b>P</b>	Print out the current result
<b>T</b>	Tare the scale
<b>Z</b>	Zero the display
<b>C</b>	Clear the current value
<b>S</b>	Define reference piece number

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## Notice

Certi ed balances and balances used for legal applications have the EU type approval. The year of the initial verification is shown next to the CE mark. Such balances are verified in the factory and carry the „M“ mark on the actual balance and the packaging. The year of initial verification is shown next to the CE mark. The GEO value of verified balances explains for which location of use the balance has been verified. This GEO value is shown on the balance itself and on the packing. Further details see GEO value table.

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