D	Montageanleitung Wägebrücken	Seite 2
E	Installation Instructions Weighing Platforms	Page 2
F	Notice de Montage Plate-formes de Pesage	Page 3
ESP	Instrucciones de Montaje Plataformas de Pesada	Página 3
I	Istruzioni di Montaggio Basamenti di Pesata	Pagina 3

KERN KTP

Version 1.0 04/12

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1. Important information

Select the **location of the weighing platform** using the following criteria:

- Never use weighing platforms in an environment subject to explosion hazards or in wet areas.
- Firm, vibration-free position as horizontal as possible.
- Temperature range from -10°C to $+40^{\circ}\text{C}$.
- No powerful drafts.

Never open the weighing platform, it contains no parts which can be serviced, repaired or changed by the user. On no account attempt to clean the inside of the weighing platform with solid objects.

2. Attaching the weighing platform

- First ensure that the terminal is disconnected from the power supply.

The connection cable of the weighing platform contains 6 lines plus shield:

Signal	Colour code
+OUT (2mV/V)	red
-OUT (2mV/V)	white
+IN (0-15V)	green
+SENSE	blue
-SENSE	brown
-IN (0-15V)	black
SHIELD	shield

- Connect the signal lines and the shield to the corresponding terminals of your weighing terminal. Consult the instructions enclosed with your terminal.

3. Levelling the weighing platform




The platform must be levelled on initial installation and after every location change:

- Turn the levelling feet of the weighing platform until the air bubble is in the centre of the level.
- To safeguard this adjustment setting, tighten the lock nuts of the levelling feet using the wrench supplied.

4. Matching the weighing platform to your operating conditions

In the factory the level is mounted on the longer side (A) of the weighing platform, however it can also be attached to the shorter side if need be:

- Turn the weighing platform over.
- Undo the two fastener screws of the level.
- Screw the level to the position provided on the shorter side of the weighing platform. The illustration opposite shows (when viewed from above) the original position "A" of the level and the alternative fastening possibility "B".

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