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# KERN MRS A01N

Version 1.0 11/2011

## Operating Instructions Temperature calibration set

---

---

### Calibrating/adjusting the temperature

We recommend periodic checks of device temperature using the optional temperature calibrating set MRS-A01N. Temperature is measured in two points (default setting 100°C and 160°C) – correction is possible in both of them. Earlier, after the heating phase is complete, the device should be left too cool to surrounding temperature.

#### 1. Calibrating temperature in Configuration menu

- ⇒ In order to enable temperature calibration, the option “Temp. CAL.” Must be set to “ON” (see chapter. 4.4.6 „Temperature calibration”):
- ⇒ Turn on the humidity meter pressing «**ON/OFF**».
- ⇒ When turning the device on hold «**MENU**» (for about 10 seconds), until information “SET DATA PRINT” is displayed erscheint.
- ⇒ Use the «» button to select the “Temp. CAL” option. Current “ON” or “OFF” setting will be displayed.

• **TEMP. CAL.**

TEMP. CAL. ON/OFF

Temperature calibration activation

In order to change the setting from “OFF” to “ON”, use the following settings:

TEMP. CAL. OFF

- If the setting is “OFF”, press the «» button, the “OFF” indication will blink.

TEMP. CAL. ON

- If the setting is “ON”, press the «» button, the “ON” indication will blink.

TEMP. CAL. ON

- Press the «» button to confirm
- Press «**esc**» to leave the menu

## 1. Temperature calibration/adjustment

- Make sure that drying is not performed – that is that the device is in weighing mode.
- Insert the temperature sensor bar into the sample chamber. Close the sample chamber.
- Press «**TARE**» until „TEMP. CALIBRATION” setting is shown.

TEMP. CALIBRATION

T100 T 160 START

- Calibrate according to the default values to temperatures of 100°C and 160°C. To set temperature calibration value to individual demands, use the  (100°C) button or the  (160°C) button. Use  and  to set the requested temperature. Use the  button to pass to the next number and enter the selected value. Then confirm the entered values using the «**START STOP**» button.

PRESS START KEY

TEMP 25°C 20.00MIN

TEMP 100°C 20.00MIN

- To start temperature calibration press «**START STOP**».
- Temperature calibration has been activated. The dryer is heated to the first temperature value. The temperature and the remaining temperature calibration time are displayed in the information strip. Heating to the first temperature value takes approx. 20 minutes.
- After the calibration temperature is reached, the indication blinks. In order to adjust, read the current temperature using a thermometer and enter it. Select the digit using the  button and increase/decrease the value using the   buttons. Afterwards press «**START STOP**».
- The device will then be heated to the second temperature value.
- After the calibration temperature is reached, the indicator blinks again. Read the current temperature using a thermometer, enter it and confirm it using the «**START STOP**» button.
- If an optional printer is connected, a protocol will be printed after the calibration is complete .

## Temperature calibration report printout

Temperature Calibration -----	Temperature calibration report
Date 16.10.2003 Time 12:51:36 Name : MRS 120-3 Software : N50-0000 P00 Serialno : 3300-1	Time of the calibration and instrument data
Temp. Reference ID : ...	Identifier of the temperature adjustment set
Temperature 100 C : 100 C Temperature 160 C : 160 C	Status of the temperature calibration
Temperature Calibration o.k.	Instrument has been correctly temperature-calibrated
Operator :	Operator ID, if activated under Set Print Format.

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