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Ziegelei 1
D-72336 Balingen
E-Mail: info@kern-sohn.com

Tel: +49-[0]7433- 9933-0
Fax: +49-[0]7433-9933-149
Internet: www.kern-sohn.com

**Zusatzbeschreibung Justierung geeichter Wägesysteme
Supplementary instructions verified weighing systems
Instructions supplémentaires systèmes de pesage étalonnés**

KERN KFB-TM

Version 2.0
09/2011

Deutsch

English

Français

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KERN KFB-TM

Version 2.0 09/2011

Supplementary instructions verified weighing systems


1 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each display unit with connected weighing plate must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the weighing system has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the display unit periodically in weighing operation.

- The adjustment is locked for verified weighing systems. To disable the access lock, destroy the seal and actuate the adjustment switch.
Attention:
After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.
- The applied adjustment weight depends on the weighing system's capacity. Carry out adjustment as near as possible to the maximum weight of the weighing system. Info about test weights can be found on the Internet at: <http://www.kern-sohn.com>
- Observe stable environmental conditions. Warm-up time required for stabilisation.

Invoke menu

<p>⇒ Switch-on balance and during the selftest press .</p>	
<p>⇒ Press subsequently , , , the first menu block „PO CHK“ will be displayed.</p>	
<p>⇒ Press  repeatedly until „P2 mode“ will be displayed. ⇒ Press the adjustment switch.</p>	
<p>⇒ Confirm by .</p> <p>⇒ Press  repeatedly until the currently balance typ will be displayed.</p> <p><i>S₁Gr</i> = Singel range <i>dUAL 1</i> = Multi range <i>dUAL 2</i> = Multi interval</p>	    
<p>⇒ Confirm by .</p>	
<p>⇒ Press  repeatedly until „CAL“ will be displayed.</p>	
<p>⇒ Confirm by .</p> <p>⇒ Press  until <i>nOLin</i> will be displayed.</p>	

How to carry out adjustments

<p>⇒ Confirm by . Ensure that there are no objects on the weighing plate.</p>	  
<p>⇒ Wait for stability display, then press .</p>	
<p>⇒ The currently set adjustment weight will be displayed.</p>	
<p>⇒ To change by using the navigation buttons (see chap. 2) select the desired setting, the active digit is flashing.</p> <p>⇒ Acknowledge with .</p>	
<p>⇒ Carefully place adjusting weight in the centre of the weighing plate. Wait for stability display, then press .</p>	
<p>⇒ After the adjustment the balance will carry out a self-test. Remove adjusting weight during selftest, balance will return into weighing mode automatically. An adjusting error or incorrect adjusting weight will be indicated by the error message; repeat adjustment procedure.</p>	

2 Numeric input via navigation keys

<p>⇒ Press ; current setting appears. The first digit is flashing and can be changed.</p>
<p>⇒ If the first digit is not to be changed, press  and the second digit will start flashing. Each time you press , the display unit jumps to the subsequent digit, returning to the first digit after the last digit has been pressed.</p>
<p>⇒ To change the selected (flashing) digit, press  repeatedly until the desired value appears. Then select by  using <i>additional digits</i> and change these by using .</p>
<p>⇒ Finish entry with .</p>

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