
Betriebsanleitung / Operating Instructions / Mode d'emploi /  
Instrucciones de funcionamiento / Manuale d'istruzioni

## KERN AGB/ PGB

Version 2.1 12/2005

Elektronische Analysen- und Präzisionswaagen/ Electronic Analytical  
and Precision Balances/ Balances électroniques d'analyse et  
de précision/ Balanzas electrónicas de precisión y analíticas/  
Bilance elettroniche analitiche e di precisione

3	Deutsch
32	English
61	Français
90	Español
119	Italiano

# Inhalt

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The non-automatic weighing instruments with the EC type-approval certificate no.:

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correspond to the production model described in the EC type-approval certificate and to the requirements of the Council Directive 90/384/EEC as amended.

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The weighing instruments were initially verified under number 0111 by the manufacturer at the place of manufacture.

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# Contents

32	Contents
32	Warnings and Safety Precautions
33	Initial Verification
34	Getting Started
38	Operation
38	Basic Weighing Function
39	Calibration/Adjustment
41	Application Programs
41	Net-Total Formulation/ Second Tare Memory
42	Counting
43	Weighing in Percent
44	Weigh Averaging
45	Toggling between Weight Units
46	Configuring the Balance (Menu Codes)
46	Setting the Parameters (Menu Codes)
47	Balance Operating Menu (Overview)
49	ISO/GLP-compliant Printout
51	Data Interface
52	Error Codes
53	Care and Maintenance
54	Instructions for Recycling
55	Overview
55	Specifications
59	Declarations of Conformity

# Warnings and Safety Precautions

## Safety

- To prevent damage to the equipment, please read these operating instructions carefully before using your balance.
- △ Do not use this balance in a hazardous area/location.
- △ Use only a commercially available non-rechargeable or rechargeable 9-volt battery
- △ Make absolutely sure to unplug the balance from AC power before you connect or disconnect a peripheral device.
- △ Exposure to excessive electromagnetic disturbance can cause the readout value to change. Once the disturbance has ceased, the instrument can be used again in accordance with its intended use.

## Setting up the Balance

- △ Warning when using pre-wired RS-232 connecting cables: The pin assignments in RS-232 cables purchased from other manufacturers may be incompatible with KERN balances. Be sure to check the pin assignment against the chart on page 51 before connecting the cable, and disconnect any lines that do not match.
- Connect only KERN accessories and options, as these are optimally designed for use with your KERN balance. Do not try to solve any problems on your own. The operator shall be responsible for any modifications to KERN equipment and for any connections of cables or equipment not supplied by KERN and must check and, if necessary, correct these modifications and connections. On request, KERN will provide information on the minimum operating specifications (in accordance with the standards for defined immunity to interference).
- Do not open the balance housing. If the seal is broken, this will result in forfeiture of all claims under the manufacturer's warranty.

# Initial Verification

## Model: PGB510-2M

- These models are verified ex factory for use in legal metrology within the European Economic Area upon delivery.
- Before using the balance as a legal measuring instrument, calibrate and adjust it at the place of use using the built-in motorized calibration weight; for details, see "Calibration/Adjustment."
- The temperature range (°C) indicated on the verification label may not be exceeded during operation.

Example:  
DT BH 210  
+10°C/+30°C

Ⓜ

## Using Verified Balances as Legal Measuring Instruments in the EU\*:

The type-approval certificate for verification applies only to non-automatic weighing instruments. For automatic operation with or without auxiliary measuring devices, you must comply with the regulations applicable to the place of installation.

- \* Including the Signatories of the Agreement on the European Economic Area

# Getting Started

## Storage and Shipping Conditions

- Do not expose the balance to extreme temperatures, blows, shocks, vibration or moisture.

## Unpacking the Balance

- After unpacking the balance, check it immediately for any visible damage
- If you see any sign of damage, proceed as directed in the chapter entitled "Care and Maintenance," under the section on "Safety Inspection"
- Save the box and all parts of the packaging until you have successfully installed your balance in case you need to return it. Before packing your balance, unplug all connected cables to prevent damage.

## Equipment Supplied

- Balance
- Weighing pan
- Pan support (only on models with a round weighing pan)
- AC adapter, plug type

Additionally supplied with AGB210-4, AGB120-4 models:

- Shield ring
- Shield plate
- Dust cover

Additionally supplied with PGB310-3, PGB150-3 models:

- Glass draft shield with cover
- PGB150-3: Test weight (100 g, F2)
- PGB310-3: Test weight (200 g, F1)
- PGB1500-2: Test weight (1,000 g, F2)
- PGB3100-2: Test weight (2,000 g, F1)

## Installation Instructions

When choosing a location to set up your balance, observe the following:

- Avoid placing the balance in close proximity to a heater or otherwise exposing the balance to heat or direct sunlight
- Protect the balance from drafts that come from open windows or doors
- Avoid exposing the balance to extreme vibrations during weighing
- Do not expose the balance to extreme moisture over long periods


## Conditioning the Balance

Moisture in the air can condense on the surfaces of a cold balance whenever it is brought into a substantially warmer place. If you transfer the balance to a warmer area, make sure to condition it for about 2 hours at room temperature, leaving it unplugged from AC power.

## Seal on Balances Verified for Use in Legal Metrology in the EU\*:

EU legislation requires that a control seal be affixed to the verified balance. The control seal consists of a sticker with the "KERN" logo. This seal will be irreparably damaged if you attempt to remove it. If the seal is broken, the validity of the verification becomes null and void, and you must have your balance re-verified.


- \* Including the Signatories of the Agreement on the European Economic Area

### Setting up the Balance


Balances with an Analytical Draft Shield

- Place the components listed below inside the chamber in the order given:
  - Shield plate
  - Shield ring
  - Pan support
  - Weighing pan

Balances with a Glass Draft Shield

- Place the components listed below inside the chamber in the order given:
  - Draft shield base – place it on the balance so that the edge for fitting the glass draft shield faces upwards and turn it until it is firmly in place
  - Pan support
  - Weighing pan
  - Glass draft shield
  - Draft shield cover – place it on the balance so that the edge faces downwards


Balances with a Round Weighing Pan

- Place the components listed below inside the chamber in the order given:
  - Pan support
  - Weighing pan

Balances with a Rectangular Weighing Pan

- Place the weighing pan on the balance


### Connecting the Balance to AC Power/ Safety Precautions

Use only original AC adapters


- Insert the right-angle plug into the jack
- The AC adapter rated to Class 2 can be plugged into any wall outlet without requiring any additional safety precautions

The ground is connected to the balance housing, which can be additionally grounded for operation.


### Using a Non-Rechargeable/Rechargeable Battery (only for models PGB610-2 and PGB6100-1)

- A non-rechargeable or rechargeable battery is not included with the equipment supplied
- ⚠ Use only a commercially available non-rechargeable or rechargeable 9-volt battery
- ⚠ When using a rechargeable battery, always use an external charger to recharge the battery
- Lay the balance on its side
- Open the battery compartment:  
Lift the compartment cover

- Insert the 9-volt (rechargeable) battery in the compartment
- Make sure to connect the positive and negative poles correctly
- ⚠ All used batteries are classified as waste that requires special handling (not "household" waste). Dispose of rechargeable batteries in accordance with the applicable special waste disposal regulations.
- Close the battery compartment:  
Press down on the cover until it clicks into place


### Leveling the Balance

(not for models PGB610-2 and PGB6100-1)

Level the balance any time you set it up in a new location. Use only the 2 front feet of the balance for leveling.

- Turn the 2 rear feet until they are in position (only on models PGB3100-2 and PGB1500-2)
  - Turn the 2 front feet as shown here in the illustration until the air bubble is centered in the level indicator
- > In most cases, this will require several adjustment steps

### Anti-theft Locking Device

To protect against theft, use the mounting lug on the rear panel of the balance.

- Secure the balance at the place of installation, for example with a chain or a lock

# Operation


## Basic Weighing Function

### Preparation

- Turn on the balance: Press [ON/OFF]
- To change configurations: See the chapter entitled “Configuring the Balance”
- To tare the balance: Press [TARE]

### Additional Functions

- To turn off the balance: Press [ON/OFF]

### Warmup Time

- To ensure accurate results, the balance must warm up for 30 minutes before operation. Only after this time will the balance have reached the required operating temperature.

### PGB510-2M

- To ensure accurate results, the balance must warm up for at least 24 hours before operation the first time it is connected to power. Only after this time will the balance have reached the required operating temperature.

### Example

#### Basic weighing

Step	Key (or instruction)	Display/Printout
1. Turn on the balance Self-test is performed	[ON/OFF]	
2. Place container on balance (here: 52 g)		+ 52.0 g
3. Tare the balance	[TARE]	+ 0.0 g
4. Place sample in container on balance (here: 150.2 g)		+ 150.2 g

# Calibration/Adjustment

## Available Features

Calibration/adjustment can only be performed when

- there is no load on the balance,
- the balance is tared,
- the internal signal is stable.

If these conditions are not met, an error message is displayed.

Otherwise, the weight required for calibration/adjustment is displayed.

It is possible to use the test weight provided to calibrate/adjust the balance; however, this is not optimal for some models from a metrological point of view.

### We recommend:

	PGB150-3	PGB310-3	PGB 1500-2	PGB 3100-2
Test weight included	100 g	200 g	1,000 g	2,000 g
Calibration/adjustment tolerance	F2	F1	F2	F1
Recommended calibration/ adjustment weight acc. to class	100 g F1	200 g E2	1,000 g F1	2,000 g E2

You can use any of the following weight units to calibrate/adjust:

g, kg\*, lb (menu code 1.4.x)



You can block calibration/adjustment of the balance:

- Select menu code 1.5.3)

Calibration/adjustment of verified balances of accuracy class  $\text{\textcircled{II}}$ :

- Calibration/adjustment is blocked (protective cap over menu access switch is sealed)
- > Calibration/adjustment is possible only after the seal has been removed. In this case, the verification becomes null and void and the balance must be re-verified.

## External calibration/adjustment

Step	Key (or instruction)	Display/Printout
1. Tare the balance		0.0 g
2. Begin calibration Calibration weight is displayed without weight unit (here: 1000 g)	[CAL]	+ 1000.0
3. Place the indicated calibration weight on the balance		1000.0
After calibration, the calibration weight is displayed with wt. unit		+ 1000.0 g
4. Remove the calibration weight		0.0 g

\* = not on models with a readability of 0.1 mg

**Calibration/adjustment with internal weight: Model PGB510-2M**

Menu setting 1 5 2 required.

Step	Key (or instruction)	Display/Printout
1. Calibrate/adjust the unloaded balance		
2. Tare the unloaded balance	TARE	0.0 g
3. Start calibration/adjustment	[CAL]	☼
4. Move the selector for the motorized weight into the "calibration/adjustment" position		☼
5. Calibration/adjustment is performed ("CC" flashes on the display)		☼☼
6. Remove weight: Return the selector to the initial position		☼☼
7. The internal weight is removed from the weighing system		0.0 g

# Application Programs

All application programs can also be used in legal metrology. Calculated weight values are indicated by one of the following symbols displayed to the right of the numerical value:

- Percentage = %
- Piece count = pcs
- Other calculated value = o

The application programs for calibrated balances not intended for use in legal metrology are labeled with a triangle  $\triangle$  in the right side of the display (e.g. averaging).




## Net-Total Formulation/Second Tare Memory

With this application program you can weigh in components for formulation of a mixture.

### Preparation

Configure the Net-Total Formulation/Second Tare Memory application in the operating menu: See "Configuring the Balance." Menu code: **2.13**

### Example

Step	Key (or instruction)	Display/Printout
1. Place an empty container on the balance		+ 65.0 g
2. Tare the balance	TARE	+ 0.0 g
3. Add the first component		+ 120.5 g
4. Store the first component weight. If the print format is set to include data ID codes, the following is printed	F	0.0 g <sub>NET</sub>
5. Add the next component		N1 + 120.5 g
6. Store the 2nd component weight	F	+ 70.5 g
7. Add further components, if desired	As described for steps 5 and 6	0.0 g <sub>NET</sub>
8. Display total weight and fill to desired final weight	CF	+ 191.0 g

# Counting

## Purpose

With the Counting program you can determine the number of parts that each have approximately equal weight.




## Preparation

- Configure the Counting application in the operating menu:  
See "Configuring the Balance"  
Menu code: **2. 1. 4**

- Reference sample quantity:  
Code **3. 3. 1** 5 pcs  
Code **3. 3. 2** 10 pcs (factory setting)  
Code **3. 3. 3** 20 pcs  
Code **3. 3. 4** 50 pcs  
Code **3. 3. 5** 100 pcs
- Storage parameter (display accuracy for counting)  
Code **3. 4. 1** Standard resolution (factory setting)  
Code **3. 4. 2** With 10 times higher resolution than standard  
See also "Configuring the Balance"

## Example

Determine an unknown piece count; weigh the preset reference sample quantity  
Menu: Application program: Counting (menu code **2. 1. 4**);  
Reference sample quantity: 20 pcs (menu code **3. 3. 3**)

Step	Key (or instruction)	Display/Printout
1. Tare the balance	<b>TARE</b>	<b>0.0</b> g
2. Display the reference sample quantity (here: 20 pcs)	<b>F</b> >2 sec.	<b>rEF 20</b> (briefly)
3. Place the reference sample quantity (20 pcs) on the balance (here: 66 g)		<b>+ 66.0</b> g
4. Start the application; if the print format is set to include data ID codes, the following piece weight is printed	<b>F</b>	<b>+ 20</b> pcs
5. Weigh uncounted parts (here: 174 pcs)		<b>wRef + 3.300</b> g <b>+ 174</b> pcs
6. Display weight	<b>F</b>	<b>+ 574.2</b> g
7. Display quantity	<b>F</b>	<b>+ 174</b> pcs
8. Unload the balance		<b>0</b> pcs
9. Delete the reference value	<b>CF</b>	

10. Repeat the procedure starting from step 5, if desired

# Weighing in Percent

## Purpose

This application program allows you to obtain weight readouts in percent which are in proportion to a reference weight.

○ Reference percentage:

- Code 3.3.1 5%
- Code 3.3.2 10% (factory setting)
- Code 3.3.3 20%
- Code 3.3.4 50%
- Code 3.3.5 100%

## Preparation

- Configure the Weighing in Percent application in the operating menu: See "Configuring the Balance."  
Menu code: 2.1.5

○ Storage parameter









- (display accuracy for counting)
  - Code 3.4.1 Standard resolution: With stability (factory setting)
  - Code 3.4.2 With 10 times higher stability than standard
- See also "Configuring the Balance"

## Example

Determine an unknown percentage: store the weight on the balance as a reference percentage

Menu: Application program: Weighing in Percent (menu code 2.1.5)

Menu: Reference percentage: 100% (menu code 3.3.5)

Step	Key (or instruction)	Display/Printout
1. Tare the balance		0.0 g
2. Display the reference percentage:	 >2 sec.	rEF 100
3. Place the reference weight for 100 % on the balance (here: 222.5 g)		+ 222.5 g
4. Start the application; if the print format is set to include data ID codes the following is printed		+ 100.00 % Wxx% + 222.500 g
5. Place an unknown weight on the balance (here: 322.5 g)		+ 144.94 %
6. Display weight		+ 322.5 g
7. Display percentage		+ 144.94 %
8. Unload the balance		0.00 %
9. Delete the reference percentage		
10. Repeat the procedure starting from step 5, if desired.		

# Weigh Averaging

## Purpose

Use this program to determine weights under unstable ambient conditions. In this program, the balance calculates the weight as the average value from a defined number of individual weighing operations, called "subweighs" in this manual.

○ Number of subweighs for weigh averaging:

- 3.3.1 5 subweighs
- 3.3.2 10 subweighs (factory setting)
- 3.3.3 20 subweighs
- 3.3.4 50 subweighs
- 3.3.5 100 subweighs

## Preparation



- Configure the Weigh Averaging application in the operating menu: See "Configuring the Balance."  
Menu code: 2.1.12

See also "Configuring the Balance"

## Example

Determine the weight of a sample in extremely unstable ambient conditions by calculating the average of 10 subweighing operations.

Menu: Application program: Weigh Averaging (menu code 2.1.12)

Step	Key (or instruction)	Display/Printout
1. Tare the balance		0.0 g
2. Display the number of subweighs (here: 10)	 >2 sec.	rEF 10 (briefly)
3. Place sample on the balance (weight readout fluctuates)		8888
4. Start measurement		8888 10 8 8 ...
After 10 subweighs		+ 275.5 g Δ
If the print format is set to include data ID codes, the following is printed		Res + 275.5 g
5. Unload the balance		+ 275.5 g Δ (stable display)
6. Delete the result		
7. Repeat the procedure starting from step 3, if desired.		

# Toggle between Weight Units

With this application program you can toggle the display of a weight value back and forth between two weight units.

Configure the "Toggle Weight Units" application in the operating menu:  
See "Configuring the Balance." Menu code 2.1.2

Menu code	Unit	Conversion factor	Abbr. on printout
1.7.1 3.1.1	Grams <sup>2)</sup>	1	o
1.7.2 o 3.1.2 o	Grams	1	g
1.7.3 3.1.3	Kilograms <sup>1)</sup>	0.001000000000	kg
1.7.4 3.1.4	Carats	5	ct
1.7.5 3.1.5	Pounds <sup>2)</sup>	0.00220462260	lb
1.7.6 3.1.6	Ounces <sup>2)</sup>	0.03527396200	oz
1.7.7 3.1.7	Troy ounces <sup>2)</sup>	0.03215074700	ozt
1.7.8 3.1.8	Hong Kong taels <sup>2)</sup>	0.02671725000	tlh
1.7.9 3.1.9	Singapore taels <sup>2)</sup>	0.02645544638	tls
1.7.10 3.1.10	Taiwanese taels <sup>2)</sup>	0.02666666000	tlt
1.7.11 3.1.11	Grains <sup>2)</sup>	15.43235835000	GN
1.7.12 3.1.12	Pennyweights <sup>2)</sup>	0.64301493100	dwt
1.7.13 3.1.13	Milligrams <sup>2)</sup>	1000	mg
1.7.14 3.1.14	Parts per pound <sup>2)</sup>	1.12876677120	/lb
1.7.15 3.1.15	Chinese taels <sup>2)</sup>	0.02645547175	tlc
1.7.16 3.1.16	Mommes <sup>2)</sup>	0.26670000000	mom
1.7.17 3.1.17	Austrian carats <sup>2)</sup>	5	K
1.7.18 3.1.18	Tola <sup>2)</sup>	0.08573333810	tol
1.7.19 3.1.19	Baht <sup>2)</sup>	0.06578947437	bat
1.7.20 3.1.20	Mesghal <sup>2)</sup>	0.21700000000	MS

o = Factory setting

<sup>1)</sup> = not for models with a readability of  $\leq 0.1$  mg

<sup>2)</sup> = not for model PGB510-2M

## Function

- To toggle the display between the 1st and 2nd weight units:  
Press the  F key

# Configuring the Balance

## Setting the Parameters (Menu Codes)

You can configure your balance to meet individual requirements by selecting from the parameters available in the menu.

Example: Adapt the balance to unstable ambient conditions


Menu code **1. 1. 4**

Step	Key (or instruction)	Display
1. Turn off the balance	[ON/OFF]	
2. Turn the balance back on; while all segments are displayed	[ON/OFF]	<b>1.</b>
○ To navigate within a menu level; the last menu option is followed by the first option	[TARE] repeatedly	<b>2.</b> ⋮ <b>9.</b> <b>1.</b>
3. Select the 2nd menu level	[PRINT]	<b>1. 1.</b>
4. Select the 3rd menu level	[PRINT]	<b>1. 1. 2 o</b>
5. In Menu Level 3: Select the desired option	[TARE] repeatedly	<b>1. 1. 4</b>
6. Confirm new setting; the "o" indicates the currently set option	[PRINT] >2 sec.	<b>1. 1. 4 o</b>
○ Select the next menu level (here: move from the 3rd to the 1st level)	[PRINT]	<b>1.</b>
○ Set other menu codes, if desired	[PRINT], [TARE]	
7. Store parameter settings and exit operating menu or	[TARE] >2 sec.	
○ Exit operating menu without storing changes	[ON/OFF]	
> Restart the application		<b>0.0 g</b>


## Balance Operating Menu (Overview)

o Factory setting

√ User setting

<sup>1)</sup>= not for models with a readability of 0.1 mg

# ISO/GLP-compliant Printout

## Features

You can have the parameters pertaining to the ambient weighing conditions printed before (GLP header) and after (GLP footer) the values of a weighing series. These parameters include:

GLP header:

- Date
- Time at beginning of measurement
- Balance manufacturer
- Balance model
- Balance serial number
- Software version number
- Identification number of the current sampling operation

GLP footer:

- Date
- Time at end of measurement
- Field for operator signature

△ The record can only be output to a KERN data printer 911-017

## Settings

- Set the following menu codes (see "Configuring the Balance"):
  - GLP-compliant record: menu code **5 5 2**
  - ISO/GLP-compliant record after calibration/adjustment only: menu code **7 2 2** or ISO/GLP-compliant record always on: menu code **7 2 3**
  - Line format for printout: With data ID codes – 22 characters: menu code **7 1 2**
  - Printout date/time:
    - 24-h format: menu code **7 3 1**
    - 12-h format: menu code **7 3 2**
- △ No ISO/GLP-compliant record is output if any of the following settings are configured: menu codes **5 1 3**, **5 1 4** (automatic printout) and **7 1 1**

## Function Keys

Press [PRINT] to output header and first measured value.

- > Header is output the first time [PRINT] is pressed

To output header and reference data automatically with an application program active: Press **[F]**

End application program:  
End application program and output GLP footer: Press **[CF]**

The ISO/GLP-compliant record can contain the following lines:

-----	Dotted line
17-Jan-2005 10:15	Date/time (beginning of measurement)
KERN	Balance manufacturer
Mod. TE6100	Balance model
Ser. no. 10105355	Balance serial number
Ver. no. 00-19-41	Software version
ID	ID
-----	Dotted line
L ID	Measurement series no.
wRef + 21.14 g	Counting: Reference weight
Qnt + 235 pcs	Counting result
Qnt + 567 pcs	Counting result
-----	Dotted line
17-Jan-2005 10:20	Date/time (end of measurement)
Name:	Field for operator signature
-----	Blank line
-----	Dotted line

ISO/GLP-compliant printout for external calibration/adjustment


-----	Dotted line
17-Jan-2005 10:30	Date/time (beginning of measurement)
KERN	Balance manufacturer
Mod. TE6100	Balance model
Ser. no. 10105355	Balance serial number
Ver. no. 00-19-41	Software version
ID	ID
-----	Dotted line
Cal. Extern	Calibration/adjustment mode
Set + 5000.0 g	Calibration weight
-----	Dotted line
17-Jan-2005 10:32	Date/time (end of measurement)
Name:	Field for operator signature
-----	Blank line
-----	Dotted line

# Data Interface

## Purpose

Your balance comes equipped with an interface port for connection to a computer or other peripheral device. You can use an on-line computer to change, start and/or monitor the functions of the balance and the application programs.

## Female interface connector


Pin Assignment Chart, 25-pin female interface connector, RS-232:

- Pin 1: Shield
- Pin 2: Data output (TxD)
- Pin 3: Data input (Rx/D)
- Pin 4: Internal ground (GND)
- Pin 5: Clear to send (CTS)
- Pin 6: Not connected
- Pin 7: Internal ground (GND)
- Pin 8: Internal ground (GND)
- Pin 9: Not connected
- Pin 10: Not connected
- Pin 11: Charging voltage for rechargeable battery pack +12 ... +20 V (I<sub>out</sub> 25mA)
- Pin 12: Reset \_ Out \*)
- Pin 13: +5 V output
- Pin 14: Internal ground (GND)
- Pin 15: Universal remote switch
- Pin 16: Not connected
- Pin 17: Not connected
- Pin 18: Not connected
- Pin 19: Not connected
- Pin 20: Data terminal ready (DTR)
- Pin 21: Internal ground (GND)
- Pin 22: Not connected
- Pin 23: Not connected
- Pin 24: Not connected
- Pin 25: +5 V output

## Preparation

You can set these parameters for other devices in the Setup menu (see the chapter entitled "Configuring the Balance"). You will also find a detailed description of the available data interface commands in the file "Data Interface Descriptions for AGB/PGB Models", which you can download from the KERN website, [www.kern-sohn.com](http://www.kern-sohn.com).

The many and versatile properties of these balances can be fully utilized for printing out records of the results when you connect your balance to a KERN data printer. The recording capability for printouts makes it easy for you to work in compliance with ISO/GLP.

\*) = Hardware restart

# Error Codes

Error codes are shown on the main display for 2 seconds. The program then returns automatically to the previous mode (e.g., weighing).

Display	Cause	Solution
No segments appear on the display	No AC power is available The AC adapter is not plugged in Battery or rechargeable battery pack is discharged	Check the AC power supply Plug in the AC adapter Replace the battery or recharge the battery pack using an external charger
H	The load exceeds the balance capacity	Unload the balance
L and E 54	The weighing pan is not in place Something is touching the weighing pan	Place the weighing pan on the balance Move that object that is touching the weighing pan
E 02	Calibration parameter not met, e.g.: – balance not zeroed – balance is loaded	Unload the balance Press <b>TARE</b> to tare the balance Calibrate only when zero is displayed
E 09	When gross value $\leq$ zero; no tare	Tare the balance
E 10	The <b>TARE</b> key is blocked when there is data in the second tare memory (net-total). Only 1 tare function can be used at a time	Press <b>CF</b> to clear the tare memory and release the tare key
E 11	Value input is not allowed for second tare memory	Press <b>TARE</b>
E 22	Weight is too light or there is no sample on the balance	Increase the weight on the balance
E 30	Interface port for printer output is blocked	Contact your local KERN Service Center
Max. weighing capacity is less than indicated under "Specifications"	The balance was turned on without the weighing pan in place	Place the weighing pan on the balance and press <b>[ON/OFF]</b> to turn the balance back on
The weight readout is obviously wrong	The balance has not been calibrated/adjusted The balance was not tared before weighing	Calibrate/adjust the balance Tare the balance

**If any other errors occur, contact your local KERN Service Center!**


# Care and Maintenance

## Service

Regular servicing by a trained KERN dealer will extend the service life of your balance and ensure its continued weighing accuracy. KERN can offer you service contracts, with your choice of regular maintenance intervals. The optimum maintenance interval depends on the operating conditions at the place of installation and on the individual tolerance requirements.

## Repairs

Repair work must be performed by trained service technicians. Any attempt by untrained persons to perform repairs may lead to hazards for the user.

## Cleaning

- Unplug the AC adapter from the wall outlet (mains supply). If you have an interface cable connected to the balance port, unplug it from the port,
- Clean the balance using a piece of cloth which has been wet with a mild detergent (soap)
- After cleaning, wipe down the balance with a soft, dry cloth.

Removing and Cleaning the Weighing Pan:

- Lift up and remove the weighing pan together with the pan support by gripping them from under the shield ring. Make sure that you do not damage the weighing system in doing so.
- ⚠ Make sure that no liquid enters the balance housing.  
⚠ Do not use any aggressive cleaning agents (solvents or similar agents).

## Cleaning Stainless Steel Surfaces

Clean all stainless steel parts regularly. Remove the stainless steel weighing pan and thoroughly clean it separately. Use a damp cloth or sponge to clean any stainless steel parts on the balance. Only use commercially available household cleaning agent that is suitable for use on stainless steel. Clean stainless steel surfaces by wiping them down. Then rinse thoroughly, making sure to remove all residues. Afterwards, allow the balance to dry. If desired, you can apply oil to the cleaned surfaces as additional protection. Solvents are permitted for use only on stainless steel parts.

### **Safety Inspection**

If there is any indication that safe operation of the balance with the AC adapter is no longer warranted:

- Turn off the power and disconnect the equipment from AC power immediately
- > Lock the equipment in a secure place to ensure that it cannot be used for the time being

In any of these cases, notify your KERN service center. Maintenance and repair work may only be performed by service technicians who are authorized by KERN and who

- have access to the required maintenance manuals
- have attended the relevant service training courses

### **Instructions for Recycling**

To ensure adequate protection for safe shipment, your balance has been packaged to the extent necessary using environmentally friendly materials. After successful installation of the balance, you should return this packaging for recycling because it is a valuable source of secondary raw material. For information on recycling options, including recycling of old weighing equipment, contact your municipal waste disposal center or local recycling depot.

# Overview

## Specifications

### AGB Series

Model		AGB120-4	AGB210-4
Weighing capacity	g	120	210
Readability	mg	0.1	0.1
Tare range (subtractive)	g	120	210
Repeatability	<±mg	0.1	0.1
Linearity	<±mg	0.2	0.2
Allowable ambient operating temperature		+10...+30°C (50°F to 86°F)	
Sensitivity drift within +10...+30°C	<±/K	2·10 <sup>-6</sup>	2·10 <sup>-6</sup>
Response time (average)	s	3	3
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels	
Display update (depends on the filter level selected)	s	0.2–0.4	
Test weight with recommended calibration/adjustment tolerance acc. to accuracy class	g	100 (E2)	200 (E2)
	lb	0.2	0.4
Net weight, approx.	kg/lb	3.0/6.6	3.0/6.6
Pan size	mm	80 Ø	80 Ø
	inches	3.3 Ø	3.3 Ø
Weighing chamber height	mm	200	200
	inches	7.9	7.9
Dimensions (WxDxH)	mm	189×251×299	
	inches	7.4×9.9×11.8	
AC power source/power requirements		AC adapter, 230 V or 115 V, +15%...–20%	
Frequency		48–60 Hz	
Power requirements, direct current	V	10 to 20	
Power consumption (average)	W	1	1
Hours of operation with the rechargeable battery pack (see "Accessories")	h	20	20

## PGB Series

Model		PGB150-3	PGB310-3	PGB1500-2	PGB3100-2
Weighing capacity	g	150	310	1,500	3,100
Readability	g	0.001	0.001	0.01	0.01
Tare range (subtractive)	g	150	310	1,500	3,100
Repeatability	±g	0.0015	0.001	0.015	0.01
Linearity	±g	0.003	0.002	0.04	0.02
Allowable ambient operating temperature		+10...+30 °C (50° to 86°F)			
Sensitivity drift within +10...+30 °C	±%/K	3.3·10 <sup>-6</sup>	4·10 <sup>-6</sup>	3.3·10 <sup>-6</sup>	2.5·10 <sup>-6</sup>
Response time (average)	s	3	3	3	3
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels			
Display update (depends on the filter level selected)	s	0.2–0.8			
Test weight included		100 g	200 g	1000 g	2000 g
Details are given in the chapter "Calibration/Adjustment"					
Calibration/adjustment tolerance		(F2)	(F1)	(F2)	(F1)
Test weight with recommended calibration/adjustment tolerance acc. to accuracy class	g lb	100 (F1) 0.2	200 (E2) 0.4	1000 (F1) –	2000 (E2) 4
Net weight, approx.	kg	1.7	2.2	1.7	2.2
Pan size	mm	100 Ø	100 Ø	174x143	174x143
Dimensions (B×T×H)	mm	189x251x120		189x251x70	
AC power source/power requirements		AC adapter 230 V or 115 V, +15% to -20%			
Frequency		48–60 Hz			
AC power source/direct current	V	10–20			
Power consumption (average)	W	0.75			
Hours of operation with the rechargeable battery pack (see "Accessories")	h	25	20	25	20

## PGB Series

Model		PGB610-2	PGB6100-1
Weighing capacity	g	610	6100
Readability	g	0.01	0.1
Tare range (subtractive)	g	610	6100
Repeatability	<±g	0.01	0.1
Linearity	<±g	0.02	0.2
Allowable ambient operating temperature		+10...+30°C (50°F to 86°F)	
Sensitivity drift within +10...+30°C	<±/K	5 · 10 <sup>-6</sup>	
Response time (average)	s	2	2
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels	
Display update (depends on the filter level selected)	s	0.2–0.8	
Test weight with recommended calibration/adjustment tolerance acc. to accuracy class	g lb	500 (F2) 1	5000 (F2) 10
Net weight, approx.	kg/lb	1.1/2.4	1.4/3.1
Pan size	mm inches	116 Ø 4.6 Ø	174×143 6.9×5.6
Dimensions (W×D×H)	mm	189×251×70 (7.4×9.9×2.8 inches)	
AC power source/power requirements		AC adapter, 230 V or 115 V, +15%...– 20%	
Frequency		48–60 Hz	
Power requirements, direct current	V	10 to 20	
Power consumption (average)	W	0.75	
Approx. hours of operation with: – 9-volt battery, approx.	h	11	
– fully charged rechargeable battery	h	2.5	
Hours of operation with the rechargeable battery pack (see "Accessories")	h	25	

Model	PGB510-2M	
Type designation	DT BH 210	
Accuracy class*	Ⓔ	
Maximum weighing capacity, Max* g	510	
Scale interval d*	g	0.01
Verification scale interval e*	g	0.1
Minimum capacity, Min*	g	0.5
Application range acc. to CD*	g	0.5 – 510
Tare range (subtractive)	<100% of maximum weighing capacity	
Operating temperature range	+10 ... +30 °C (50°F to 86°F)	
Response time (average)	s	2
Adaptation to application and ambient conditions	By selection of 1 of 4 optimized filter levels	
Display update interval (depends on filter level selected)	s	0.1 – 0.4
Pan/load plate size	mm	116 Ø
Dimensions (WxDxH)	mm	189x251x70
Net weight, approx	kg	1.4
AC power source/ power requirements	AC adapter, 230 V or 115 V, +15%...-20%	
Frequency	48 – 60 Hz	
AC power/constant voltage	V	10 to 20
Power consumption (average)	W	0.75
Hours of operation using the rechargeable battery pack (see "Accessories")	h	25

\* CD = Council Directive 90/384/EEC on non-automatic weighing instruments used in the EU (including the Signatories of the Agreement on the European Economic Area)

# Declarations of Conformity

The CE marking affixed to the equipment indicates that the equipment meets the requirements of the following Directive(s) issued by the Council of the European Union:

**Council Directive 89/336/EU  
“Electromagnetic compatibility (EMC)”**

**Important Note:**

The operator shall be responsible for any modifications to KERN equipment and for any connections of cables or equipment not supplied by KERN and must check and, if necessary, correct these modifications and connections. On request, KERN will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

**73/23/EU “Electrical equipment designed for use within certain voltage limits”**

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installations in your country. For information on the current legal requirements of your country, please contact your local KERN customer service office.

**Additional Directive for Weighing Instruments Used in Legal Metrology: Directive 90/384/EU “Non-automatic weighing instruments”**


This Directive also regulates the performance of the EC verification by the manufacturer, provided that an EC Type-Approval Certificate has been issued and the manufacturer has been accredited by an officer of a Notified Body registered at the Commission of the European Community for performing such verification. On February 15, 1993, the manufacturers received accreditation for performing verification.

**Subsequent Verifications within the European Countries**

The validity of the verification will become void in accordance with the national regulations of the country in which the balance is used. For information on verification and legal regulations currently applicable in your country, and to obtain the names of the persons to contact, please contact your local KERN office, dealer or service center.

**Konformitätserklärung zur Richtlinie 90/384/EWG**  
**Declaration of Type Conformity to Directive No. 90/384/EEC**  
**Déclaration de Conformite au Type selon la Directive 90/384/CEE**  
**Declaración de Conformidad de tipo según Directiva N° 90/384/CEE**  
**Dichiarazione di Conformità alla Direttiva 90/384/CEE**

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E-Mail: [info@kern-sohn.de](mailto:info@kern-sohn.de)

Internet: [www.kern-sohn.de](http://www.kern-sohn.de)

Copyright by Kern, Balingen, Germany.

All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Kern.

The status of the information, specifications and illustrations in this manual is indicated by the date given below. Kern reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Status: December 2005, Kern, Balingen, Germany

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## **KERN & Sohn GmbH**