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Ziegelei 1
D-72336 Balingen
E-Mail: info@kern-sohn.com

Tel: +49-[0]7433-9933-0
Fax: +49-[0]7433-9933-149
Internet: www.kern-sohn.com

**Betriebsanleitung Temperatur-Kalibrierset
Operating instruction temperature calibrating set
Mode d'emploi kit de calibrage de température**


KERN DLB_A01

Version 1.1
10/2011

Français English Deutsch

DLB_A01-BA-def-1111

KERN DLB_A01

Version 1.1 10/2011

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Temperatur-Kalibrierset

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
2.1 Overview of display

Display	Description
MAX	Means the maximum value of a measurement
HOLD	Means a recorded display value
!	Measuring range exceeded
	Capacity of battery exhausted

2.2 Keyboard overview

		Description
1		Switch unit on and off
2		Press button, the display background will be illuminated for approx. 5 s.
3		Not documented
4		Select unit, °C or °F or K
5	0.1°	Switch-over readability 1°C ⇔ 0.1°C / 1°F ⇔ 0.1°F.
	1°C	
	MAX	Press button, the highest measured value will be displayed, identified by the indicator „MAX“
	HOLD	Press the button to record the current displayed value, identified by the indicator „HOLD“. It remains displayed until it is deleted by pressing button again.

3 Change battery

If the battery voltage is not more sufficient, the battery symbol  will appear in the display.

- ☞ Switch off appliance and remove cable from temperature sensor plate.
- ☞ Open the lid of the battery compartment loosening the screw.
- ☞ Replace battery.

Ensure correct polarisation.


- ☞ Close again the lid.

4 Adjust/calibrate the temperature value of the KERN DLB-A

4.1 Preparation

We recommend sometimes to check the temperature value of the device using the optional temperature calibrating set DLB-A01. Before you do this, allow the device to cool down for at least 3 hours after the last heating phase.

- ⇒ Remove the separate component parts „sample dish“ according to their sequence

⇒ Install the temperature-calibration set acc. to fig.

- ⇒ Switch-on the heating top on the rear side.
- ⇒ Close cover of the heating top
- ⇒ Switch-on the temperature calibration set and carry out the following adjustments, see chap. 2.2. Select unit °C, readability 1°C

Invoke service function:

- ⇒ In weighing mode press the **MENU** button, the menu item „Therm“ will be displayed.

ThErM

- ⇒ Confirm by pressing the **PRINT** button. Press repeatedly the **MENU** button until „Service“ will be displayed.

SERvICE

- ⇒ Keep pressed the **PRINT** button for **2 sec.** Select the desired settings by pressing the **MENU** button

tNP tEST

Temperature calibration

tNP AdJ

Temperature adjustment


HEAt SEL

Settings of lamp

4.2 Calibrate temperature


During temperature calibration only a check is carried out, i.e. no values are changed.

- ⇒ Invoke service function „Tmp test“, see chap. 4.1.

The image shows a rectangular LCD display with a black border. Inside the display, the text 'TMP TEST' is shown in a white, monospaced, all-caps font.

- ⇒ Confirm by pressing the **PRINT** button. The current setting will be displayed. Use the navigation buttons **↓** **↑** to select the temperature during the check, available 35- 160°C.


The image shows a rectangular LCD display with a black border. Inside the display, the text '80 C' is shown in a white, monospaced, all-caps font.

- ⇒ Confirm using the **PRINT** button, the appliance heats up to the set temperature. After approx. 15 minutes compare the temperature displayed on the thermometer with that of the moisture analyzer. If the two values do not match, we recommend a temperature adjustment, see chap. 3.
- ⇒ Use the **ON/OFF** button to finish the calibration.
- ⇒ Use the **ON/OFF** button to return into weighing mode.

4.3 Adjust temperature


The temperature is measured at two selectable points and it is possible to correct it there.

- ⇒ Invoke service function „Tmp Adj“, see chap. 4.1.


The LCD display shows the text "TMP ADJ" in a monospaced font.

- ⇒ Confirm by pressing the **PRINT** button. The current setting for the first temperature point is displayed. Or leave it on factory setting 80°C or select a desired temperature with the help of the navigation buttons ↓ ↑ (available 50-130°C).


The LCD display shows the text "P1 80 C" in a monospaced font.

- ⇒ Confirm with the **PRINT** button, the first heat-up phase is started. The remaining time is displayed.


The LCD display shows the text "15 min" in a monospaced font.

- ⇒ Temperature calibration of first point will take 15 min. Compare the temperature displayed on the thermometer with that of the moisture analyzer. If the two values do not match, correct them using the navigation buttons ↓ ↑ and confirm with the **PRINT** button


The LCD display shows the text "t1 80 C" in a monospaced font.

- ⇒ The current setting for the second temperature point is displayed. Or leave it on factory setting 150°C or select the desired temperature with the help of the navigation buttons ↓ ↑. The second temperature point has to be at least 30°C over the first one, max. 160°C.


The LCD display shows the text "P2 150 C" in a monospaced font.

- ⇒ Confirm with the **PRINT** button, the second heat-up phase is started. The remaining time is displayed.

The LCD display shows the text "15 min" in a monospaced font.

- ⇒ Temperature calibration of the second point will take 15 min. Compare the temperature displayed on the thermometer with that of the moisture analyzer. If the two values do not match, correct with the help of the navigation buttons ↓ ↑.

The LCD display shows the text "t2 150 C" in a monospaced font.

- ⇒ Confirm by pressing the **PRINT** button. The adjustment is completed, the appliance returns into the menu. Use the **ON/OFF** button to return into weighing mode.

Due to the AUTO-POWER OFF function the display unit of the temperature calibration set switches off automatically during adjustment. Turn-on appliance again, the current temperature adjustment will not be disturbed.

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