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KERN MH
Version 1.0 12/99
Operating Instructions
Electronic Handy Balances

1 Technical Data

1.1 Balance

<i>KERN</i>	<i>MH5K5</i>	<i>MH10K10</i>
<i>Readout</i>	<i>5 g</i>	<i>10 g</i>
<i>Weighing range</i>	<i>5 kg</i>	<i>10 kg</i>

1.2 Quartz Clock

Indicating hours, minutes, a.m. and p.m.

1.3 Thermometer

To measure ambient temperature Range from – 10°C to 40°C.

2 Parts of the Balance

There are three function buttons to operate the three functions.

3 Operation

3.1 First time use

Push the lid of the battery housing in the direction indicated. Insert two #7(AAA) batteries in the polarity shown. The power source is then connected. After a brief self –test, lasting for 1 to 2 seconds, the operation is in clock mode and starts registering time. Close the lid.

3.2 Set Time

Switch to Clock Mode:

Set Minute

Press the **MODE** key once. The minute figure blinks. Press the **Set** key to adjust minutes.

Set Hour

Press the **MODE** key once again. The hour figure blinks. Press the **Set** key to adjust hours.

Press the **MODE** key once again to return to normal clock mode. If the **MODE** key is not pressed after setting, the display will return to normal clock mode automatically after a few seconds.

3.3 Scale Mode

Under the clock mode, press the **Scale** key once. The display will switch to scale mode.

Weighing

Press the **TARE** key once. The display registers zero. Hook the object being weighed. The steady reading is the weight of the object.

Net weight

Hook the container. Press the **TARE** key to zero the display. Put the object inside the container. The reading is the net weight of the object.

Memory

The product has the feature to retain reading. When the reading (over 0.05 kg) reaches a steady state the memory is activated. the symbol “j” for seeing the reading easily is displayed and the buzzer sounds. Now the load could be released. This is very helpful in weighing heavy objects. When the weight varies the reading also varies. When the load is released the reading indicates the last steady weight. After seven seconds the display disappears.


3.4 Measuring Temperature

When under the scale mode press **Temp** button once the ambient temperature is displayed. Pressing the button again returns the operation to scale mode. Please note. the thermometer works only under steady state. When the temperature varies it takes time to reach a steady state before accurately displaying the ambient temperature.

4 Important information

These electronic balances are precision instruments. Avoid disruptive environmental conditions such as drafts, vibrations, and electrical and magnetic fields. If considerable electromagnetic fields occur, display deviations are possible. In that case change the place. Equally avoid rapid temperature changes. Avoid high humidity, vapour, and dust, since the balance is not hermetically sealed. Do not allow the balance to come into direct contact with liquids, which could penetrate the measuring system. Therefore clean only with dry or damp material. Do not use solvents because paint or plastic parts could be damaged as a result. Place the item to be weighed carefully. Do not hang a constant load. Shocks and severe impacts must be avoided. Remove any spilled weighing substances immediately. Never operate the balance in areas where there is a danger of explosion. The standard unit is not explosion-proof. Check the balance regularly with external check weights. Opening the balance voids the guarantee. When not in use for a long time, remove the batteries to avoid damages from possible battery leakage.

Always make sure that there are no living beings or materials below the load that could be injured or damaged!

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