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Ziegelei 1

D-72336 Balingen

E-Mail: info@kern-sohn.com

Tel: +49-[0]7433-9933-0

Fax: +49-[0]7433-9933-149

Internet: www.kern-sohn.com

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The diagram illustrates the navigation through the P2.2.6 menu to assign functions to F1 and F2. The sequence of screens and actions is as follows:

- Screen 1:** TARA 0000 G. Action: MENU button.
- Screen 2:** P1 JUSTAGE. Action: F1 (up arrow) button.
- Screen 3:** P2 BETRIEBSMODUL. Action: F2 (right arrow) button.
- Screen 4:** P2.1 ZUGANG. Action: F1 (up arrow) button.
- Screen 5:** P2.2 WÄGUNG. Action: F2 (right arrow) button.
- Screen 6:** P2.2.1 ABLESUNG. Action: F1 (up arrow) button, 5 x.
- Screen 7:** P2.2.6 ABKÜRZUNGEN. Action: F2 (right arrow) button.
- Screen 8:** P2.2.6 F1 KEIN. Action: F2 (right arrow) button.

Two callouts show the 'KEIN' selection process:

- Callout 1:** From the 'P2.2.6' screen, pressing F1 (up arrow) leads to 'KEIN', and pressing F3 (down arrow) leads to 'TARA EINGEBEN'. Pressing F1 again leads to 'EIN-LOGGEN'.
- Callout 2:** From the 'P2.2.6' screen, pressing F1 (up arrow) leads to 'KEIN', and pressing F3 (down arrow) leads to 'TARA EINGEBEN'. Pressing F1 again leads to 'EIN-LOGGEN'.

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The diagram illustrates the sequence of menu screens and key presses to activate a weighing unit. It starts with the scale in 'TARA 0000 G' mode. Pressing the MENU key leads to the 'P8.1.1' menu. Pressing F3 down leads to 'G' and 'NEIN'. Pressing F2 right leads to 'P8.1.2' and 'MG' and 'NEIN'. Pressing F2 right leads to 'P8.1.2' and 'JA'. Pressing F2 right leads to 'P8.1' and 'G' and 'JA'. Pressing F2 right leads to 'P8.1.2' and 'G' and 'JA'. Navigation arrows and a PRINT key are also shown.

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Analysis and precision balances

Contents

1	Technical data	6
2	Declaration of Conformity	11
3	Device overview	12
3.1	Display and Keyboard Summary	15
3.1.1	Navigation keys / numeric input	16
4	Basic Information (General)	17
4.1	Utilisation in accordance with specification	17
4.2	Improper Use	17
4.3	Warranty	17
4.4	Monitoring of Test Resources	18
5	Basic Safety Precautions	18
5.1	Pay attention to the instructions in the Operation Manual	18
5.2	Personnel training	18
6	Transport and storage	18
6.1	Testing upon acceptance	18
6.2	Packaging / return transport	18
7	Unpacking, Setup and Commissioning	21
7.1	Installation Site, Location of Use	21
7.2	Unpacking and checking	21
7.3	Placing	22
7.4	Levelling	25
7.5	Mains connection	26
7.6	Commissioning	26
7.6.1	Connecting to power supply	26
7.6.2	Switching on and off	27
7.6.3	Selecting user language	27
7.6.4	Log on and off (activated user profile only)	27
7.7	Peripheral equipment	30
8	Menu	31
8.1	Menu structure	31
8.2	Navigation in the menu	33
9	Weighing mode	34
9.1	Weighing	34
9.1.1	Simple weighing	34
9.1.2	Zeroing	34
9.1.3	Weighing with tare	35
9.2	Underfloor weighing	38
9.3	Settings <P2.2 WEIGHING>	39
9.3.1	<P2.2.1.1> Filter Matching weighing type to ambient conditions	41
9.3.2	<P2.2.1.2> Rest position control display	42
9.3.3	<P2.2.1.3> Auto zero function	43
9.3.4	<P2.2.1.4> Hiding last digit of display	44
9.3.5	<P2.2.1.5> Environment	45
9.3.6	<P2.2.2 > AUTOTARE / <P2.2.3.2> AUTO THRESHOLD	46
9.3.7	<P2.2.3.> Manual / automatic data output / <P2.2.3.2> AUTO THRESHOLD	47

9.3.8	<P2.2.4> Informational line	49
9.3.9	<P2.2.5> Special information	50
9.3.10	<P2.2.6> Assigning function keys	51
10	< P 8 UNTIS >	52
10.1	Change weighing unit < P8.1>	52
10.2	Defining start-up weighing unit < P8.2>	53
10.3	Defining freely programmable weighing units < P8.3 / P8.4>	54
11	General parameter < P 6 OTHER>	55
12	Show device info < P 7 INFO>	57
13	Adjustment < P1 JUSTAGE >	58
13.1	Automatic adjustment using internal weight (models AEJ, PLJ)	59
13.2	Settings in menu <P1 JUSTAGE>	61
13.2.1	<P1.1 INT. JUST.>	62
13.2.2	<P1.2 EXT. JUST.> (AES models)	63
13.2.3	< P1.3 USER ADJUSTMENT > Adjustment by using user-defined external weight	65
13.2.4	< P1.4 JUST.TEST > Adjustment test	67
14	Define logs < P5 PRINTOUTS >	68
14.1	Adjustment log <P5.1>	68
14.2	Weighing log	70
14.2.1	Define header<P5.2>	70
14.2.2	Print weighing results / define GLP log<P5.3>	71
14.2.3	Define footer<P5.4>	72
14.3	User-defined data output < P5.5> - <P5.8 SPECIAL PRINTOUT >	73
14.4	Universal variables	75
14.5	Save logs to USB storage medium / print	76
15	Data bases	77
15.1	User database	78
15.1.1	Create new user	79
15.1.2	Delete user	80
15.1.3	For invoking a user	81
15.2	Product database	81
15.2.1	Create new product	82
15.3	Container weights	84
15.3.1	Create new container weight	84
15.4	Weighing results	85
15.5	Alibi memory	87
15.6	Export / import databases to USB storage medium	89
16	Application settings	91
16.1	Selecting an application	91
17	Parts counting	93
17.1	Settings <P2.3> PARTS COUNTING>	93
17.2	Carry out parts counting	94
17.2.1	Determination of the reference piece weight by weighing	94
17.2.2	Enter known reference single weight as numeric value	97
17.2.3	Invoking reference single weight from the product database	98
18	Checkweighing	99
18.1	Settings <P2.4> CHECKWEIGHING>	99
18.1.1	Entering limits with a function key	100
18.1.2	Invoke limits from the product database	102
19	Dosing	103
19.1	Settings <P2.5> DOSING	103
19.1.1	Enter target weight with the function key	104
19.1.2	Invoke target weight from product database	106
20	Percent weighing	107

20.1	Settings <P2.6> DEVIATION>	107
20.2	Perform percentage weighing	108
20.2.1	Determination of the reference weight by weighing	108
20.2.2	Enter known reference weight as numeric value	109
20.2.3	Invoke the reference weight from the product database.....	110
21	Animal weighing	111
21.1	Animal weighing by manual start	111
21.2	Animal weighing with automatic start.....	113
22	Determining the density of solid matter and fluids	116
22.1	Settings <P2.7> or <P2.8> DENSITY DETERMINATION>	117
23	Statistics	118
24	Totalization	122
25	Peak value function	126
26	Pipette calibration	128
26.1	Settings <P2.13 PIPETTE CALIBRATION >.....	128
26.2	Procedure.....	128
27	Communication / setup menu < P3 >	135
27.1	Configuring a <P3.1 / P3.2> RS-232 interface.....	136
27.2	<P3.3> WLAN interface (factory option KERN PLJ-A06).....	138
27.2.1	P3.3.2 - SELECT NET	139
27.2.2	<P3.3.3> - NETWORK PARAMETER	141
27.3	USB interfaces	143
27.3.1	USB type A.....	143
27.3.2	USB type B (PC connection).....	143
28	Devices / setup menu < P4 >	145
28.1	Continuous data transmission for remote control commands < P4.1 Computer >	145
28.2	To start data transmission, press the PRINT-key < P4.2 Printer >	149
28.3	< P4.3 Barcode Scanner >	152
29	Communication with external devices (printer / PC)	153
29.1	Pin assignment of RS 232 interface (frontal view):.....	153
29.2	Interface cable	154
29.3	Data transmission format	155
29.4	Interface commands.....	157
29.5	General reply format	158
29.6	Detailed information to communication logs	159
29.6.1	Zeroing	159
29.6.2	Taring	159
29.6.3	Query tare value.....	159
29.6.4	Set tare value	160
29.6.5	Send stable weighing value in standard weighing unit	160
29.6.6	Send weighing value immediately in standard weighing unit.....	161
29.6.7	Send stable weighing value in current weighing value	162
29.6.8	Send weighing value immediately in current weighing unit	163
29.6.9	Start continuous output in standard weighing unit	163
29.6.10	Stop continuous output in standard weighing unit	163
29.6.11	Start continuous output in current weighing unit.....	164
29.6.12	Stop continuous output in current weighing unit	164
29.6.13	Enter lower tolerance limit.....	164
29.6.14	Enter upper tolerance limit	164
29.6.15	Query lower tolerance unit	165
29.6.16	Query upper tolerance unit.....	165
29.6.17	Enter reference single weight	165
29.6.18	Enter target weight (dosing)	165
29.6.19	Enter reference weight (percentage weighing)	166
29.6.20	Confirm / simulation "Press PRINT-key"	166
29.6.21	Start internal adjustment	166

29.6.22	Disable automatic internal adjustment	166
29.6.23	Enable automatic internal adjustment.....	167
29.6.24	Enter serial number	167
29.6.25	Disable weighing balance keyboard	167
29.6.26	Enable weighing balance keyboard	167
29.6.27	List of available operating modes	168
29.6.28	Set operating mode	168
29.6.29	Query current operating mode	168
29.6.30	Confirm sound signal	169
29.6.31	Send all messages implemented	169
30	 Servicing, maintenance, disposal	170
31	 Troubleshooting / error messages	171

1 Technical data

KERN	AEJ 100-4CM	AEJ 200-4CM	AEJ 200-5CM
Readability (d)	0.1 mg	0.1 mg	0.01 mg; 0.1 mg
Weighing range (max)	160 g	220 g	82 g; 220 g
Minimum load (Min)	10 mg	10 mg	1 mg
Verification value (e)	1 mg	1 mg	1 mg; 1 mg
Verification class	I	I	I
Reproducibility	0.2 mg	0.2 mg	0.04 mg; 0.1 mg
Linearity	± 0.3 mg	± 0.3 mg	± 0.1 mg; 0.2 mg
Stabilization time	4 sec.	4 sec.	6 sec; 3.5sec.
Adjustment weight	internal	internal	internal
Warm-up time	8 h	8 h	8 h; 8 h
Weighing Units	g, mg, ct	g, mg, ct	g, mg, ct
Smallest part weight for piece counting	0.1 mg	0.1 mg	0.1 mg
Reference quantities at piece counting	10, 20, 50, freely selectable	10, 20, 50, freely selectable	10, 20, 50, freely selectable
Weighing plate, stainless steel [cm]	Ø 10	Ø 10	Ø 8.5
Net weight (kg)	5.4	5.4	5.4
Permissible ambient condition	+10° C to +40° C	+10° C to +40° C	+10° C to +40° C
Humidity of air	80 % relative (not condensing)	80 % relative (not condensing)	80 % relative (not condensing)
Input Voltage	110 - 230 V AC,	110 - 230 V AC,	110 - 230 V AC,
Power pack secondary voltage	13.5 V – 16 V	13.5 V – 16 V	13.5 V – 16 V
Interface	RS 232C; USB	RS 232C; USB	RS 232C; USB
Underfloor weighing device	Hooks	Hooks	Hooks


KERN	AES 100-4C	AES 200-4C
Readability (d)	0.1 mg	0.1 mg
Weighing range (max)	160 g	220 g
Reproducibility	0.2 mg	0.2 mg
Linearity	± 0.2 mg	± 0.2 mg
Stabilization time	4 sec.	4 sec.
Recommended adjusting weight not supplied (class)	150 g (E2)	200 g (E2)
Warm-up time	8 h	8 h
Weighing Units	ct, dwt, g, gn, mg, oz	ct, dwt, g, gn, mg, oz
Smallest part weight for piece counting	0.5 mg	0.5 mg
Reference quantities at piece counting	10, 20, 50, freely selectable	10, 20, 50, freely selectable
Weighing plate, stainless steel [cm]	Ø 10	Ø 10
Net weight (kg)	5.4	5.4
Permissible ambient condition	+18° C to +30° C	+18° C to +30° C
Humidity of air	80 % relative (not condensing)	80 % relative (not condensing)
Input Voltage	110 - 230 V AC,	110 - 230 V AC,
Power pack secondary voltage	13.8 V	13.8 V
Interface	RS 232C; USB	RS 232C; USB
Underfloor weighing device	Hooks	Hooks

KERN	PLJ 300-3CM	PLJ 600-3CM
Readability (d)	1 mg	1 mg
Weighing range (max)	360 g	600 g
Minimum load (Min)	20 mg	20 mg
Verification value (e)	10 mg	10 mg
Verification class	II	II
Reproducibility	1 mg	2 mg
Linearity	± 4 mg	± 4 mg
Stabilization time	4 sec.	4 sec.
Adjustment weight	internal	internal
Warm-up time	2 h	4 h
Weighing Units	ct, g, kg	ct, g, kg
Smallest part weight for piece counting	5 mg	5 mg
Reference quantities at piece counting	10, 20, 50, freely selectable	10, 20, 50, freely selectable
Weighing plate, stainless steel [mm]	128 x 128	128 x 128
Net weight (kg)	4.6	4
Permissible ambient condition	+10° C to +40° C	+10° C to +40° C
Humidity of air	80 % relative (not condensing)	80 % relative (not condensing)
Input Voltage	100 - 240 V	100 - 240 V
Power pack secondary voltage	13.5 V – 16 V	13.5 V – 16 V
Interface	RS 232C; USB	RS 232C; USB
Underfloor weighing device	Hooks	Hooks


KERN	PLJ 700-3CM	PLJ 3000-2CM
Readability (d)	1 mg	10 mg
Weighing range (max)	750 g	3500 g
Minimum load (Min)	20 mg	500 mg
Verification value (e)	10 mg	100 mg
Verification class	II	II
Reproducibility	2 mg	10 mg
Linearity	± 4 mg	± 40 mg
Stabilization time	4 sec.	4 sec.
Adjustment weight	internal	
Warm-up time	4 h	2 h
Weighing Units	ct, g, kg	
Smallest part weight for piece counting	5 mg	
Reference quantities at piece counting	10, 20, 50, freely selectable	
Weighing plate, stainless steel [mm]	128 x 128	195 x 195
Net weight (kg)	3.5	4.4
Permissible ambient condition	+10° C to +40° C	
Humidity of air	80 % relative (not condensing)	
Input Voltage	100 - 240 V	
Power pack secondary voltage	13.5 V – 16 V	
Interface	RS 232C; USB	
Underfloor weighing device	Hooks	

Dimensions

AEJ-C/AES-C:

PLJ-C:

English

2 Declaration of Conformity

To view the current EC/EU Declaration of Conformity go to:

www.kern-sohn.com/ce

- i** The scope of delivery for verified weighing balances (= conformity-rated weighing balances) includes a Declaration of Conformity.

3 Device overview

KERN AEJ/AES-C:


Pos.	Designation
1	Glass wind screen
2	Weighing pan
3	Display
4	Keyboard
5	Levelling screw

KERN PLJ-C:

Pos.	Designation
1	Glass wind screen
2	Weighing pan
3	Display
4	Keyboard
5	Levelling screw

Connections (all models):


Pos.	Designation
1	Mains connection
2	COM 2 port
3	COM 1 port
4	USB 2 port (for devices such as printer, PC)
5	USB 1 port (for devices such as keyboard)

English

3.1 Display and Keyboard Summary


Example AEJ-C:


Pos. Designation

- 1 Stability display
- 2 Indicator “zero display”
- 3 Information line
- 4 The non-verified value is given in brackets in verified scales.
- 5 Weighing unit “g”

Button	Function
	<ul style="list-style-type: none"> Select operating mode
	<ul style="list-style-type: none"> Start adjustment
	<ul style="list-style-type: none"> Invoke setup menu
	<ul style="list-style-type: none"> Switch machine on/off
	<ul style="list-style-type: none"> Invoke database menu

	<ul style="list-style-type: none"> Function key / shortcut to basic functions
	<ul style="list-style-type: none"> Taring Add data record
	<ul style="list-style-type: none"> Switch-over weighing unit Exit menu
	<ul style="list-style-type: none"> Shortcut key for most frequently occurring functions and settings (For configuration see chap. 9.3.10) Navigation keys, see chap. 2.1.1
	<ul style="list-style-type: none"> Data output to external device (weighing mode) Confirm / store settings
	<ul style="list-style-type: none"> Zeroing Data record deletion

3.1.1 Navigation keys / numeric input

Button	Designation	Function
	Navigation button ↑	Increase flashing digit Move through numbers and alphabet in a forward sequence
	Navigation button ↓	Decrease flashing digit Move through numbers and alphabet in a backward sequence
	Navigation button →	Digit selection to the right
	Navigation button ←	Digit selection to the left
	Navigation button ↵	Confirm entry
		Cancel input

The connection of a PC keyboard via USB interface makes numerical entry easier.

4 Basic Information (General)

4.1 Utilisation in accordance with specification

The balance you purchased is intended to determine the weighing value of material to be weighed. It is intended to be used as a “non-automatic balance”, i.e. the material to be weighed is manually and carefully placed in the centre of the weighing pan. As soon as a stable weighing value is reached the weighing value can be read.

4.2 Improper Use

Do not use balance for dynamic add-on weighing procedures, if small amounts of goods to be weighed are removed or added. The “stability compensation“ installed in the balance may result in displaying an incorrect measuring value! (Example: Slowly draining fluids from a container on the balance.)

Do not leave permanent load on the weighing pan. This may damage the measuring system.

--

Impacts and overloading exceeding the stated maximum load (max) of the balance, minus a possibly existing tare load, must be strictly avoided. Balance may be damaged by this.

Never operate balance in explosive environment. The serial version is not explosion protected.

The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.

The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.

4.3 Warranty

Warranty claims shall be voided in case

- Our conditions in the operation manual are ignored
- The appliance is used outside the described uses
- The appliance is modified or opened
- Mechanical damage or damage by media, liquids, natural wear and tear
- The appliance is improperly set up or incorrectly electrically connected
- The measuring system is overloaded

4.4 Monitoring of Test Resources

In the framework of quality assurance the measuring-related properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (www.kern-sohn.com) with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited DKD calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

5 Basic Safety Precautions

5.1 Pay attention to the instructions in the Operation Manual

- Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.

5.2 Personnel training

The weighing balance may only be operated and maintained by trained staff.


6 Transport and storage

6.1 Testing upon acceptance


When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

6.2 Packaging / return transport

- Keep all parts of the original packaging for a possibly required return.
- Only use original packaging for returning.





- Prior to dispatch disconnect all cables and remove loose/mobile parts.

See example shown for KERN AEJ/AES-C

- Reattach possibly supplied transport securing devices.
- Secure all parts such as the glass wind screen, the weighing platform, power unit etc. against shifting and damage.

7 Unpacking, Setup and Commissioning

7.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

Therefore, observe the following for the installation site:

- Place the balance on a firm, level surface;
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight;
- Protect the balance against direct draughts due to open windows and doors;
- Avoid jarring during weighing;
- Protect the balance against high humidity, vapours and dust;
- Do not expose the device to extreme dampness for longer periods of time. Non-permitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for ca. 2 hours at room temperature.
- Avoid static charge of goods to be weighed or weighing container.

If electro-magnetic fields or static charge occur, or if the power supply is unstable major deviations on the display (incorrect weighing results) are possible. In that case, the location must be changed.

7.2 Unpacking and checking

Remove device and accessories carefully from packaging, remove packaging material and place device at the planned work place. Verify that there has been no damage and that all packing items are present.

Scope of delivery / serial accessories:

- Balance, see chap. 2
- Mains adapter
- Operating manual
- Protective cover
- Flush-mounted hook

7.3 Placing

- i** The right place is decisive for the accuracy of the weighing results of high-resolution analytical and precision balances (see chap. 7.1).

Model AEJ-CM


Turn the transport locking device in the direction of the arrow <OPEN> and remove.

Keep transport locking device in a safe place for possible future transportation.


⇒ Installation of weighing plate

Models AEJ/AES-C (except AEJ-5CM)


- Open glass draft shields on all sides
- Place protection sheet (1) into the weighing chamber
- Place centring ring (2) on the protection sheet
- Put weighing plate (3) on the centring ring
- Place wind protection ring (4) into the weighing chamber

Model AEJ-5CM


- Open glass draft shields on all sides
- Remove transport securing (1).
- Place protection sheet (2) into the weighing chamber
- Place centring ring (3) on the protection sheet
- Put weighing plate (4) on the centring ring
- Place wind protection ring (5) into the weighing chamber

**Models PLJ-C:
10 mg models**

- Remove adhesive tape from the rubber element (1)
- Place weighing pan (2) on the rubber elements

1 mg models



- Remove adhesive tape from the rubber element (1)
- Place weighing pan (2) on the rubber elements
- Place wind shield (3)

7.4 Levelling

Exact alignment and stable installation are an essential requirement for repeatable results. You can compensate for minor unevenness or tilts of the footprint by levelling the weighing balance.

- ⇒ Level balance with foot screws until the air bubble of the water balance is in the prescribed circle.

- ⇒ Check levelling regularly

7.5 Mains connection

Select a country-specific power plug and insert it in the power unit.

Check, whether the voltage acceptance on the scales is set correctly. Do not connect the scales to the power grid unless the information on the scales (sticker) matches the local mains voltage.

Only use KERN original mains adapter. Using other makes requires consent by KERN.

Important:

- Before starting your weighing balance, check the mains cable for damage.
- Ensure that the power unit does not come into contact with liquids.
- Ensure access to power plug at all times.

In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap. 1). During this warming up time the balance must be connected to the power supply (mains, accumulator or battery).

The accuracy of the balance depends on the local acceleration of gravity. Strictly observe hints in chapter Adjustment.

7.6 Commissioning

7.6.1 Connecting to power supply

Connect the weighing balance by plugging the mains cable into the power supply. The weighing balance will perform a self-test, briefly showing software version.

In models with integrated adjustment weight, the noise of the charging system's motor may be audible. An internal adjustment will be taking place.

Wait until the basic display appears:


7.6.2 Switching on and off

The weighing balance will remain switched on as long as it is connected to the power supply. Operate the ON/OFF-key to turn on/off the display.

To turn on the display, press the **ON/OFF**-key.
Wait for the basic display to appear.

To **switch off**, press the **ON/OFF** button again. During this process the weighing balance will be moving into ready mode.

The weighing balance will be ready for operation immediately after activation without requiring any time to warm up.

7.6.3 Selecting user language

The display is configured to German upon delivery. For setting additional languages see chap.11.

7.6.4 Log on and off (activated user profile only)

The weighing balance is delivered without user profile and you can configure your own settings. If the settings are to be protected against unauthorised access, create a password protected user profile.

The user database (<operator / user> database allows you to create up to ten user profiles with specific settings and user rights, see chap.15.1.

The weighing balance is also designed for use without user management.


An enabled user profile will be indicated by the  icon on the display.

Login:


Several keys are available for logging in: You can use either, the **F**-key, **DATA**-key or one of the shortcut keys (provided it is assigned to <LOG IN>).

⇒ Use the navigation keys to select the desired user such as ADMIN.





⇒ If necessary use the navigation keys to enter the password and confirm by pressing the PRINT-key; for numeric entry see chap. 3.1.1.




-key


⇒ Use the navigation keys to select <ADMIN> and confirm by pressing the PRINT-key.

⇒ If necessary use the navigation keys to enter the password and confirm by pressing the PRINT-key, see chap. 3.1.1 Numeric entry.

Logout:

7.7 Peripheral equipment

Before connecting or disconnecting of additional devices (printer, PC) to the data interface, always disconnect the balance from the power supply.

With your balance, only use accessories and printers by KERN, as they are ideally tuned to your balance.

8 Menu

8.1 Menu structure

The menu is subdivided into seven menu blocks (P1 – P7).

P1	Adjustment P1.1 Internal adjustment P1.2 External adjustment P1.3 User adjustment P1.4 Adjustment test	Available settings see chap. 13.2.1. s. Kap. 13.2.2 only AES models s. Kap. 13.2.3 only AES models see chap. 13.2.4.
P2	Operation mode P2.1 Access P2.2 Weighing P2.3 Parts counting P2.4 Checkweighing P2.5 Dosing P2.6 Percent weighing P2.7 Density determination „solids“ P2.8 Density determination „liquids“ P2.9 Animal weighing P2.10 Statistics P2.11 Adding up P2.12 Peak value function P2.13 Pipette calibration	see chap. 16.1. see chap. 9. see chap. 17. see chap. 18. see chap. 19. see chap. 20. see chap. 22. see chap. 21. see chap. 23. see chap. 24. see chap. 25. see chap. 26.
P3	Communication P3.1 COM 1 P3.2 COM 2 P3.3 WIFI	see chap. 27. see chap. 27.1. see chap. 27.2.
P4	Equipment P4.1 Computer P4.2 Printer P4.3 Barcode reader P4.4 Additional display P4.5 External keys	see chap. 28. see chap. 28.1. see chap. 28.2. see chap. 28.3. see chap. 28.3. Not documented

P5	Printouts	
	P5.1 Adjustment log	see chap. 14.1.
	P5.2 Header	see chap. 14.2.1.
	P5.3 GLP log	see chap. 14.2.2.
	P5.4 Footer	see chap. 14.2.3.
	P5.5 Special printout 1	see chap.
	P5.6 Special printout 2	
	P5.7 Special printout 3	see chap. 14.3.
	P5.8 Special printout 4	
	P5.10 Variable 1	see chap. 14.4.
	P5.11 Variable 2	
P6	General parameters	see chap. 11.
	P6.1 User language	
	P6.2 Access right level for menu access	
	P6.3 Key sound	
	P6.4 Display backlight	
	P6.5 Sleep mode	
	P6.6 Auto off	
	P6.7 Setting date	
	P6.8 Setting time	
	P6.9 Date format	
	P6.10 Time format	
P7	Info	see chap. 12.
	P7.1 Balance identification no.	
	P7.2 Model designation	
	P7.3 Software version status	
	P7.4 Temperature environment	
	P7.5 Printout device settings	
P8	Units	see chap. 10.
	P8.1 Access.	
	P8.2 Standard weighing unit	
	P8.3 Free weighing unit 1	
	P8.4 Free weighing unit 2	

8.2 Navigation in the menu

Call up menu	Press  , and the first menu block 1 CAL will be displayed.
Select menu block	Use the navigation buttons to select the individual menu blocks one by one. Use the navigation key  to scroll down. Use the navigation key  to scroll up.
Select menu item	Confirm selected menu block by pressing  . The first menu item of the selected menu block will be shown. Use the navigation buttons to select the individual menu items one by one. Use the navigation key  to scroll down. Use the navigation key  to scroll up.
Select setting	Confirm selected menu item with  and the current setting will be shown.
Change settings	Use the navigation keys to switch over into the available settings. Use the navigation key  to scroll down. Use the navigation key  to scroll up.
Confirm setting / back to menu	Either confirm by pressing  or cancel by pressing  .
Exit menu	Press 
Return to weighing mode	Press  repeatedly.

9 Weighing mode

The sequence of a simple weighing process is described in the chapter on Weighing below. Apart from the work steps described there (simple weighing, reset to zero, taring) the weighing balance offers additional options for adapting the “Weighing” application to your requirements. For the settings available see chap. 9.3

9.1 Weighing

 For warm-up time required for stabilisation, see chpt 1.

9.1.1 Simple weighing

- Check zero display [] and set to zero with the help of the  key, as required.
- Place load on sample dish
- Wait until the stability display appears ().
- Read weighing result.

After pressing  you can save the weighing value and if an optional printer is connected, print out the result.

9.1.2 Zeroing

In order to obtain optimal weighing results, reset to zero the balance before weighing. Reset to zero is restricted to a range of $\pm 2\%$ Max. An error message Err2 will appear in case of values exceeding $\pm 2\%$ Max.

- Unload the balance
- Press  and zero display as well as indicator  will appear.

9.1.3 Weighing with tare

➤ Taring

The dead weight of any weighing container may be tared away by pressing a button, so that the following weighing procedures show the net weight of the goods to be weighed.

- Place the weighing container on the weighing pan and close the draft shields, as required.
- Wait until the stability display appears (☐), then press . Zero display and indicator (**Net**) appear. The weight of the container is now internally saved.
- Weigh the material.
- Wait until the stability display appears (☐).
- Read net weight.

- When the balance is unloaded the saved taring value is displayed with negative sign.
- Taring negative values is not permitted. Any attempt to tare negative values will result in error message “Err 3”
- The taring process can be repeated any number of times. The limit is reached when the whole weighing range is exhausted.
- You can assign the tare weight to a product in the database. Selecting the product will automatically load the tare weight, see chap. 15.2.1.

➤ Automatic tare function

When the automatic tare function is enabled <P2.2.2 AUTOTARA YES > the weight on the pan will first be saved as tare, see chap. 9.3.6

➤ **Numerical input of tare (PRE-TARE)**

Condition:

Shortcut keys F1, F2, F3 or F4 will be assigned to the <ENTER TARE> function, see chap. 9.3.10.

- Press the shortcut key and the display for entering the tare value will appear.

- Use the navigation keys to enter known tare weight (See chap. 3.1.1) and

confirm by pressing .

- The display will return to weighing mode. The entered weight will be saved as tare weight, the indicator (**net**) and the tare weight showing the minus-sign appear.

- Place the filled weighing container.
- Wait until the stability display appears (▢).
- Read net weight.

➤ **Load container weight from the database.**

- Press shortcut key (F1, F2, F3, F4) assigned to the <SELECT TARE> function (See chap. 9.3.10.).

or

- After pressing , select <SELECT TARE> function.

The display for selecting a tare memory will appear.

- Use the navigation keys to select the desired data record and confirm by pressing .


- The display will return to weighing mode, the indicator (net) and the selected container weight with minus-sign will appear.

- For creating a container weight in the database see chap. 15.3.1

➤ **Delete tare**

- Remove the load from the weighing pan and press . The (Net) indicator turns off, the zero display shows.

9.2 Underfloor weighing

Objects not fitting on the weighing pan due to their weight or shape can be weighed with the help of the underfloor weighing option.

Proceed as follows:

- Switch off the balance.
- Open closing cover at the balance bottom.
- Suspend hook for underfloor weighing **carefully and completely**.
- Do not turn lug.
- Place weighing balance over an opening.
- Attach load to hook, set to zero and carry out weighing procedure.


Fig.1: Setup of balance for underfloor weighing

- Always ensure that all suspended objects are stable enough to hold the desired goods to be weighed safely (danger of breaking).
- Never suspend loads that exceed the stated maximum load (max) (danger of breaking)
- Always ensure that there are **no persons, animals or objects** that might be damaged underneath the load.

After completing the underfloor weighing the opening on the bottom of the balance must always be closed (dust protection).

9.3 Settings <P2.2 WEIGHING>

i Navigation in menu see chap. 8.2

Menu structure:


Menu item	Parameters	Available settings / explanation
P2.2.1 READING	P2.2.1.1 Filter	EXTRA FAST Weighing balance reacts very sensitively and very fast, very stable place of installation.
		FAST Weighing balance reacts sensitively and fast, stable place of installation.
		MEDIUM Standard – normal ambient conditions
		SLOW Weighing balance reacts insensitively but slowly, unstable place of installation.
		VERY SLOW Weighing balance reacts very insensitively, but very slowly, very unstable place of installation.
	P2.2.1.2 Weighing speed	PRECISE Rest position control precise, unstable place of installation
		FAST Rest position control fast, very stable place of installation
		FAST + PRECISE Rest position fast + precise, stable place of installation
	P2.2.1.3 Autozero	Selectable yes / no
	P2.2.1.4 Last digit	ALWAYS Last decimal place will be shown
		NEVER Last decimal place will be hidden
		WHEN STABLE Last decimal place will only be shown for stable weighing values
	P2.2.1.5 Environment	STABLE
		UNSTABLE

P2.2.2 AUTOTARE	Selectable yes / no	
P2.2.3 PRINT MOD. (specifications PRINT-key)	P2.2.3.1 MODE	ALL Weighing log will be generated immediately after the PRINT-key has been pressed, regardless whether the weighing value is stable or not
		WHEN STABLE Output after the PRINT-key has been pressed will not take place until the weighing value has stabilised
		AUTO Automatic data output after exceeding set minimum weight <P2.1.3.2 Auto threshold>
	P2.2.3.2 AUTO THRESHOLD for automatic tare and automatic operation	Input minimum weight for automatic output: Weighing value issued automatically if current weighing value exceeds entered minimum value. The next weighing value will not be issued unless the weighing value has meanwhile dropped below the entered weighing value.
P2.2.4 INFORM.	Define info line - options: TARE, NET, GROSS, USER, PRODUCT, STNFO, NONE, BARGRAPH	
P2.2.5 NST INF	Define special information (arbitrary text, max. 19 characters)	
P2.2.6 HOT KEYS	F1	Options: NONE / LOG ON / ENTER TARE / SELECT TARE / PRINT HEADER / PRINT FOOTER / VARIABLE 1, 2 /m change mode
	F2	
	F3	
	F3	

i When working with user management, please remember that all settings for the weighing application will be saved to the enabled user profile. Each user can select his/her individual settings for this application. For that reason ensure that the correct user profile was selected before making any changes to settings.


9.3.1 <P2.2.1.1> Filter Matching weighing type to ambient conditions

Please note that data processing geared to higher stability will always bring about a deceleration of reaction times while an acceleration of reaction times will have a negative impact on stability, that is, the higher the filter level, the longer the reaction time.

9.3.2 <P2.2.1.2> Rest position control display

This setting is used to define how quickly the weighing balance considers a measured value as stable for release.


9.3.3 <P2.2.1.3> Auto zero function

This function is used to automatically tare minor weight fluctuations (such as minor soiling on weighing pan). By enabling this function you automatically ensure that each weighing process starts with zero display.

In the event that small quantities are removed or added to the material to be weighed, incorrect weighing results can be displayed due to the “stability compensation”. (e.g. slow flow of liquids from a container placed on the balance, evaporating processes).

When apportioning involves small variations of weight, it is advisable to switch off this function.

9.3.4 <P2.2.1.4> Hiding last digit of display

Readability may be reduced by 1 digit on the instrument, as required. The last decimal place will be rounded and removed from the display.


			4 x
			3 x

English

9.3.5 <P2.2.1.5> Environment

This function allows you to optimally adapt the weighing balance to ambient conditions. In a stable environment (for instance free of draft and vibrations) select <STABLE> as a setting.

In an instable environment exposed to constantly changing conditions select <INSTABLE>.


9.3.6 <P2.2.2 > AUTOTARE / <P2.2.3.2> AUTO THRESHOLD

The automatic tare function automatically saves the first stable weighing value as tare value.

The parameter <P2.2.3.2> AUTO THRESHOLD is used to define which minimum weight must be placed before automatic taring takes place.

Further automatic taring will only take place when the display drops below the set value for the auto threshold.

**9.3.7 <P2.2.3.> Manual / automatic data output /
<P2.2.3.2> AUTO THRESHOLD**

This function is used whether and how data output is to take place.


Available settings:

Menu item	Available settings / explanation	
<P2.2.3.1> MODE	<IF STABLE>	Manual output of stable weighing values after pressing PRINT key
	<ANY>	Manual output for stable and instable weighing values after pressing PRINT key
	<AUTO>	Sequence of operations: <ul style="list-style-type: none"> • Set balance to zero • Place goods to be weighed on balance. The first stable weighing value will automatically be issued when the current weighing value exceeds the set threshold value. • A new output requires that the weighing value first falls short and then exceeds the set threshold.
<P2.2.3.2> AUTO THRESHOLD	Setting for threshold in [g], from which automatic data output is to start.	

0.0000g
TARE 0000 G

P2.2.3.1
EACH

P1
CALIBRATION

P2.2.3.1
WHEN STABIL

P2
WORKING MODE

P2.2.3.1
MODES WHEN STABIL

P2.1
ACCESSIBILITY


P2.2.3.2
AUTO THRES. 0000

P2.2
WEIGHING

P2.2.3.2
0000000

P2.2.1
READING

P2.2.3.2
000 10000

P2.2.3
PRINT MOD

P2.2.3.2
AUTO THRES. 10000

4 x


P2.2.3.1
MODE EACH

0.0000g
TARE 0000 G

9.3.8 <P2.2.4> Informational line


This function is used to define which additional information is to be shown in the line of text at the lower edge of the display.

9.3.9 <P2.2.5> Special information

This function is used to define user-specific information.

When this function <P2.2.5 SINFO> is enabled, it will be shown in the line of text at the lower edge of the display.


9.3.10 <P2.2.6> Assigning function keys

The function keys F1, F2, F3 and F4 allow direct access to frequently required functions and settings. The available settings depend on the operating mode selected. The function assigned to a menu item of the two keys will be executed when the respective key is tapped.

Selectable functions in weighing mode:

Adjustment	Funktion, welche die Taste ausführen soll
<NONE>	No function
<SELECT PRODUCT>	Select product from product data base
<LOG ON>	Select user
<ENTER TARE>	Enter tare weight
<SELECT TARE>	Select tare weight from product database
<PRINT HEADER>	Print created header
<PRINT FOOTER>	Print created footer
<VARIABLE 1>	Select / edit variable 1
<VARIABLE 2>	Select / edit variable 2

10 < P 8 UNTIS >

This function is used to define which weighing units are to be used in the weighing balance. Available settings:

Menu item	Available settings / explanation
P8.1 ACCESS	This function is used to define the units into which the weighing result is to be converted upon pressing the UNIT-key.
P8.2 START WEIGHING UNIT	This function is used to define the unit to be shown on starting up the weighing balance.
P8.3 / P8.4 B. UNIT 1 / 2	These functions are used to define two independent weighing units.

10.1 Change weighing unit < P8.1>

Enabling units to which conversion is to take place:

- Use the navigation buttons (↓ ↑) to select more units and to enable/disable them as described before.

- Press back to weighing mode  several times.


How to change weighing units

- To return to the enabled weighing units go to weighing mode by .

- The weighing units available depend on the model and calibration of the weighing balance.
- When switching-on the balance, the unit in which the balance has been switched off, will be displayed.


10.2 Defining start-up weighing unit < P8.2 >

- The weighing units available depend on the model and calibration of the weighing balance.
- Disconnect scale shortly from power supply. The selected start-up weighing unit will be shown upon reconnection.

10.3 Defining freely programmable weighing units < P8.3 / P8.4>

i Freely programmable weighing units are only available for weighing balances without type approval.

- Use the navigation keys to enter the conversion factor and confirm by pressing the PRINT-key.
- Press  several times to return to weighing mode .

11 General parameter < P 6 OTHER>

This allows you to configure parameters that influence the operation of the weighing balance such as user language, date / time display, key sounds, display brightness etc.

Menu structure:

Navigation in menu see chap. 8.2

Menu item	Available settings / explanation
P6.1 LANGUAGE	User language selectable: Deutsch/Espanol/Francais/Turk/Cesky/Italiano/Magyar/Polски/English
P6.2 ACCESS	Access levels selectable: ADMIN / USER / MASTER USER
P6.3 KEY SOUND	Options: YES/NO
P6.4 BACKLIGHT	Backlighting of display, brightness selectable: NONE / 10 ⇆ 100 (dark ⇆ bright)
P6.5 STANDBY MODE	When sleep mode is enabled, the weighing balance will turn off backlighting after a set length of time provided that during this time the weighing pan was not loaded nor a key was pressed. As soon as any key is pressed or there is a change of weight the display and backlighting will once again be turned on. Options: none, 0.5, 1, 2, 3, 5 minutes
P6.6 AUTO. TURN OFF	When this function is enabled the display will automatically be switched off and the time displayed, after a defined length of time of inactivity has passed. During this process the weighing balance will be moving into ready mode. The weighing balance will be ready for operation immediately after activation without requiring any time to warm up. To restart the weighing balance, press  . Options: none, 1, 2, 3, 5, 10 minutes
P6.7 DATE	Setting date
P6.8 TIME	Time setting



<p>P6.9 DATE FORMAT (on logs)</p>	<p>Options: YYYY.MM.DD / YYYY.DD.MM / DD.MM.YYYY / MM. DD.YYYY</p>
<p>P6.10 TIME FORMAT (on logs)</p>	<p>Options: 12H / 24H The selected format AM or PM will be printed on the log next to the time.</p>
<p>P6.11 GLP AUTO TEST</p>	<p>Repeatability test activated by placing the internal weight Sequence of operations:</p> <ul style="list-style-type: none"> ➤ To start the test, confirm the function <P6.9 AUTOTEST GLP> by pressing . ➤ The internal weight will be applied twice, followed by 10 repeat measurements. ➤ Adjustment ➤ Calculation of standard deviation ➤ Automatic log printout <p>Printout example KERN YKB-01N</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <pre> ----- GLP Autotest Rep.----- Scale Type PL Balance ID 150854 User Admin Soft. Ver. r3.2.6 Date 05.01.2016 Time 11:34:17 ----- Measur. No. 10 Division 0.01 g Int. Weight mass 765.35 Filter Means Result cong. Fast+rel. ----- Max. Dev. 0.01 g Repeatab. 0.003 g Signature ----- </pre> </div>


i




- Additional log output by pressing .
- You can cancel the test by pressing .
- Press back to weighing mode  several times.

12 Show device info < P 7 INFO >

This function is used to show and print information about the weighing balance.

- Use  to scroll down the list of weighing balance information.

- To print the current weighing balance settings confirm parameter <P7.5 PRINT SETTINGS> by pressing 
- Press back to weighing mode  several times.

13 Adjustment < P1 JUSTAGE >

As the acceleration value due to gravity is not the same at every location on earth, each balance must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the balance has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the balance periodically in weighing operation.

Possible adjustments include:

- Automatic adjustment using internal weight (models AEJ, PLJ) time and temperature controlled, see chap. 13.1.
- Manual internal adjustment (models AEJ, PLJ)
To activate, press the CAL-key, see chap. 13.2.1.
- Adjustment using external weight (model AES), see chap. 13.2.2.

- Observe stable environmental conditions. A warming up time (see chapter 1) is required for stabilization.
- Internal adjustment can only be carried out when the weighing pan is free of objects. If the weighing pan is loaded, the error message <Range exceeded> will appear.
- You can cancel the adjustment by pressing the **TARE**-key.
- Carry out external adjustment as near as possible to the balance's maximum weight (recommended adjustment weight see chap. 1). Weights of different nominal values (>30% max) may be used for adjustment but are not optimal for technical measuring. Info about test weights can be found on the Internet at: <http://www.kern-sohn.com>.


13.1 Automatic adjustment using internal weight (models AEJ, PLJ)

With the internal adjustment weight, the weighing accuracy can be checked and re-adjusted at any time.

The automatic adjustment function is always enabled. You can start adjustment at any time by pressing the **CAL**-key manually.

Temperature / time controlled adjustment is taking place:

Two minutes prior to automatic adjustment „CAL“ in combination with a clock and thermometer icon will be shown at the upper edge of the display.

The user must complete his/her weighing process within this time. After two minutes [**CAL 30**] appears on the display.

This starts a “count down” of 30 seconds [**CAL 30**] → [**CAL 0**]. You may cancel adjustment during these 30 seconds by pressing the **TARE**-key. This makes the balance return to weighing mode in order to e. g. complete an unfinished measurement.

After approx. further 5 minutes the automatic adjustment will restart and [**CAL 30**] appears.

You may interrupt the adjustment process several times. However, keep in mind that an interruption of the adjustment may lead to incorrect weighing results.

The adjustment log will be printed automatically if an optional printer is connected.

For defining an adjustment log see chap. 14.1 / <P5.1 Adjustment log>.

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Automatic adjustment is carried out:

- When balance was disconnected from the mains.
- When the temperature is changing
The moment a change in temperature of 3°C is detected by the temperature sensor, the weighing scale will carry out a fully automatic adjustment. The message „CAL“ in combination with a thermometer icon indicates that an adjustment is about to take place.

- at the end of a defined time interval
After a time interval set in the menu **<P1.5 AU.TIME.JUST>** has lapsed (selectable 1 – 12 h) the weighing balance will announce that adjustment is about to take place by sending message „CAL“ in combination with a clock icon:

In weighing scales with type approval certificate

- Menu item **<P1.5 AU.TIME.JUST>** disabled
- Adjustment takes place every 3 hours


13.2 Settings in menu <P1 JUSTAGE>

i Navigation in menu see chap. 8.2

Menu item	Explanation
P1.1 INT JUST	Internal adjustment (AEJ, PLJ models only)
P1.2 EXT JUST	External adjustment (AES models only)
P1.3 USER ADJUSTMENT	External adjustment with user-defined adjustment weight (AES models only)
P1.4 ADJUSTMENT TEST	Adjustment test (AEJ, PLJ models only)
P1.4 AU.THERM.JUST	Define start criteria for automatic adjustment (disabled for weighing balances with type approval)
P1.5 AU.TIME.JUST	Define time interval for automatic adjustment (disabled for weighing balances with type approval)

13.2.1 <P1.1 INT. JUST.>

- ⇒ For automatic internal adjustment see chap. 13.1
- ⇒ Starting internal adjustment manually



- Ensure that there are no objects on the weighing pan.
- After successful adjustment the balance automatically returns to weighing mode.

13.2.2 <P1.2 EXT. JUST.> (AES models)

- Carry out adjustment with the help of the recommended adjustment weight (see chapter 1 “Technical Specifications”).
- Observe stable environmental conditions. A warming up time (see chapter 1) is required for stabilization.


- Invoke menu item < P1.2 EXT. JUST. >

- Ensure that there are no objects on the sample dish. Press the **PRINT**-key and wait until the weighing value for the required adjustment weight is shown.

- Place the required adjustment weight carefully in the centre of the sample dish and press the **PRINT** key. Wait until the “MASS REMOVAL” is shown.

- Take away adjustment weight. If the adjustment was successful, this will be followed by automatic return to the menu.

P 1.2
EXT CALIB

0.000g
TARE 0000 G


- The adjustment log will automatically be printed if an optional printer is connected.
For defining an adjustment log see chap. 14.1 / <P5.1 Adjustment log>.
- Press the UNIT-key to return to weighing mode

Printout example KERN YKB-01N:

```
-----Cal. Report-----  
Calib. type           External  
User  
Project  
Date                 04.05.2015  
Time                 9:49:09  
Balance ID           461677  
Cal. differ.         0.001 g  
-----  
Signature  
.....
```

13.2.3 < P1.3 USER ADJUSTMENT > Adjustment by using user-defined external weight


- Invoke menu item < P1.3 USER JUSTAGE >


- Use the navigation keys to enter the nominal value for the desired adjustment weight (>30% max) (See chap.3.1.1) and confirm by pressing the PRINT-key.


- Ensure that there are no objects on the sample dish. Press the **PRINT**-key and wait until the weighing value for the required adjustment weight is shown.

- Place the required adjustment weight carefully in the centre of the sample dish and press the **PRINT** key. Wait until the “MASS REMOVAL” is shown.


- Take away adjustment weight. If the adjustment was successful, this will be followed by automatic return to the menu.

- The adjustment log will be printed automatically if an optional printer is connected. For defining an adjustment log see chap. 14.1 / <P5.1 Adjustment log>.
- Press the UNIT-key to return to weighing mode

13.2.4 < P1.4 JUST.TEST > Adjustment test

This function is applied to compare the internal adjustment weight with the saved value. The test is carried out automatically. The result is shown on the display and, if an optional printer is connected, a printout will be issued via the RS 232 interface.

14 Define logs < P5 PRINTOUTS >

14.1 Adjustment log <P5.1>

This menu is used to define which information is to be printed in the adjustment logs.

Settings default log <P5.1 ADJUSTMENT LOG>:

- ⇒ Navigation in menu see chap. 8.2
- ⇒ All information confirmed by <yes> will be printed.

Menu item	Available settings / explanation
<P5.1.1>DRAFT	Enter name of project, max. 16 characters Numerical input, see chap. 3.1.1.
<P5.1.2> ADJUSTMENT TYPE	Type of adjustment executed, yes / no
<P5.1.3> USER	Print name of logged on user, yes / no
<P5.1.4> DRAFT	Print name of draft yes / no
<P5.1.5> DATE	Print date of adjustment, yes / no
<P5.1.6> TIME	Print time of adjustment, yes / no
<P5.1.7> ID WEIGHING BALANCE	Print serial number, yes / no
<P5.1.8> DIFFERENCE ADJUSTMENT	Print deviation from last adjustment, yes / no
<P5.1.9> LINES	Print lines to separate signature, yes / no
<P5.1.10> SIGNATURE	Insert line block for signing log, yes / no

Printout example 1 (KERN YKB-01N):

Menu item	Adjustment	Protocol
P5.1.2	yes	<pre> -----Cal. Report----- Calib. type External User ADMIN Project PROJEKT1 Date 05.05.2015 Time 10:22:11 Balance ID 461677 Cal. differ. 0.000 g ----- Signature </pre>
P5.1.3	yes	
P5.1.4	yes	
P5.1.5	yes	
P5.1.6	yes	
P5.1.7	yes	
P5.1.8	yes	
P5.1.9	yes	
P5.1.10	yes	

Printout example 2 (KERN YKB-01N):

Menu item	Adjustment	Protocol
-----------	------------	----------

P5.1.2	no
P5.1.3	yes
P5.1.4	no
P5.1.5	yes
P5.1.6	no
P5.1.7	yes
P5.1.8	yes
P5.1.9	no
P5.1.10	no

```
-----Cal. Report-----
User                ADMIN
Date                05.05.2015
Balance ID          461677
Cal. differ.        0.000 g
Signature

.....
```

14.2 Weighing log

This menu item is used to define which information is to be printed for default weighing logs.

- ⇒ Navigation in menu see chap. 8.2
- ⇒ The menu is subdivided into three submenus where you can define the settings for the header, the logging of weighing results and the footer.
- ⇒ All information confirmed by **<yes>** will be printed.
- ⇒ These settings only apply to weighing mode.

14.2.1 Define header<P5.2>

This submenu is used to define which information is to be printed on the header of weighing logs. The header will be printed on pressing the shortcut keys **F1** or **F2** (provided they are assigned to <Print header>) or on printing the **F**-key and selecting <Print header>.

Menu item	Available settings / explanation
<P5.2.1> LINES	Insert broken lines, yes / no
<P5.2.2> OPERATING MODE	Print name of application, yes / no
<P5.2.3> DATE	Print current date, yes / no
<P5.2.4> TIME	Print current time, yes / no
<P5.2.5> WEIGHING BALANCE TYPE	Print weighing balance type, yes / no
<P5.2.6> ID WEIGHING BALANCE	Print serial number, yes / no
<P5.2.7> USER	Print name of logged on user, yes / no
<P5.2.8> PRODUCT	Print name of selected product, yes / no
<P5.2.9> VARIABLE 1	Print variable 1, yes / no
<P5.2.10> VARIABLE 2	Print variable 2, yes / no For defining variable see chap. 14.4
<P5.2.11> EMPTY LINE	Print empty line, yes / no
<P5.2.12>ADJUSTMENT LOG	PRINT log of last adjustment, yes / no. For settings of adjustment log see chap. 14.1
<P5.2.13> SPECIAL PRINTOUT	Place user-defined printout in header. Options: NO / SPECIAL PRINTOUT 1 / SPECIAL PRINTOUT 2 / SPECIAL PRINTOUT 3 / SPECIAL PRINTOUT 4 / For creating special printouts see chap. 14.3

Printout example (KERN YKB-01N):

Menu item	Adjustment	Protocol
-----------	------------	----------

P5.2.1	yes
P5.2.2	yes
P5.2.3	yes
P5.2.4	yes
P5.2.5	no
P5.2.6	yes
P5.2.7	yes
P5.2.8	no
P5.2.9	no
P5.2.10	no
P5.2.11	yes
P5.2.12	no
P5.2.13	none

```

-----
Work. mode           Weighing
Date                07.05.2015
Time                9:14:52
Balance ID         461677
User                ADMIN
  
```


14.2.2 Print weighing results / define GLP log<P5.3>

Quality assurance systems require printouts of weighing results as well as of correct adjustment of the balance stating date and time and balance identification.

To print the GLP log, press the PRINT-key.

Menu item	Available settings / explanation
<P5.3.1> DATE	Print current date, yes / no
<P5.3.2> TIME	Print current time, yes / no
<P5.3.3> USER	Print name of logged on user, yes / no
<P5.3.4> PRODUCT	Print name of selected product, yes / no
<P5.3.5> VARIABLE 1	Print variable 1, yes / no For defining variable see chap. 14.4
<P5.3.6> VARIABLE 2	Print variable 2, yes / no For defining variable see chap. 14.4
<P5.3.7> NET	Print net weight value of current weighing
<P5.3.8> TARE	Print tare value of current weighing
<P5.3.9> GROSS	Print gross weight value of current weighing
<P5.3.10> CURRENT RESULT	Print weighing value in current unit
<P5.3.11> ADJUSTMENT LOG	PRINT log of last adjustment, yes / no. For settings of adjustment log see chap. 14.1
<P5.3.12> SPECIAL PRINTOUT	Options: NO / SPECIAL PRINTOUT 1 / SPECIAL PRINTOUT 2 / SPECIAL PRINTOUT 3 / SPECIAL PRINTOUT 4 / For creating a special printout see chap. 14.3

Printout example (KERN YKB-01N):

Menu item	Adjustment	Protocol
P5.3.1	yes	
P5.3.2	yes	
P5.3.3	no	
P5.3.4	no	
P5.3.5	no	
P5.3.6	no	
P5.3.7	yes	
P5.3.8	yes	
P5.3.9	yes	
P5.3.10	none	

14.2.3 Define footer<P5.4>

This submenu is used to define which information is to be printed on the footer of weighing logs. The footer will be printed on pressing the shortcut keys **F1** or **F2** (provided they are assigned to <Print footer>) or on pressing the **F**-key and selecting <Print footer>.

Menu item	Available settings / explanation
<P5.4.1> OPERATING MODE	Print name of application, yes / no
<P5.4.2> DATE	Print current date, yes / no
<P5.4.3> TIME	Print current time, yes / no
<P5.4.4> BALANCE TYPE	Print weighing balance type, yes / no
<P5.4.5> ID BALANCE	Print serial number, yes / no
<P5.4.6> USER	Print name of logged on user, yes / no
<P5.4.7> PRODUCT	Print name of selected product, yes / no
<P5.4.8> VARIABLE 1	Print variable 1, yes / no; For defining variable see chap. 14.4
<P5.4.9> VARIABLE 2	Print variable 2, yes / no For defining variable see chap. 14.4
<P5.4.10> LINES	Insert broken lines, yes / no
<P5.4.11> EMPTY LINE	Print empty line, yes / no
<P5.4.12> ADJUSTMENT LOG:	PRINT log of last adjustment, yes / no. For settings of adjustment log see chap. 14.1
<P5.4.12> SIGNATURE	Insert line block for signing log, yes / no
<P5.4.13> SPECIAL PRINTOUT	Place user-defined printout in footer. Options: No / special printout 1 / special printout 2 / special printout 3 / special printout 4 / For creating special printouts see chap. 14.3

14.3 User-defined data output < P5.5> - <P5.8 SPECIAL PRINTOUT >

Menu <P5> provides options that allow you to define up to four different special printouts of maximal 160 characters each (<P5.5> - <P5.8>).

In addition to the option for entering running text (For numeric entry see chap.) 3.1.1) the menu also provides the following variables:

General variables:

%%	Printout of 1 character “%” (i.e. in order to print out a character %, two %% must be entered)
%V;	Net weight in current unit
%N	Net weight in default unit
%G	Gross weight in current unit
%T	Gross weight in default unit
%D	Current date
%M	Current time
%I	Balance ID-no.
%R	Program no.
%P	Project number.
%U	User no.
%F	Current working mode
%C	Date and time of last adjustment
%K	Type of last adjustment
%S	Currently selected product
%Y	Deviation from last adjustment
%1	Variable 1; for defining variable see chap. 14.4
%2	Variable 2; for defining variable see chap. 14.4

Function based variables:

	Description	Function
%W	Reference quantity	Parts counting
%H	Upper limit	Tolerance weighing
%L	Lower limit	
%A	Target weight	Dosing
%B	Reference weight	Percent weighing

Print related variables:

\\	Printout of 1 character “/” i.e. in order to print out 1 character /, two // must be entered)
\C	CRLF (carriage return line feed) line beginning next line
\R	CR (carriage return) line beginning
\N	LF (line feed) next line
\T	tabulator
\F	Paper feed (for PCL printers)
\0	End character string
%E	Cut paper for EPSON printer

These variables must be entered in UPPER CASE.

The connection of a PC keyboard via USB interface makes the entry easier.

Example of text input in weighing mode:

KERN*SOHN

DATE: <current date,>


TIME: <current time>

MASS: <current net weight in default unit>

*****SIGNATURE:

Input:

- Invoke desired menu item P5.5 – P5.8 (special printout 1 -4) and enter text; for numeric entry see chap. 3.1.1

- Confirm entry by pressing PRINT key and the weighing balance will return to the menu. Press UNITS key to return to weighing mode.

14.4 Universal variables

Menu items P5.9 and P5.10 are used to define the content (such as batch no.) of variable1 and variable2 (max. 32 characters).


The principle of entry is identical to that of special printouts, see chap. 14.3.

14.5 Save logs to USB storage medium / print

This option allows you to save weighing results as text file to a USB storage medium and to print these files on a printer connected to your PC.


Procedure:

- ⇒ You will see a prompt asking you to connect a USB storage medium to the USB PORT, <IEI IMPORT / EXPORT>. Press the UNITS-key to return to work mode.
- ⇒ Go to menu <P4.2.1 Devices / Printer / Port> to select the setting <USB-STICK >, see chap.28.2.
- ⇒ From then onwards weighing results will be automatically exported as text file (printout.txt.) to the USB storage medium the moment the PRINT-key is pressed.
- ⇒ Press the ON/OFF-key to save your data as file. Only then remove the USB storage medium and connect to the PC.

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- ⇒ Open printout.txt., then edit or print on a printer connected to the PC, as required.

```
printout.txt - Editor
Datei Bearbeiten Format Ansicht ?
Datum 22.02.2016
Uhrzeit 13:07:14
Bediener
Produkt test
Netto 19.999[44] g
Tara 4.99941 g
Brutto 24.998[85] g
19.999[44] g
Datum 22.02.2016
Uhrzeit 13:08:11
Bediener
Produkt test
Netto 24.999[77] g
Tara 0.00000 g
Brutto 24.999[77] g
24.999[77] g
Datum 22.02.2016
Uhrzeit 13:09:00
Bediener
Produkt test
Netto 21.999[82] g
Tara 3.00000 g
Brutto 24.999[82] g
21.999[82] g
```

- ⇒ Any additional logs will be added to the same folder on re-inserting the USB medium.

15 Data bases

Available databases:

Databank	Data	editable	
➤ User <b1>	10 users	yes	see chap. 15.1.
➤ Product <b2>	1000 products	yes	see chap. 15.2.
➤ Tare values <b3>	100 container weights	yes	see chap. 15.3.
➤ Weighing results <b4.4>	1000 weighing results	no	see chap. 15.4.
➤ Alibi <b4.5>	100 000 weighing results (alibi memory)	no	

General data record operation for databases <b1> - <b3>:

Invoke database menu (For navigation in menu see chap. 8.2)

Add new data record to database

Delete selected data record in database

15.1 User database

The user database < **b1 Operator / User** > allows you to create up to 100 user profiles with specific settings and user rights.

The user management provides for users with different rights, the administrator and several users.

The **Administrator** has full access to the menu and database. He alone can create new users and grant individual user rights.

A user on the other hand cannot use all the functions and has the following limited rights that are specified in the user profile.

Authorisation level	Available rights and functions
USER	<ul style="list-style-type: none"> ⇒ Starting and implementing weighing processes ⇒ Defining universal variables ⇒ Exporting weighing data ⇒ Viewing data from the database ⇒ Access to the following functions: <ul style="list-style-type: none"> < P2.1.1 Read > < P6 Other> except <date / time>
MASTER USER (ADVANCED USER)	<ul style="list-style-type: none"> ⇒ Starting and implementing weighing processes ⇒ Access to the following functions: <ul style="list-style-type: none"> < P2 Operating modes> < P2.1.1 Read > < P3 Communication> < P4 Devices> < P6 Other> except <date / time>
ADMIN	<p>The administrator can use all the functions and has all rights.</p> <p>There is only one administrator.</p>

15.1.1 Create new user

The following data can be entered for each user:

NAME	max. 30 characters, CODE max. 6 characters
PASSWORD	max. 8 digits
ACCESS RIGHTS	User, master user, administrator
LANGUAGE	Deutsch/Espanol/Francais/Turk/Cesky/Italiano/Magyar/ Polski/English

The connection of a PC keyboard via USB interface makes the entry of user parameters easier.

⇒ Specify user name

⇒ After entering a user name specify user code, password, access right and language

64.1.1.1
NAME MILLER


64.1.14
ACCESS LEV USER


64.1.1.2
CODE

64.1.14
USER

64.1.1.2
1234

64.1.14
USER ADVANCED

64.1.1.2
CODE 1234

64.1.15
LANGUAGE ENGLISH


64.1.1.3
PASSWORD *****

64.1.15
DEUTSCH

64.1.1.3

64.1.15
LANGUAGE DEUTSCH

4 x

64.1.1.3
PASSWORD *****

0.000g
TARE 0000 g

15.1.2 Delete user

0.000g
TARE 0000 g

61
USER

3 x

64
VIEW/EDIT



64.1
USERS

64.1.1
ADMIN

- Use the navigation keys to select and delete a user by .

- Confirm the query “DELETE“ by pressing the PRINT-key and the display will return to the menu. Press UNITS key to return to weighing mode.

15.1.3 For invoking a user

See chap. 7.6.4 „Login“

15.2 Product database

For up to 1000 products in the product database < **b2 Produkt / Product** >: you can enter the following data











NAME	max. 30 digits
CODE	max. 6 digits
EAN	max. 16 digits
WEIGHT	Reference weight / reference single weight
TARA	Container weight
MIN	Define lower limit for check weighing
MAX	Define upper limit for check weighing
TOLERANCE	Define ± tolerance for target weight in dosing mode


For deleting product see chap. 15.1.2

15.2.1 Create new product

⇒ Define product name

⇒ After entering a product name define all parameters one by one


15.3 Container weights

The product database < **b3 TARA** > allows for up to 10 container weights to be created

you can enter the following data

NAME	max. 30 digits
TARA	Container weight

15.3.1 Create new container weight

- For downloading container weight from database see chap. 9.1.3
- For deleting container weight see chap. 15.1.2

15.4 Weighing results

The database < **b4.4 Weighing results** > is not editable. Data will be saved automatically when you press the PRINT-key. These data records can be viewed, printed or exported to a USB storage medium by the user.

The database allows for up to 1000 weighing results to be saved. When the capacity of the memory is exhausted, data records will be overwritten in turn, starting with the first data record.

The following additional data are saved together with the weighing result:


- Date of measurement
- Time of measurement
- Weighing results
- Container weight
- Product Name
- Logged-in user
- Enabled work mode
- Value of variables one and two


Data records of weighing results cannot be deleted.

Invoke data record:

- ⇒ The first weighing result <4.4.1> including date / time will be shown. Each weighing result will be saved under a consecutive number on form <b4.4.n>. Use the navigations keys to select weighing results.

- ⇒ To access further information on the selected weighing result, press the F2-key. Use the navigation keys to select and view the relevant data. Selecting <Print> followed by confirmation with the PRINT-key will trigger a printout at the optional printer. The content of the data output depends on the settings in the menu <P5.3 GLP>, see chap. 14.2.2.

Printout example (KERN YKB-01N):

Date	18.05.2015
Time	11:12:30
User	Admin
Product	Test
Net	23.998 g
Tare	3.000 g
Gross	26.998 g
	23.998 g

- ⇒ Back to menu / press weighing mode **UNITS**-key repeatedly.

15.5 Alibi memory

Within the framework of consumer protection the Weights and Measures Act stipulates electronic archiving in the form of a data memory that is appropriate for verification and cannot be manipulated with respect to weighing results subject to verification. Alibi memories by KERN meet this requirement.

They are designed for paperless recording of up to 100 000 weighing results.

All data transmitted to the PC will be saved including date, time and all the important weighing values. These saved data records are available for viewing on the weighing balance at any time.

The database < **b4.5 Weighing results** > is not editable. Data will be saved automatically when you press the PRINT-key. The user can view and print these data records. When the capacity of the memory is exhausted, data records will be overwritten in turn, starting with the first data record.

Weighing results are saved including the following data:


- Date of measurement
- Time of measurement
- Weighing results
- Container weight
- Product Name
- Logged-in user


Data records of weighing results cannot be deleted.

Invoke data record:

- ⇒ The first weighing result <4.4.1> including date / time will be shown. Each weighing result will be saved under a consecutive number on form <b4.4.n>. Use the navigations keys to select weighing results.

- ⇒ To access further information on the selected weighing result, press the F2-key. Use the navigation keys to select and view the relevant data. Selecting <Print> followed by confirmation with the PRINT-key will trigger a printout at the optional printer. The content of the data output depends on the settings in the menu <P5.3 GLP>, see chap. 14.2.2.

Printout example (KERN YKB-01N):

Date	18.05.2015
Time	11:12:30
Result	23.998 g
Tare	3.000 g
User	Admin
Product	Test


- ⇒ Back to menu / press weighing mode **UNITS**-key repeatedly.

15.6 Export / import databases to USB storage medium


This option allows to:

- Copy / transmit data of the databases “Products”, “User” and “Tare” from one weighing balance to another.
- Archive the data of the databases “Weighing Results” and “Alibi”

⇒ Connect a USB storage medium to an USB port

⇒ Use the navigation keys to select the options <Import> or <Export>


⇒ Use the navigation keys to select the database to be imported / exported and to confirm the selection by pressing the **PRINT**-KEY, available options:

OPERATOR
 PRODUCTS
 TARE
 WEIGHINGS*
 ALIBI*
 USER PARAMETERS
 ALL DATABASES

⇒ Back to menu / press weighing mode **UNITS**-key repeatedly.

* = It is not possible to import data records from the “weighing results” or the “alibi memory”

⇒ Connect the USB storage medium to the PC.

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- To open the “alibi” file you need the optionally available KERN PET-A16 software.
- It is not possible to open the „NC“-files as they merely serve the purpose of import to additional devices.


16 Application settings


Overview of available applications:

Weighing	see chap. 9.
Parts counting	see chap. 17.
Checkweighing	see chap. 18.
Dosing	see chap. 19.
Percent determination	see chap. 20.
Animal weighing	see chap. 21.
Density determination	see chap. 22.
Statistics	see chap. 23.
Totalization	see chap. 24.
Peak value function	see chap. 25.
Pipette calibration	see chap. 26.

16.1 Selecting an application


⇒ Press the **MODE**-key and the first function <F1 WEIGHING> will be shown.

⇒ Use the navigation keys   to select the required application and confirm by .

The factory-setting for the MODE-key makes all applications available. However, the list of available applications can be personalised by the user. That means, only applications set to “YES” will be shown when the MODE-key is pressed.

Enable / disable applications

17 Parts counting

Before the balance can count parts, it must know the average part weight (i.e. reference). The balance determines the total weight and divides it by the number of parts (the so-called reference quantity). Based on the calculated reference, single weight counting will follow next.

As a rule:

The higher the reference quantity the higher the counting exactness.

Two different methods are available for calculating the reference (See chap. 17.2.1 and chap.17.2.2).

17.1 Settings <P2.3> PARTS COUNTING>

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.

A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2

Menu item	Parameters	Explanation
P2.3.4 INFORM	There is an additional info box called “REFERENCE WEIGHT” available for Parts Counting	
P2.3.6 SHORTCUTS Functions keys F1; F2; F3; F4	The following additional function keys are available for Parts Counting	
	➤ ENTER SAMPLE	Enter reference single weight as numeric value
	➤ DETERMINE SAMPLE	Determine reference single weight by weighing
P2.3.7 ACAI	YES	Enable / disable automatic reference optimisation This function allows you to define whether the reference single weight is to be optimised automatically during counting. At every reference optimisation, the reference piece weight is calculated anew. As the additional pieces increase the base for the calculation, the reference also becomes more exact.
	NO:	Notes: Reference optimisation will only take place when: <ul style="list-style-type: none"> • The number of placed parts is greater than the previous reference number of parts • The number of placed parts is not greater than double the last reference number of parts • The weighing value is stable

17.2 Carry out parts counting

Preparation

For assigning two function keys of your choice with the functions <ENTER SAMPLE> and <DETERMINE SAMPLE> see chap. 17.1 or chap. 9.3.10.

17.2.1 Determination of the reference piece weight by weighing

- ⇒ Select application <PARTS COUNTING>, see chap.16.1. The reference single weight last used will be shown in the informational line of text. At the first start-up as value 0.0000g.

- ⇒ Press function key assigned to < DETERMINE SAMPLE > function, (Display example F1)

- ⇒ Use the navigation keys to select the desired reference number of parts, selectable 10, 20, 50 or user-defined.

Selecting reference number of parts 10, 20, 50 or 100:


- ⇒ Confirm selection (such as 100 items) by pressing the PRINT-key.

- ⇒ If necessary, tare the empty balance container.

- ⇒ Place as many pieces to add-up as required by the set reference piece number. Wait for steady state display and confirm by pressing the PRINT key.

- ⇒ The weighing balance will calculate the reference single weight and show it in the informational line of text. Remove reference weight.


- ⇒ The balance is now in parts counting mode counting all units on the weighing plate.


i The reference single weight will remain saved until you define a new reference.

Return to weighing mode:

Selecting reference number of parts user-defined (such as 30 items):

- ⇒ If necessary, tare the empty balance container.
- ⇒ Place as many pieces to add-up as required by the set reference piece number. Wait for steady state display and confirm by pressing the PRINT key.

- ⇒ The weighing balance will calculate the reference single weight and show it in the informational line of text. Remove reference weight.


- ⇒ The balance is now in parts counting mode counting all units on the weighing plate.


i The reference single weight will remain saved until you define a new reference.

17.2.2 Enter known reference single weight as numeric value

- ⇒ Select application <PARTS COUNTING>, see chap.16.1. The reference single weight last used will be shown in the informational line of text. At the first start-up as value 0.0000g.

- ⇒ Press the function key assigned to the <ENTER REFERENCE WEIGHT> function (Display example F2). Enter known reference single weight with the help of the navigation keys and the active digit will be flashing.

- ⇒ Confirm input by pressing the **PRINT** button.

- ⇒ The saved reference single weight will be shown in the informational line of text. Tare if necessary, place weighing good and read off the number of items.

- The reference single weight will remain saved until you define a new reference.
- Return to weighing mode, see chap. 17.2.1

17.2.3 Invoking reference single weight from the product database

Preparation

For assigning a function key of your choice to the < SELECT PRODUCT> function see chap. 17.1

- ⇒ Press the function key assigned to the <SELECT PRODUCT> function (Display example F1)

- ⇒ Use the navigation keys to select the desired product and confirm by pressing the PRINT-key.

- ⇒ The saved reference single weight will be shown in the informational line of text.
- ⇒ Tare if necessary, place weighing good and read off the number of items.

- i**
- For saving a product to the product database see chap. 15.2.1.
 - Return to weighing mode, see chap. 17.2.1

18 Checkweighing

In checkweighing an upper and a lower limit value (lower / upper threshold) can be established to ensure that the weighed goods are exactly within the established tolerance limits.

Any falling short of or exceeding the limits will be indicated on the screen by the indicators <MIN>, <OK> or <MAX>.

18.1 Settings <P2.4> CHECKWEIGHING

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.

A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2


Menu item	Parameters	Explanation
P2.4.4 INFORM		For check weighing, the additional info box <THRESHOLD> is available
P2.4.6 SHORTCUTS Functions keys F1; F2; F3; F4		For check weighing, the additional configuration <ENTER THRESHOLD> is available

18.1.1 Entering limits with a function key


- ⇒ Select application <CHECKWEIGHING>, see chap.16.1. The most recently used limits will be shown in the informational line of text. At the first start-up as value 0.0000g.

- ⇒ Press the function key assigned to the <ENTER THRESHOLD> function (display example F1). The display for entering the lower limit value will appear.

- ⇒ Use the navigation keys to enter the lower limit and confirm by pressing the PRINT-key. The display for entering the upper limit value will appear.

- ⇒ Use the navigation keys to enter the upper limit and confirm by pressing the PRINT-key.

The saved limits will be shown in the informational line of text (provided this option is set for the informational line of text, see chap. 18.1).

After setting the limit values, the weighing balance will be ready for check weighing.

- ⇒ Tare, as required, place the load on the weighing pan and start tolerance control.
The indicators <MIN>, <OK> or <MAX> will show the sample weight in relation to the tolerance limits.

Sample display	State
	Load below set tolerance
	Load within set tolerance
	Load above set tolerance

- The limit values will be saved until new values are defined
- To delete limit values, set value to [0.000 g]
- Return to weighing mode, see chap. 17.2.1

18.1.2 Invoke limits from the product database

Preparation

For assigning the function < SELECT PRODUCT > to a function key see chap. 17.1

- ⇒ Press the function key assigned to the <SELECT PRODUCT> function (Display example F1)

- ⇒ Use the navigation keys to select the desired product and confirm by pressing the PRINT-key.

- ⇒ The saved limits will be shown in the informational line of text (provided this option is set for the informational line of text, see chap. 18.1).
- ⇒ Tare, as required, place the load on the weighing pan and start tolerance control.
The indicators <MIN>, <OK> or <MAX> show the sample weight in relation to the tolerance limits, see chap. 18.1.1.

- i**
- For saving a product to the product database see chap. 15.2.1.
 - Return to weighing mode, see chap. 17.2.1

19 Dosing

During dosing you may for instance fill a container up to a specified target weight or a sample to a specified target weight.

The indicators <MIN>, <OK> or < MAX> show the filling status.

The additional entry defining \pm tolerance can be used to specify the degree of accuracy to be applied.

19.1 Settings <P2.5> DOSING

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.

A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2


Menu item	Parameters	Explanation
P2.5.4 INFORM		For dosing the additional info box <REFERENCE WEIGHT> is available
P2.5.6 SHORTCUTS Functions keys F1; F2; F3; F4		For dosing the additional configuration < ENTER REFERENCE WEIGHT> is available

19.1.1 Enter target weight with the function key

- ⇒ Select application <DOSING>, see chap.16.1. The most recently used target weight will be shown in the informational line of text. At the first start-up as value 0.0000g.

- ⇒ Press the function key assigned to the <ENTER REFERENCE WEIGHT> function (Display example F1). The display to enter the target weight appears.

- ⇒ Use the navigation keys to enter target weight and confirm by pressing the PRINT-key. The display to enter the tolerance interval.

- ⇒ Use the navigation keys to enter tolerance and confirm by pressing the PRINT-key.

The negative value of the target weight will be shown.

After setting the target weight / \pm tolerance the weighing balance will be ready for dosing.

- ⇒ Tare, as required and dose weighing load.
The indicators <MIN>, <OK> or <MAX> show the filling status.

Sample display	State
	Weighing load below the set target weight / tolerance
	Target weight reached and within the set tolerance
	Load above set target weight / tolerance

- The target weight / \pm tolerance will be saved until new values are defined.
- To delete, set value to [0.000 g].
- Return to weighing mode, see chap. 17.2.1

19.1.2 Invoke target weight from product database

Preparation

For assigning function keys of your choice to the < SELECT PRODUCT > function see chap. 17.1

- ⇒ Press the function key assigned to the <SELECT PRODUCT> function (Display example F1)

- ⇒ Use the navigation keys to select the desired product and confirm by pressing the PRINT-key.

- ⇒ The negative value of the target weight will be shown
- ⇒ Tare, as required and dose weighing load.
The indicators <MIN>, <OK> or <MAX> show the filling status, see chap.19.1.1.

- For saving a product to the product database see chap. 15.2.1.
- Return to weighing mode, see chap. 17.2.1

20 Percent weighing

Percent weighing allows to display weight in percent, in relation to a reference weight (100%).

Two different methods are available for calculating the reference (See chap. 20.2.1 and chap.20.2.2).

20.1 Settings <P2.6> DEVIATION>

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.

A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2

Menu item	Parameters	Explanation
P2.6.4 INFORM	For percentage weighing, the additional info box <REFERENCE WEIGHT> is available	
P2.6.6 SHORTCUTS Functions keys F1; F2; F3; F4	The following additional function keys are available for Parts Counting	
	➤ ENTER REFERENCE WEIGHT	Enter reference weight as numeric value
	➤ DEFINE REFERENCE WEIGHT	Determine reference weight by weighing

20.2 Perform percentage weighing

Preparation


For assigning two function keys of your choice with the functions <ENTER REFERENCE> and <DEFINE REFERENCE WEIGHT> see chap. 20.1 or chap. 9.3.10.

20.2.1 Determination of the reference weight by weighing


- ⇒ For selecting application <DEVIATION> see chap.16.1. The reference weight used will be shown in the informational line of text. At the first start-up as value 0.0000g.

- ⇒ Press the function key assigned to the <DEFINE REFERENCE WEIGHT> function (display example F1)


- ⇒ Put down a reference weight equalling 100 %, wait until the stability sign is shown on the weighing balance screen and confirm by pressing the PRINT-key.

The weight will be adopted as reference (100%) and shown in the informational line of text.

- ⇒ Tare, as required and place load on pan. The weight of the item to be weighed is displayed in percent, referring to the reference weight.


- The reference weight will remain saved until a new reference is defined.
- Return to weighing mode, see chap. 17.2.1

20.2.2 Enter known reference weight as numeric value


- ⇒ For selecting application <DEVIATION> see chap.16.1. The reference weight used will be shown in the informational line of text. At the first start-up as value 0.0000g.


- ⇒ Press the function key assigned to the <ENTER REFERENCE WEIGHT> function (Display example F2). Use the navigation keys to enter known reference weight and the active digit will be flashing.

- ⇒ Confirm input by pressing the **PRINT** button. The value will be accepted as reference (100%) and shown in the informational line of text.

- ⇒ Tare, as required and place load on pan. The weight of the item to be weighed is displayed in percent, referring to the reference weight.

- The reference weight will remain saved until a new reference is defined.
- Return to weighing mode, see chap. 17.2.1

20.2.3 Invoke the reference weight from the product database


Preparation

For assigning function keys of your choice to the < SELECT PRODUCT > function see chap. 17.1

- ⇒ Press the function key assigned to the <SELECT PRODUCT> function (Display example F1)

- ⇒ Use the navigation keys to select the desired product and confirm by pressing the PRINT-key.

The saved reference single weight will be shown in the informational line of text.

- ⇒ Tare, as required and place load on pan. The weight of the item to be weighed is displayed in percent, referring to the reference weight.

- For saving a product to the product database see chap. 15.2.1.
- Return to weighing mode, see chap. 17.2.1

21 Animal weighing

The Animal weighing application facilitates weighing of unstable loads or weighing under unstable ambient conditions. The weighing balance calculates the mean value of several weighing results obtained over a set time.

Two different start / reset modes are available.

- **Manual** (via keystroke, see chap. 21.1)
or
- **Automatic** (automatic start after exceeding a minimum weight, see chap. 21.2).

21.1 Animal weighing by manual start

Preferences:

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application. A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”). The description below merely deals with settings deviating from the default settings.

Menu item	Adjustment	Description
<P2.9.7 AUTOSTART>	<NO>	The measuring cycle is started manually, that is, placing of load must be confirmed by keystroke.
<P2.9.7 MEAN. Z.>	selectable 1 – 90 s	Duration after which the weighing balance calculates the mean value of several weighing results.
<P2.9.4> SHORTCUTS Functions keys F1; F2; F3; F4	<START>	Assign key with “Measuring cycle manual start” function

0.0000g
TARE 0000 G


P2.9.5
AVER TIME 5

P1
CALIBRATION

P2.9.5
05

P2
WORKING MODE

P2.9.5
10

P2.1
ACCESSIBILITY

P2.9.5
AVER TIME 10

P2.9
ANIMAL WEIGH

P2.9.4
HOZ KEYS

P2.9.1
READING

P2.9.4.1
F1 NONE

P2.9.7
AUTOSTART YES

P2.9.4.1
NONE

P2.9.7
YES

P2.9.4.1
START

P2.9.7
NO

P2.9.4.1
F1 START

4x

P2.9.7
AUTOSTART NO

2x



0.0000g
TARE 0000 G

The weighing balance will now be ready for animal weighing:


- ⇒ Select application < ANIMAL WEIGHING >, see chap.16.1.
When using a weighing container tare with **TARE**-key.

- ⇒ Place load on pan and start the measuring cycle with the function key assigned to the <START> function (display example F1)

- ⇒ Wait until the result is displayed. Remove load and press the UNITS-key. The weighing log will be printed automatically if an optional printer is connected (printout examples see chap. 14.2.2).
Then the balance is ready for further measurements.

21.2 Animal weighing with automatic start

Previous adjustments

Menu item	Adjustment	Description
<P2.9.7 AUTOSTART>	<YES>	The measuring cycle will start automatically as soon as the load is placed on the pan and the weight exceeds the set threshold (minimum weight).
<P2.9.7 MEAN. Z.>	selectable 1 – 90 s	Duration after which the weighing balance calculates the mean value of several weighing results.
<P2.9.7 THRESHOLD>	User-defined	Definition of minimum weight that starts the measuring cycle automatically when exceeded.

 0.000g
TARE 0000 G


P2.9.6
THRESHOLD 10.000

P1
CALIBRATION

P2.9.6
00 10.000


P2
WORKING MODE

P2.9.6
0020.000

P2.1
ACCESSIBILITY

P2.9.6
THRESHOLD 20.000

P2.9
ANIMAL WEIGH


P2.9.5
AVER TIME 5

P2.9.1
READING

P2.9.5
05

P2.9.7
AUTOSTART NO

P2.9.5
10

P2.9.7
NO

P2.9.5
AVER TIME 10

3x

P2.9.7
JA


 0.000g
TARE 0000 G

P2.9.7
AUTOSTART JA



The weighing balance will now be ready for animal weighing:


- ⇒ Select application < ANIMAL WEIGHING >, see chap.16.1.
When using a weighing container tare with **TARE**-key.

- ⇒ Place load on pan and the measuring cycle will be started automatically when the set threshold is reached.

- ⇒ Wait until the result is displayed. Remove load and press the UNITS-key. The weighing log will be printed automatically if an optional printer is connected (printout examples see chap. 14.2.2).
Then the balance is ready for further measurements.

- Return to weighing mode, see chap. 17.2.1

22 Determining the density of solid matter and fluids


For density determination we recommend working with the optionally available density determination sets.

1. Analytical balances

Density determination set KERN YDB-03

2. Precision balances

Density determination set KERN PLS-A01

The sets contain all the accessories and aids required for easy and precise density determination. For instructions please see the operating instructions enclosed with the density determination set.

22.1 Settings <P2.7> or <P2.8> DENSITY DETERMINATION

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.

A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2

Menu item	Parameters	Explanation
<P2.7> or <P2.8> SHORTCUTS Function keys F1; F2; F3; F4	For density determination, the following additional function key is available	
	➤ START	Start density determination

23 Statistics

The Statistics application facilitates the statistical evaluation of weighing values.

Settings <P2.10> STATISTICS >

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application. A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”). The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2

Menu item	Parameters	Explanation
P2.10.4 INFORM	For statistics, the following info boxes are available	
	N	Number of samples
	SUM	The added up weight of all single weights
	AVG	Average weight of all samples
	MIN	Smallest weight value
	MAX	Greatest weight value
	SDV	Standard Deviation
	DIF	Difference between smallest and greatest weight value
	SDV	Standard Deviation
	RDV	Relative standard deviation


Selecting an application

⇒ Select application <F 9 STATISTICS>, see chap.16.1.

Weight first sample

- ⇒ If required, set weighing balance to zero or tare. Place first sample on pan, wait for stability sign, then press the PRINT-key.


Number of sample <1> will appear on the display; the weight will be saved and then printed if an optional printer is connected.

- ⇒ Remove the first sample from the weighing pan.


Weigh additional samples

- ⇒ Follow the same procedure as for the first sample


Display / print result

- ⇒ After weighing the last sample press the F-key. The informational line of text will show <RESULT>.

⇒ After having pressed the F2-key you can display the next statistics value each time you press the navigation keys.


⇒ To activate log output, confirm option <PRINT> by pressing the PRINT-key.

Sample log KERN YKB-01N

-----Statistics-----		
N	9	Number of samples
Sum	455.60 g	Sum
Avg	50.6222 g	Mean Value
Min	49.939 g	Smallest weight value
Max	51.380 g	Greatest weight value
Dif	1.441 g	Difference smallest / greatest weight value
Sdv	0.39605 g	Standard deviation
Rdv	0.78 %	Relative standard deviation

Deleting statistical values

24 Totalization

The Statistics application allows you to weigh various samples and to add up the weighing values. The total number of samples that can be weighed is limited to 30.

Settings <P2.11> ADD UP >


Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application. A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”). The description below merely deals with settings deviating from the default settings.

i Navigation in menu see chap. 8.2

Menu item	Parameters	Explanation
P2.11.2 INFORM	For adding up, the following additional info box is available	
	➤ RESULT	N = Number samples S = Total
P2.11.4 ABBREVIATIONS Functions keys F1; F2; F3; F4	For adding up, the following additional function keys are available	
	➤ FINISH	Finish adding up process
	➤ DELETE LAST	Undo last weighing

Selecting an application

- ⇒ Select application <F10 ADD UP>, see chap.16.1. The informational line of text will show the sample counter “N” and the total “S”.

Weight first sample


- ⇒ If required, set weighing balance to zero or tare. Place first sample on pan, wait for stability sign, then press the PRINT-key.

The informational line of text will show the number of the sample < 1 >; the weight will be saved. The result on the display will be tared automatically.

Weigh additional samples

- ⇒ Follow the same procedure as for the first sample


When the sample number is displayed you can use the function key assigned to the <DELETE LAST> to undo the last weighing process.

Displaying / printing total

Either

⇒ After weighing the last sample press the F-key. Confirm query <FINISH> by pressing the F2-key.

The total will be displayed and then printed if an optional printer is connected.

or

⇒ Press the function key assigned to the function <FINISH>. The total will be shown and then printed if an optional printer is connected.

Sample log KERN YKB-01N


-----Adding-up-----	
1.	38.000 g
2.	100.00 g
3.	50.000 g
4.	10.000 g
5.	125.000 g
6.	15.100 g
7.	148.000 g
8.	6.000 g
9.	41.000 g
10.	15.000 g

Sum	548.100 g
Tare	100.00 g

Delete totals memory / start new adding up process

- ⇒ Remove sample, press UNIT-key. The sample counter “N” and the total “S” will be reset to zero.

The weighing scale is ready for a new adding up process.

25 Peak value function

This function displays the highest load value (peak value) of a weighing. The peak value remains in the display until it will be deleted.












Preferences:

Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application.


A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”).

The description below merely deals with settings deviating from the default settings.

Menu item	Adjustment	Description
<P2.12.5 THRESHOLD>	User-defined	The measuring cycle will start automatically as soon as the load is placed on the pan and the weight exceeds the set threshold (minimum weight).

3x


English

The weighing balance is ready to calculate a peak value:

- ⇒ Select application < F11 MAX. HOLD >, see chap.16.1.
From that point onwards the weighing balance will register every weighing value greater than the set threshold.


- ⇒ Place load on pan and weighing will start automatically on reaching the set threshold.

The highest load value will be indicated by the “Max” indicator.

- ⇒ The peak value remains in the display until it is deleted by the UNIT button.
Then the balance is ready for further measurements.

- Return to weighing mode, see chap. 17.2.1

26 Pipette calibration

- ⇒ Before starting the weighing process, acclimatise the sample liquid, weighing container, pipette and pipette tip for at least 2 hours.
- ⇒ Use distilled water as a sample fluid.
- Ambient temperature 20 – 25°C, temperature fluctuations during calibration \pm 0.5°C
- Relative humidity 50 -75 %

26.1 Settings <P2.13 PIPETTE CALIBRATION >





Applications can be adapted to your requirements by creating application-related settings. The settings available depend on the currently enabled application. A part of these settings is global, that is, they are – apart from a few exceptions – identical to the weighing application (See chap 9.3 „Settings <P2.1 WEIGHING>”. The description below merely deals with settings deviating from the default settings.

Navigation in menu see chap. 8.2


Menu item	Parameters	Explanation
P2.13.4 SHORTCUTS Functions keys F1; F2; F3; F4	For pipette calibration, the following additional function key is available.	
	➤ START	Start pipette calibration
P2.13.5 AMOUNT TEST VOLUME	For devices with fixed volume select setting <1> (test volume = nominal volume). For devices with variable volume two or three volumes can be tested:	
	1	As a rule 10 % of the nominal volume [μl]
	2	As a rule 50 % of the nominal volume [μl]
	3	As a rule 100 % of the nominal volume [μl]
P2.13.6 AMOUNT MEASUREMENTS	Number of measurements per test volume, selectable 6 -20 measurements	
P2.13.7 TARE AUTO	YES	Automatic taring after test cycle Turn on / off
	NO:	

26.2 Procedure


Select application:


Enter test parameters:

- ⇒ Use the navigation keys to enter the current test temperature; the active digit will be flashing.
- ⇒ Confirm by pressing the PRINT-key and the display for entering humidity will be shown. Use the navigations keys to enter the current humidity; the active digit will be flashing.


- ⇒ Confirm by pressing the PRINT-key and the display for entering air pressure will be shown. Use the navigations keys to enter the current barometric air pressure; the active digit will be flashing.


- ⇒ Confirm by pressing the PRINT-key and the entry for the first test volume <V1> will be shown.

- ⇒ Use the navigation keys to enter the test volume <V1> (such as 100 µl); the active digit will be flashing.


This completes the entry for pipettes with fixed volume. For pipettes with variable volume select volume V2 and V3 as described below.


- ⇒ Confirm your entry by pressing the PRINT-key and the display for entering the second test volume <V2> will be shown. Use the navigation keys to enter the test volume <V2> (such as 500 µl), the active digit will be flashing.

- ⇒ Confirm your entry by pressing the PRINT-key and the display for entering the third test volume <V3> will be shown. Use the navigation keys to enter the test volume <V3> (such as 1000 µl), the active digit will be flashing.


- ⇒ Confirm your entry by pressing the PRINT-key and the display for performing pipette calibration will be shown. Follow the instructions in the informational line of text.


Performing a pipette test:

1. Measuring sequence for test volume V1:



- ⇒ Place weighing container on pan, wait for stability sign, then tare by pressing the **TARE**-key.


1. Draw up selected test volume <V1> into pipette and dispense liquid into weighing container.
2. Wait for stability sign and accept weighing value by pressing the PRINT-key.
3. For menu setting P2.13.7 AUTOTARA <YES> automatic taring will follow. For menu setting tare <No> before the next dosing, process by pressing the TARE-key.

- ⇒ Repeat this test cycle (steps 1 -3) as often as set in the menu <P2.13.6> (see chap. 9.3).
- ⇒ After the confirmation of the last test cycle, the systematic measurement deviation [**E_s**] and the default deviation [**S_r**] will be shown.

- ⇒ To start the second measuring sequence <V2 >, press the **PRINT** -key. To finish, press the ESC-key.


2. Measuring sequence for test volume V2:

- ⇒ Follow the same sequence of operations for the second measuring sequence but apply test volume <V2>wit

3. Measuring sequence for test volume V3:

- ⇒ Follow the same sequence of operations for the third measuring sequence but apply test volume <V3>
- ⇒ After pressing the **PRINT**-key and the display of the systematic measurement deviation [E_s] / default deviation [S_r] printing of log will take place via RS232 interface if an optional printer is connected (See sample printout below).
- ⇒ To start a new pipette test, press the F-key.
- ⇒ Return to weighing mode:

Printout example KERN YKB-01N


-----Pipette calibration-----		
Meas. No.	10	
Date	05.01.2016	
Time	11:34:17	
Temp.	20.0°C	
Humidity	50 %	
Pressure	1017hPa	
-----Tested volume 1000 µl-----		
1	1003 µl	
2	993 µl	
3	1013 µl	
4	1003 µl	
5	993 µl	
6	1003 µl	
7	1013 µl	
8	1053 µl	
9	1003 µl	
10	1023 µl	
Average Volume [Va]	1010 µl	<i>Arithmetic average</i>
Systematic error [ES]	1.03 %	<i>Systematic measurement deviation</i>
Random error [Sr]	17.7 µl	<i>Standard deviation</i>
----- Tested volume 5000 µl-----		
1	4966 µl	
2	4966 µl	
3	4966 µl	
4	4986 µl	
5	4976 µl	
6	4966µl	
7	4966 µl	
8	4976 µl	
9	4976 µl	
10	4976 µl	
Average Volume [Va]	4972 µl	<i>Arithmetic average</i>
Systematic error [ES]	0.56 %	<i>Systematic measurement deviation</i>
Random error [Sr]	7.0 µl	<i>Standard deviation</i>
----- Tested volume 10000 µl-----		

1	10033 µl	
2	10033 µl	
3	10033 µl	
4	10033 µl	
5	10043 µl	
6	10043 µl	
7	10043 µl	
8	10043 µl	
9	10043 µl	
10	10043 µl	
Average Volume [Va]	10039 µl	<i>Arithmetic average</i>
Systematic error [ES]	0.39 %	<i>Systematic measurement deviation</i>
Random error [Sr]	5.2 µl	<i>Standard deviation</i>

Signature	

English

27 Communication / setup menu < P3 >



Available interfaces:

1. COM 2 (RS-232)
2. COM 1 (RS-232)
3. USB 2 / type B
(connection to printer or PC).
4. USB 1 / type A
(for instance connection to USB storage medium or barcode scanner)
5. WLAN (factory option KERN PLJ-A06)




The interface has to be configured according to peripheral and desired function, see chap. 0 / 27.2 (the parameters of the USB port are not configurable).

27.1 Configuring a <P3.1 / P3.2> RS-232 interface


- ⇒ Go to Communication menu <P3> and select the interface <COM 1> or <COM 2> to be configured.

- ⇒ After selecting <COM 1> or <COM 2> use the navigation keys to select baud rate under communication parameter.


⇒ After selecting baud rate, select the next communication parameter <Parity>.

27.2 <P3.3> WLAN interface (factory option KERN PLJ-A06)

Call up menu:


Available settings <P3.3>


P3.3.1 - STATUS	(Configure / connected)
P3.3.2 - SELECT NET	(Search for network)
P3.3.3 - NETWORK PARAMETER	
P3.3.3.1 - NAME	(Network name)
P3.3.3.2 - PASSWORD	(Password for access to network)
P3.3.3.3 – CHANNEL No.	(Standard <AUTO>)
P3.3.3.4- - IP	(Address of weighing balance in network)
P3.3.3.5 - MASK	(subnet mask)
P3.3.3.6 - GATE	
P3.3.3.7 - PORT	(Standard 4000)
P3.3.3.8 – MAC ADDRESS	

27.2.1 P3.3.2 - SELECT NET

⇒ For invoking menu item < P3.3.2>, see chap. 27.2.

⇒ Use the navigation keys to select the desired network and confirm by pressing the **PRINT**-key. The password query will appear.

- ⇒ Enter password and confirm by pressing the **PRINT**-key.
Remark: The connection of a PC keyboard via USB interface makes data entry easier.
- ⇒ The weighing balance program will guide you automatically through the network parameters.
Confirm displayed network parameters by pressing the **PRINT**-key or change, as required; see chap. 0.
The weighing balance will return automatically to the menu <P3.3.2 – SELECT NETWORK>


- ⇒ The weighing balance will be connected automatically to the network.
After successful network connection the display  will appear.
- ⇒ Back to weighing mode, press **UNIT**-key repeatedly.

27.2.2 <P3.3.3> - NETWORK PARAMETER


The connection of a PC keyboard via USB interface makes the entry of the network parameters easier.


⇒ For invoking menu item < P3.3.3> see chap. 27.2.


⇒ Under menu item < P3.3.3.3> select channel <AUTO (YES / NO)> and confirm by pressing the **PRINT**-key. When connecting next time to the network, the setting <YES> will trigger a check that determines whether the channel on which the router is running was changed. If it was changed, the channel will be adapted automatically. To select the next menu item, press the **F1**-key.


⇒ Under menu item < P3.3.3.4> select the IP address for the weighing balance in the network and confirm by pressing the **PRINT**-key. To select the next menu item, press the **F1**-key.

⇒ Under menu item < P3.3.3.5> select the sub-network mask <default 255.255.000.000> and confirm by pressing the **PRINT**-key. To select the next menu item, press the **F1**-key.


- ⇒ Under menu item < P3.3.3.6> select GATE and confirm by pressing the **PRINT**-key. To select the next menu item, press the **F1**-key.

- ⇒ Under menu item < P3.3.3.7> select PORT <default 4000> and confirm by pressing the **PRINT**-key. To select the next menu item, press the **F1**-key.

- ⇒ Under menu item < P3.3.3.8> select the MAC address and confirm by pressing the **PRINT**-key. Press the **F4**-key to return to menu. Back to weighing mode, press **UNIT**-key repeatedly.

27.3 USB interfaces

27.3.1 USB type A


The USB type A interface is used to connect a USB storage medium, barcode scanner or PC keyboard.

You can export / import weighing data with the help of a USB medium (See chap. 15.6) and measuring logs stored and printed by a printer connected to the PC (See chap.)


27.3.2 USB type B (PC connection)


Connecting a device via USB:

- ⇒ Disconnect the device from the power supply.
- ⇒ Install the USB driver which is available from our KERN- Homepage www.kern-sohn.com, / Downloads. Select the driver version compatible with your system and execute the exe file.

Click on Dialogue <Forward> and the driver will be installed.


To finish installation, click on <Finish>.

- ⇒ Connect the USB cable to the moisture meter and to the USB port on the computer.
Wait until the Windows Installer Assistant is shown.

- ⇒ Invoke device manager.
Click on „connections,, and the new virtual COM port for the device will be shown.
Correct entry is clearly indicated by the name of KERN PLJ series.


- ⇒ Select the shown COM Port such as COM6 according to the transmission software.
We recommend our transfer software 'Balance Connection KERN SCD 4.0' for the import of data to a PC program.
- ⇒ Go to menu **<P4.2.1 Devices / Printer / Port>** and select the setting **<USB-PC>**, see chap.28.2.

28 Devices / setup menu < P4 >

28.1 Continuous data transmission for remote control commands < P4.1 Computer >


1. Select interface <P4.1.1>

⇒ Use the navigation keys to select the interface to be applied for connecting the device to the PC.


Selectable:

- COM 1:** RS 232
- USB:** Type B
- WIFI:** Factory option KERN PLJ-A06


2. Select type of continuous data transmission <P4.1.2>

⇒ Use the navigation keys to select the display for setting continuous data transmission.


⇒ Use the navigation keys to select the desired setting and confirm by pressing the PRINT-key. Options:

- <NONE> Continuous data transmission disabled
- <ACT.UNIT> Continuous data transmission in current unit
- <CAL.UNIT> Continuous data transmission in adjustment unit


3. Select output interval <P4.1.3>

⇒ Use the navigation keys to select the display for setting the output intervals for continuous data transmission.


4. Select log type <P4.1.5>

⇒ Use the navigation keys to select the display for selecting the log type.

⇒ Use the navigation keys to select the desired setting and confirm by pressing the PRINT-key. Options:


- | | |
|---------------|---|
| <NONE> | None |
| <DEFAULT> | Standard Protocol |
| <DEF PRINT 1> | User-defined log: |
| ↓ | For special printout 1 -4 definition see chap. 14.3 |
| <DEFPRINT 4> | |

28.2 To start data transmission, press the PRINT-key < P4.2 Printer >

Go to menu item <P4.2> to select and configure the interface to be used for output after pressing the PRINT-key. The content for the data output is defined in < P5 Printout >, see chap. 14.



1. Select interface <P4.2.1>


⇒ Use the navigation keys to select the interface to be used for data transmission.

Options:

- <NONE> None
- <COM 1> Default printer interface RS 232
- <USB> USB port type B for EPSON- and PCL printer
- <WIFI> Not documented
- <USB stick> USB port type A, connection USB storage medium.
For saving measuring logs see chap. 0
- <USB PC> USB port type B, PC connection in combination with
transmission software (such as KERN SCD 4.0)



2. Enter values (hexadecimal) for prefix / suffix <P4.2.1>


⇒ Use the navigation keys to enter the control code for the start of the log and confirm by pressing the PRINT-key.

To select the display for entering the suffix, press the F1-key.


⇒ Use the navigation keys to enter the control code for the end of the log (such as paper section) and confirm by pressing the PRINT-key.


28.3 < P4.3 Barcode Scanner >

To select the interface for connection to a barcode scanner, go to menu item <P4.3>.

⇒ Use the navigation keys to select the interface to be used for data transmission. Confirm your entry by pressing the PRINT-key.

29 Communication with external devices (printer / PC)

The following options are available for data exchange via the RS 232C interface or USB with external devices.


- **Manual** - After printing the PRINT-key
- **Automatic** - After the stability sign has appeared on the display
- **Continuous** - When function is enabled or via remote control command

The following conditions must be met to provide successful communication between the balance and the external devices.

Disconnect balance from the power supply and connect to the device interface with a suitable cable. Faultless operation requires an adequate KERN interface cable.


For communication parameters of the RS 232 interface (baud rate, bits and parity) for the balance and the external device must match.

29.1 Pin assignment of RS 232 interface (frontal view):


29.2 Interface cable

⇒ Balance – PC 9-pole

⇒ Weighing balance - printer (CITIZEN, EPSON) 25-pin

⇒ Balance - printer

3 (TxD)	1 (RxD)
5 (GND)	3 (GND)
7 - 8 jumpered	

29.3 Data transmission format

1. After printing the PRINT-key

1	2	3	4-12	13	14-16	17	18
Stability	▬	Sign	Weighing Value	▬	Unit	CR	LF

Stability symbol	[▬] stable
	[?] unstable
	[^] excess load
	[v] underload
Signs	[▬] positive values
	[-] negative values
Measuring Value	9 characters aligned right
Unit	3 characters aligned left

Example (stable / positive weighing value):

<pre> ▬▬▬▬▬▬▬ 1 8 3 2 . 0 ▬g▬▬CR LF </pre>
--

2. Via remote control command

1-3	4	5	6	7	8-16	17	18-20	21	22
Command	▬	Stability sign	▬	Signs	Measuring Value	▬	Unit	CR	LF

Command	1 – 3 characters
Stability symbol	[_] stable
	[?] unstable
	[^] excess load
	[v] underload
Signs	[_] positive values
	[-] negative values
Measuring Value	9 characters aligned right
Unit	3 characters aligned left

29.4 Interface commands

The weighing balance recognises the commands listed below.

Commands	Function
Z	Zeroing
T	Taring
OT	Query tare value
UT	Set tare value
S	Send stable weighing value in standard weighing unit
SI	Send weighing value immediately in standard weighing unit
SIA	Send weighing values of all weighing pans in default weighing unit immediately
SU	Send stable weighing value in current weighing value
SUI	Send weighing value immediately in current weighing unit
C1	Start continuous output in standard weighing unit
C0	Stop continuous output in standard weighing unit
CU1	Start continuous output in current weighing unit
CU0	Stop continuous output in current weighing unit
DH	Enter lower tolerance limit (tolerance weighing)
UH	Enter upper tolerance limit (tolerance weighing)
ODH	Query lower tolerance limit (tolerance weighing)
OUH	Query upper tolerance limit (tolerance weighing)
SM	Enter reference single weight (parts counting)
TV	Enter target weight (dosing)
RM	Enter reference weight (percentage weighing)
Note	Enter serial number
K1	Disable device keyboard
K0	Enable device keyboard
SS	Confirm / Press simulation "PRINT-key"
IC	Start internal adjustment
IC1	Disable automatic internal adjustment
IC0	Enable automatic internal adjustment
OMI	List of available operating modes
OMS	Set operating mode
OMG	Query current operating mode
BP	Confirm sound signal
PC	Send all implemented commands

Finish commands with CR/LF character.

29.5 General reply format

Commands	Function
XX_A CR LF	Command accepted, command executing
XX_D CR LF	Command execution completed (occurs after XX_A only)
XX_I CR LF	Command valid but cannot be executed at the moment
XX_^ CR LF	Command valid but range exceeded
XX_v CR LF	Command valid but below range
XX_OK CR LF	Command accepted and executed
ES_ CR LF	Invalid entry
XX_E CR LF	Time limit for stabilisation of weighing balance display was exceeded

i	XX	Command such as Z = Reset to zero
	_	Space (20h, 0x20)

29.6 Detailed information to communication logs

29.6.1 Zeroing

Command: **Z CR LF**

Possible replies:

Z_A CR LF	Command accepted, command executing
Z_D CR LF	Command execution completed
Z_A CR LF	Command accepted, command executing
Z^ CR LF	Command valid but zero reset range was exceeded
Z_A CR LF	Command accepted, command executing
Z_E CR LF	Time limit for stabilisation of weighing balance display was exceeded
Z CR LF	Command valid but cannot be executed at the moment

29.6.2 Taring

Command: **T CR LF**

Possible replies:

T_A CR LF	Command accepted, command executing
T_D CR LF	Command execution completed
T_A CR LF	Command accepted, command executing
T^ CR LF	Command valid but tare range was exceeded
T_A CR LF	Command accepted, command executing
T_E CR LF	Time limit for stabilisation of weighing balance display was exceeded
T_I CR LF	Command valid but cannot be executed at the moment

29.6.3 Query tare value

Command: **OT CR LF**

Response:

1	2	3	4 - 12	13	14 - 16	17	18	19
O	T	▬	Tare	▬	Unit	▬	CR	LF

Tare value: - 9 characters right aligned

Unit: - 3 characters left aligned

i	The tare value will always be shown in the adjustment unit.
----------	---

29.6.4 Set tare value

Command: **UT_TARA CR LF**, (**TARA** = tare value)

Possible replies:

UT_OK CR LF	Command executed
UT_I CR LF	Command valid but cannot be executed at the moment
ES CR LF	Invalid entry

i	Setting tare value: <ul style="list-style-type: none"> • Mark decimal places with a point • Without weighing unit
----------	---

29.6.5 Send stable weighing value in standard weighing unit

Command: **S CR LF**

Possible replies:

S_A CR LF	Command accepted, command executing
S_E CR LF	Time limit for stabilisation of weighing balance display was exceeded
S_I CR LF	Command valid but cannot be executed at the moment

1	2-3	4	5	6	7-15	16	17-19	20	21
S	▬	Stability sign*	▬	Signs	Weight	▬	Unit	CR	LF

Example:

S CR LF	Sent command: Send stable weighing value in standard weighing unit
S _ A CR LF	Command accepted, command executing
S _ _ _ _ _ _ _ _ _ _ 8 . 5 _ g _ _ CR LF	Command executed, Weighing value is shown in default weighing unit

i	Stability symbol*	
	?	instable (3Fh, 0x3F)
	_	stable (20h, 0x20)

29.6.6 Send weighing value immediately in standard weighing unit

Command: **SI CR LF**

Possible replies:

SI _ CR LF	Command valid but cannot be executed at the moment
-------------------	--

1	2	3	4	5	6	7-15	16	17-19	20	21
S	I	_	Stability sign	_	Signs	Weight	_	Unit	CR	LF

Example:

SI CR LF	Sent command, see chap. 29.4
SI _ ? _ _ _ _ _ _ _ _ 1 8 . 5 _ k g _ CR LF	Command executed, Weighing value is shown in default weighing unit

29.6.7 Send stable weighing value in current weighing value

Command: **SU CR LF**

Possible replies:

SU_A CR LF	Command accepted, command executing
SU_E CR LF	Time limit for stabilisation of weighing balance display was exceeded
SU_I CR LF	Command valid but cannot be executed at the moment
SU_A CR LF	Command accepted, command executing

1	2	3	4	5	6	7-15	16	17-19	20	21
S	U	▬	Stability sign	▬	Signs	Weight	▬	Unit	CR	LF

Example:

S U CR LF	Sent command, see chap. 29.4
S U _ A CR LF	Command accepted, command executing
S U _ _ _ _ _ _ _ 1 7 2 . 1 3 5 _ N _ _ _ CR LF	Command executed. Weighing value is shown in current weighing unit

29.6.8 Send weighing value immediately in current weighing unit

Command: **SUI CR LF**

Possible replies:

SUI	Command valid but cannot be executed at the moment
------------	--

1	2	3	4	5	6	7-15	16	17-19	20	21
S	U	I	Stability sign	▬	Signs	Weight	▬	Unit	CR	LF

Example:

SUI CR LF	Sent command, see chap. 29.4
SUI ?_▬_▬_▬5 8 . 2 3 7_k g_▬ CR LF	Command executed. Weighing value is shown in current weighing unit

29.6.9 Start continuous output in standard weighing unit

Command: **C1 CR LF**

Possible replies:

C1	Command valid but cannot be executed at the moment
C1_A	Command accepted, command executing
	Weighing values are shown in default weighing unit

1	2	3	4	5	6	7-15	16	17-19	20	21
S	I	▬	Stability sign	▬	Signs	Weight	▬	Unit	CR	LF

29.6.10 Stop continuous output in standard weighing unit

Command: **C0 CR LF**

Possible replies:

C0	Command valid but cannot be executed at the moment
C0_A	Command accepted, command executing

29.6.11 Start continuous output in current weighing unit

Command: **CU1_***CR LF*

Possible replies:

CU1_I <i>CR LF</i>	Command valid but cannot be executed at the moment
CU1_A <i>CR LF</i>	Command accepted, command executing
	Weighing values are shown in current weighing unit

1	2	3	4	5	6	7-15	16	17-19	20	21
S	U	I	Stability sign	▬	Signs	Weight	▬	Unit	CR	LF

29.6.12 Stop continuous output in current weighing unit

Command: **CU0_***CR LF*

Possible replies:

CU0_I <i>CR LF</i>	Command valid but cannot be executed at the moment
CU0_A <i>CR LF</i>	Command accepted, command executing

29.6.13 Enter lower tolerance limit

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Possible replies:

DH_OK <i>CR LF</i>	Command executed
ES <i>CR LF</i>	Syntax error (incorrect weight format)

29.6.14 Enter upper tolerance limit

Command: **UH_XXXXX** *CR LF*

Possible replies:

UH_OK <i>CR LF</i>	Command executed
ES <i>CR LF</i>	Syntax error (incorrect weight format)

29.6.15 Query lower tolerance unit

Command: **ODH CR LF**

Response:

DH_MASA CR LF	Command executed
----------------------	------------------

1	2	3	4 - 12	13	14 - 16	17	18	19
D	H	_	Tare	_	Unit	_	<i>CR</i>	<i>LF</i>

Weight: - 9 characters right aligned

Unit: - 3 characters left aligned

29.6.16 Query upper tolerance unit

Command: **OUH CR LF**

Response:

UH_MASA CR LF	Command executed
----------------------	------------------

1	2	3	4 - 12	13	14 - 16	17	18	19
U	H	_	Tare	_	Unit	_	<i>CR</i>	<i>LF</i>

Weight: - 9 characters right aligned

Unit: - 3 characters left aligned

29.6.17 Enter reference single weight

Command: **SM_XXXXX CR LF**

Possible replies:

SM_OK CR LF	Command executed
SM_I CR LF	Command valid but cannot be executed at the moment (e. g. incorrect operating mode)
ES CR LF	Syntax error (incorrect weight format)

29.6.18 Enter target weight (dosing)

Command: **TV_XXXXX CR LF**

Possible replies:

TV_OK CR LF	Command executed
TV_I CR LF	Command valid but cannot be executed at the moment (e. g. incorrect operating mode)
ES CR LF	Syntax error (incorrect weight format)

29.6.19 Enter reference weight (percentage weighing)

Command: **RM_XXXXX CR LF**

Possible replies:

RM_OK CR LF	Command executed
RM_I CR LF	Command valid but cannot be executed at the moment (e. g. incorrect operating mode)
ES CR LF	Syntax error (incorrect weight format)

29.6.20 Confirm / simulation "Press PRINT-key"

Command: **SS CR LF**

Command **SS CR LF** is used to save the weighing result automatically to the database and if an optional printer is connected to print out the result.

29.6.21 Start internal adjustment

Command: **IC CR LF**

Possible replies:

IC_A CR LF	Adjustment running
IC_D CR LF	Adjustment finished
IC_E CR LF	Command valid but cannot be finished at the moment (e. g. range or time limit for stability display exceeded)
IC_I CR LF	Command valid but cannot be executed at the moment

29.6.22 Disable automatic internal adjustment

Command: **IC1 CR LF**

Possible replies:

IC1_OK CR LF	Command executed
IC1_I CR LF	Command valid but cannot be executed at the moment (Not available in weighing balances with type approval)
IC1_E CR LF	Command valid but not available at the moment

29.6.23 Enable automatic internal adjustment

Command: **IC0 CR LF**

Possible replies:

IC0 _OK CR LF	Command executed
IC0 _I CR LF	Command valid but cannot be executed at the moment (Not available in weighing balances with type approval)

29.6.24 Enter serial number

Command: **NB CR LF**

Possible replies:

NB_A "x" CR LF	Command accepted, stating serial number
NB_I CR LF	Command valid but cannot be executed at the moment

"x": **Serial number of device**

29.6.25 Disable weighing balance keyboard

Command: **K1 CR LF**

Possible replies:

K1_I CR LF	Command valid but cannot be executed at the moment
K1_OK CR LF	Command accepted, command executing

29.6.26 Enable weighing balance keyboard

Command: **K0 CR LF**

Possible replies:

K0_I CR LF	Command valid but cannot be executed at the moment
K0_OK CR LF	Command accepted, command executing

29.6.27 List of available operating modes

Command: **OMI** CR LF

Possible replies:

OMI _CR LF N° designation CR LF N° designation CR LF OK CR LF	Command accepted, available operating modes are shown
OMI _I CR LF	Command valid but cannot be executed at the moment

n	Operating mode
1	- Weighing
2	- Parts counting
3	- Percent determination
4	- Dosing
5	- Make up recipes
6	- Animal weighing
7	- Density determination "solid matter"
8	- Density determination "liquids"
9	- Maximum value function
10	- Totalization
11	- Tolerance weighing
12	- Statistics
13	- Pipette calibration

29.6.28 Set operating mode

Command: **OMS** CR LF

Possible replies:

OMS _OK CR LF	Command valid
OMS _E CR LF	No parameter or incorrect format
OMS _I CR LF	Command valid but cannot be executed at the moment

29.6.29 Query current operating mode

Command: **OMG** CR LF

Possible replies:

OMG _n_OK CR LF	Command accepted, number of current operating mode will be shown
OMG _I CR LF	Command valid but cannot be executed at the moment

n Number of operating mode, see chap. 29.6.27

29.6.30 Confirm sound signal

Command: **BP_CZAS** *CR LF*

Possible replies:

BP_OK <i>CR LF</i>	Command executed
BP_E <i>CR LF</i>	No parameter or incorrect format
BP_I <i>CR LF</i>	Command valid but cannot be executed at the moment

29.6.31 Send all messages implemented

Command: **PC** *CR LF*

Response:

PC_A_"Z, T, S, SI, SU, SUI, C1, C0, CU1, CU0, NB, K1, K0, SS, OMI, OMS, OMG, BP, PC "

30 Servicing, maintenance, disposal

Before any maintenance, cleaning and repair work disconnect the appliance from the operating voltage.

30.1 Cleaning

Please do not use aggressive cleaning agents (solvents or similar agents), but a cloth dampened with mild soap suds. Ensure that no liquid penetrates into the device. Polish with a dry soft cloth.

Loose residue sample/powder can be removed carefully with a brush or manual vacuum cleaner.

Spilled weighing goods must be removed immediately.

30.2 Servicing, maintenance

⇒ The appliance may only be opened by trained service technicians who are authorized by KERN.

⇒ Before opening, disconnect from power supply.

30.3 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

31 Troubleshooting / error messages

Possible causes of errors:

In case of an error in the program process, briefly turn off the balance and disconnect from power supply. The weighing process must then restart from the beginning.

Help:

Fault

Possible cause

The displayed weight does not glow.

- The balance is not switched on.
- The mains supply connection has been interrupted (mains cable not plugged in/faulty).
- Power supply interrupted.

The displayed weight is permanently changing

- Draught/air movement
- Table/floor vibrations
- Weighing pan has contact with other objects.
- Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)

The weighing result is obviously incorrect

- The display of the balance is not at zero
- Adjustment is no longer correct.
- The balance is on an uneven surface.
- Great fluctuations in temperature.
- Warm-up time was ignored.
- Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)

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