

Bench scales KERN ECE · ECB



Easy, mobile and economical

Features

- **High mobility:** Battery operation, compact, flat construction and low weight make this balance ideal for use in several locations
- **Particularly flat design**
- **Very fast display:** steady weight values within 2 sec
- **Simple and convenient 2-key operation**
- **Adjusting program CAL,** external test weights at an additional price

Technical data

- Overall dimensions WxDxH 320x300x60 mm
- Optional battery operation, 6 x 1.5 V AA not included, operating time approx. 100 h. AUTO-OFF function to preserve the battery
- Permissible ambient temperature 5°C / 35°C

ECB:

- Backlit LCD display, digit height 14 mm
- Dimensions of weighing plate (stainless steel) WxD 320x260 mm
- Net weight approx. 1,7 kg

ECE:

- Large backlit LCD display, digit height 14 mm
- Dimensions of weighing plate (plastic) WxD 320x260 mm
- Net weight approx. 1,5 kg

Options

- **Rechargeable battery pack internal,** operating time approx. 30 h, charging time approx. 10 h. AUTO-OFF function to preserve the battery, KERN PCB-A01
- **Rechargeable battery pack external,** operating time approx. 270 h, charging time approx. 10 h, KERN KS-A01

STANDARD



OPTION



STANDARD



OPTION



Model	Weighing range [Max] kg	Readout [d] g	Reproducibility g	Linearity g	Option DKD Calibr. Certificate	
					DKD KERN	
ECE 10K5	10	5	5	± 10	963-128	
ECE 20K10	20	10	10	± 20	963-128	
ECE 50K20	50	20	20	± 400	963-128	
ECB 10K5	10	5	5	± 10	963-128	
ECB 20K10	20	10	10	± 20	963-128	
ECB 50K20	50	20	20	± 400	963-128	

Improved readout

