

KERN 303-046

Set of weights, knob shape, polished stainless steel

KERN®



Piktogramme

STANDARD



OPTION



Measuring system

Nominal value: 1 mg | 200000 mg
OIML class: E1

Accessories

Aluminium box: 313-040-600

Services (optional)

DAkkS Certificate: 963-203
DAkkS Certificate (Recalibration): 962-203R

Category

Brand: KERN

Packaging & shipping

Delivery: 24 h
Dimensions packaging (WxDxH): 190 x 130 x 80 mm
Gross weight: 1,100 kg
Net weight: 900 g
Shipping method: Parcel service

Construction

Material: polished stainless steel
Shape: cylindrical