

# KERN 316-07

Single weights, compact shape, polished stainless steel



## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	22 x 22 x 35 mm
Gross weight:	100 g
Net weight:	100 g
Shipping method:	Parcel service

## Construction

Dimensions:	22   35
Material:	polished stainless steel
Shape:	compact shape

## Piktogramme

### STANDARD



### OPTION



## Measuring system

Nominal value:	100 g
OIML class:	E2
Tolerance:	0,160 mg

## Accessories

Aluminium box:	317-070-600
Plastic box:	317-070-400
Wooden box:	317-070-100

## Services (optional)

DAkkS Certificate:	962-337
DAkkS Certificate (Recalibration):	962-337R
Positive tolerance:	921-337
Verification:	952-337

## Category

Brand:	KERN
--------	------