

KERN 316-14

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	100 x 100 x 158 mm
Gross weight:	10 kg
Net weight:	10 kg
Shipping method:	Parcel service

Construction

Dimensions:	100 158
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1e+007 mg
OIML class:	E2
Tolerance:	16 mg

Accessories

Aluminium box:	317-140-600
Plastic box:	317-140-400
Wooden box:	317-140-100

Services (optional)

DAkKS Certificate:	962-344
DAkKS Certificate (Recalibration):	962-344R
Positive tolerance:	921-344
Verification:	952-344

Category

Brand:	KERN
--------	------