

# KERN 317-02

Single weights, cylindrical, polished stainless steel

**KERN**<sup>®</sup>



## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	6 x 6 x 10 mm
Gross weight:	2 g
Net weight:	2 g
Shipping method:	Parcel service

## Construction

Dimensions:	6   10
Material:	polished stainless steel
Shape:	cylindrical

## Piktogramme

### STANDARD



1 DAY

### OPTION



+3 DAYS

## Measuring system

Nominal value:	2 g
OIML class:	E2
Tolerance:	0,040 mg

## Accessories

Aluminium box:	317-020-600
Plastic box:	317-020-400
Wooden box:	317-020-100

## Services (optional)

DAkkS Certificate:	962-332
DAkkS Certificate (Recalibration):	962-332R
Positive tolerance:	921-332
Verification:	952-332

## Category

Brand:	KERN
--------	------