

# KERN 326-37

Rectangular weights, polished stainless steel



## Piktogramme

### STANDARD



### OPTION



## Measuring system

Nominal value:	10 kg
OIML class:	F1
Tolerance:	50 mg

## Accessories

Aluminium box:	346-070-600
Plastic box:	nicht verfügbar
Wooden box:	nicht verfügbar

## Services (optional)

DAkkS Certificate:	962-444
DAkkS Certificate (Recalibration):	962-444R

## Category

Brand:	KERN
--------	------

## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	193 x 96 x 117 mm
Gross weight:	10 kg
Net weight:	10 kg
Shipping method:	Parcel service

## Construction

Dimensions:	193 x 96 x 117 mm
Material:	polished stainless steel
Shape:	Rectangular weight