

KERN 326-38

Rectangular weights, polished stainless steel

KERN®



Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	20 kg
OIML class:	F1
Tolerance:	100 mg

Accessories

Aluminium box:	346-080-600
Plastic box:	nicht verfügbar
Wooden box:	nicht verfügbar

Services (optional)

DAkkS Certificate:	962-445
DAkkS Certificate (Recalibration):	962-445R

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	233 x 115 x 147 mm
Gross weight:	20 kg
Net weight:	20 kg
Shipping method:	Parcel service

Construction

Dimensions:	233 x 115 x 147 mm
Material:	polished stainless steel
Shape:	Rectangular weight