

KERN 327-01

Single weights, cylindrical, polished stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	6 x 6 x 6 mm
Gross weight:	1 g
Net weight:	1 g
Shipping method:	Parcel service

Construction

Dimensions:	6 6
Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 g
OIML class:	F1
Tolerance:	0,100 mg

Accessories

Aluminium box:	317-010-600
Plastic box:	347-030-400
Wooden box:	317-010-100

Services (optional)

DAkKS Certificate:	962-431
DAkKS Certificate (Recalibration):	962-431R
Positive tolerance:	921-431
Verification:	952-431

Category

Brand:	KERN
--------	------