

KERN 327-06

Single weights, cylindrical, polished stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	18 x 18 x 29 mm
Gross weight:	50 g
Net weight:	50 g
Shipping method:	Parcel service

Construction

Dimensions:	18 29,500
Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



1 DAY

OPTION



DAkKS

+3 DAYS

Measuring system

Nominal value:	50 g
OIML class:	F1
Tolerance:	0,300 mg

Accessories

Aluminium box:	317-060-600
Plastic box:	347-070-400
Wooden box:	317-060-100

Services (optional)

DAkKS Certificate:	962-436
DAkKS Certificate (Recalibration):	962-436R
Positive tolerance:	921-436
Verification:	952-436

Category

Brand:	KERN
--------	------