

KERN 327-12

Single weights, cylindrical, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	60 x 60 x 115 mm
Gross weight:	2 kg
Net weight:	2 kg
Shipping method:	Parcel service

Construction

Dimensions:	60 115
Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	2e+006 mg
OIML class:	F1
Tolerance:	10 mg

Accessories

Aluminium box:	317-120-600
Plastic box:	347-120-400
Wooden box:	317-120-100

Services (optional)

DAkkS Certificate:	962-442
DAkkS Certificate (Recalibration):	962-442R
Positive tolerance:	921-442
Verification:	952-442

Category

Brand:	KERN
--------	------