

KERN 327-13

Single weights, cylindrical, polished stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	80 x 80 x 147 mm
Gross weight:	5 kg
Net weight:	5 kg
Shipping method:	Parcel service

Construction

Dimensions:	80 147
Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5e+006 mg
OIML class:	F1
Tolerance:	25 mg

Accessories

Aluminium box:	317-130-600
Plastic box:	347-130-400
Wooden box:	317-130-100

Services (optional)

DAkKS Certificate:	962-443
DAkKS Certificate (Recalibration):	962-443R
Positive tolerance:	921-443
Verification:	952-443

Category

Brand:	KERN
--------	------