


KERN 328-06

Milligram weight, flat polygonal sheet, polished stainless steel



Measuring system

Nominal value:	50 mg
OIML class:	F1
Tolerance:	0,040 mg

Accessories

Aluminium box:	317-009-600
Plastic box:	347-009-400
Wooden box:	338-090-200

Services (optional)

DAkkS Certificate:	962-456
DAkkS Certificate (Recalibration):	962-456R
Positive tolerance:	921-456
Verification:	952-456

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	7,5 x 6 x 3 mm
Gross weight:	50 mg
Net weight:	50 mg
Shipping method:	Parcel service

Construction

Dimensions:	7,5 x 6 x 3 mm
Material:	polished stainless steel
Shape:	flat polygonal sheet

Piktogramme

STANDARD



OPTION

