

KERN 333-084

Set of weights, knob shape, finely turned stainless steel

KERN®



Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 mg 5000000 mg
OIML class:	F2

Accessories

Plastic box:	313-082-400
--------------	-------------

Services (optional)

DAkkS Certificate:	962-407
DAkkS Certificate (Recalibration):	962-407R
Positive tolerance:	921-407

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	350 x 320 x 170 mm
Gross weight:	12,690 kg
Net weight:	12,690 kg
Shipping method:	Parcel service

Construction

Material:	finely turned stainless steel
Shape:	cylindrical