

# KERN 334-02

Set of weights, knob shape, finely turned stainless steel

**KERN**<sup>®</sup>



## Piktogramme

### STANDARD



### OPTION



## Measuring system

Nominal value:	1 g   50 g
OIML class:	F2

## Accessories

Wooden box:	353-420-200
-------------	-------------

## Services (optional)

DAkkS Certificate:	962-415
DAkkS Certificate (Recalibration):	962-415R
Positive tolerance:	921-415

## Category

Brand:	KERN
--------	------

## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	120 x 50 x 40 mm
Gross weight:	260 g
Net weight:	260 g
Shipping method:	Parcel service

## Construction

Material:	finely turned stainless steel
Shape:	cylindrical