

# KERN 337-02

Single weights, cylindrical, finely turned stainless steel

**KERN**<sup>®</sup>



## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	6 x 6 x 11 mm
Gross weight:	2 g
Net weight:	2 g
Shipping method:	Parcel service

## Construction

Dimensions:	6   11
Material:	finely turned stainless steel
Shape:	cylindrical

## Piktogramme

### STANDARD



### OPTION



## Measuring system

Nominal value:	2 g
OIML class:	F2
Tolerance:	0,400 mg

## Accessories

Aluminium box:	317-020-600
Plastic box:	347-030-400
Wooden box:	337-020-200

## Services (optional)

DAkkS Certificate:	962-432
DAkkS Certificate (Recalibration):	962-432R
Positive tolerance:	921-432

## Category

Brand:	KERN
--------	------