

KERN 337-05

Single weights, cylindrical, finely turned stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	13 x 13 x 23 mm
Gross weight:	20 g
Net weight:	20 g
Shipping method:	Parcel service

Construction

Dimensions:	13 23,500
Material:	finely turned stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



1 DAY

OPTION



DAKkS
+3 DAYS

Measuring system

Nominal value:	20 g
OIML class:	F2
Tolerance:	0,800 mg

Accessories

Aluminium box:	317-050-600
Plastic box:	347-050-400
Wooden box:	337-050-200

Services (optional)

DAkKS Certificate:	962-435
DAkKS Certificate (Recalibration):	962-435R
Positive tolerance:	921-435

Category

Brand:	KERN
--------	------