

KERN 337-07

Single weights, cylindrical, finely turned stainless steel

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	22 x 22 x 39 mm
Gross weight:	100 g
Net weight:	100 g
Shipping method:	Parcel service

Construction

Dimensions:	22 39
Material:	finely turned stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	100 g
OIML class:	F2
Tolerance:	1,600 mg

Accessories

Aluminium box:	317-070-600
Plastic box:	347-070-400
Wooden box:	337-070-200

Services (optional)

DAkkS Certificate:	962-437
DAkkS Certificate (Recalibration):	962-437R
Positive tolerance:	921-437

Category

Brand:	KERN
--------	------