

KERN 344-08

Set of weights, knob shape, finely turned stainless steel, Wooden box

KERN®



Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 g 5000 g
OIML class:	M1

Accessories

Wooden box:	353-480-200
-------------	-------------

Services (optional)

DAkkS Certificate:	962-621
DAkkS Certificate (Recalibration):	962-621R
Verification:	952-621

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	400 x 190 x 200 mm
Gross weight:	17,810 kg
Net weight:	17,810 kg
Shipping method:	Parcel service

Construction

Material:	finely turned stainless steel
Shape:	cylindrical