


KERN 347-445-100

Beam bar, aluminium / finely turned stainless steel

KERN[®]



Measuring system

Nominal value:	10 g
OIML class:	M1
Tolerance:	2 mg

Services (optional)

DAkkS Certificate:	962-634
DAkkS Certificate (Recalibration):	962-634R

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	14 x 14 x 118 mm
Gross weight:	10 g
Net weight:	10 g
Shipping method:	Parcel service

Construction

Dimensions:	14 118
Material:	aluminium
Shape:	Beam bar

Piktogramme

STANDARD



OPTION

