

KERN 347-495-100



Beam bar, aluminium . Note: Brass weights while stocks last. Inexpensive alternative: stainless steel button weights. Stainless steel slotted and hook weights available soon.



Measuring system

Nominal value:	500 g
OIML class:	M1
Tolerance:	25 mg

Accessories

Plastic box:	347-515-400
--------------	-------------

Services (optional)

DAkKS Certificate:	962-639
DAkKS Certificate (Recalibration):	962-639R

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	40 x 40 x 645 mm
Gross weight:	500 g
Net weight:	500 g
Shipping method:	Parcel service

Construction

Dimensions:	40 645
Material:	finely turned brass
Shape:	Beam bar

Piktogramme

STANDARD



OPTION

