

# KERN 347-515-100-N00



Beam bar, aluminium . Note: Brass weights while stocks last. Inexpensive alternative: stainless steel button weights. Stainless steel slotted and hook weights available soon.



## Accessories

Plastic box: 347-515-400

## Services (optional)

DAkkS Certificate: 962-641N

DAkkS Certificate (Recalibration): 962-641RN

## Category

Brand: KERN

## Packaging & shipping

Delivery: 24 h

Dimensions packaging (WxDxH): 60 x 60 x 1030 mm

Gross weight: 1 kg

Net weight: 1 kg

Shipping method: Parcel service

## Construction

Dimensions: 60 | 1030

Material: finely turned brass

Shape: Slotted weight

## Piktogramme

### STANDARD



1 DAY

### OPTION



+3 DAYS

## Measuring system

OIML class: M1