

KERN 366-89

Rectangular weight, cast iron, lacquered

KERN®



Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5e+007 mg
OIML class:	M3
Tolerance:	25 g

Accessories

Aluminium box:	346-090-600
----------------	-------------

Services (optional)

DAkkS Certificate:	962-646
DAkkS Certificate (Recalibration):	962-646R

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	48 h
Dimensions packaging (WxDxH):	310 x 155 x 192 mm
Gross weight:	50 kg
Net weight:	50 kg
Shipping method:	Freight forwarder

Construction

Dimensions:	310 x 155 x 192 mm
Material:	cast iron lacquered
Shape:	Rectangular weight