

KERN 367-43

Single weights, cylindrical, finely turned brass/finely turned stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	8 x 8 x 14 mm
Gross weight:	10 g
Net weight:	10 g
Shipping method:	Parcel service

Construction

Dimensions:	8 14
Material:	finely turned brass
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5 g
OIML class:	M3
Tolerance:	16 mg

Accessories

Aluminium box:	317-030-600
Plastic box:	347-030-400
Wooden box:	337-030-200

Services (optional)

DAkkS Certificate:	962-633
DAkkS Certificate (Recalibration):	962-633R

Category

Brand:	KERN
--------	------