

# KERN 367-51

Single weights, cylindrical, finely turned brass/finely turned stainless steel



## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	48 x 48 x 80 mm
Gross weight:	1 kg
Net weight:	1 kg
Shipping method:	Parcel service

## Construction

Dimensions:	48   80
Material:	finely turned brass
Shape:	cylindrical

## Piktogramme

### STANDARD



### OPTION



## Measuring system

Nominal value:	1e+006 mg
OIML class:	M3
Tolerance:	500 mg

## Accessories

Aluminium box:	317-110-600
Plastic box:	347-110-400
Wooden box:	337-110-200

## Services (optional)

DAkkS Certificate:	962-641
DAkkS Certificate (Recalibration):	962-641R

## Category

Brand:	KERN
--------	------