

KERN 367-52

Single weights, cylindrical, finely turned brass

KERN[®]



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	60 x 60 x 100 mm
Gross weight:	2 kg
Net weight:	2 kg
Shipping method:	Parcel service

Construction

Dimensions:	60 100
Material:	finely turned brass
Shape:	cylindrical

Piktogramme

STANDARD



1 DAY

OPTION



+3 DAYS

Measuring system

Nominal value:	2e+006 mg
OIML class:	M3
Tolerance:	1 g

Accessories

Aluminium box:	317-120-600
Plastic box:	347-120-400
Wooden box:	337-120-200

Services (optional)

DAkkS Certificate:	962-642
DAkkS Certificate (Recalibration):	962-642R

Category

Brand:	KERN
--------	------