

# KERN ORA10BA

Area of application: Sugar



## Measuring system

ATC:	yes
Division:	0,100 %
Measuring capacity:	0 % - 10 %
Scales:	Brix

## Services (optional)

Factory calibration certificate:	961-290
----------------------------------	---------

## Category

Brand:	Optics
--------	--------

## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	215 x 94 x 65 mm
Gross weight:	320 g
Net weight:	267,720 g
Shipping method:	Parcel service

## Construction

Dimensions (WxDxH):	198 x 40 x 40 mm
---------------------	------------------

## Piktogramme

### STANDARD



### OPTION

