

KERN ORA3SA

Area of application: Salt



Measuring system

ATC:	yes
Division:	0.2 % 0.2 %
Measuring capacity:	0 % - 32 % 0 % - 28 %
Scales:	BrixSalt (NaCl)

Services (optional)

Factory calibration certificate:	961-290
----------------------------------	---------

Category

Brand:	Optics
--------	--------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	215 x 94 x 65 mm
Gross weight:	310 g
Net weight:	257,720 g
Shipping method:	Parcel service

Construction

Dimensions (WxDxH):	178 x 40 x 40 mm
---------------------	------------------

Piktogramme

STANDARD



OPTION

