

KERN ORA4RR

Area of application: Expert applications



Measuring system

Division:	0,001 nD
Measuring capacity:	1.440 nD - 1 520 nD
Scales:	Refractive index

Services (optional)

Factory calibration certificate:	961-290
----------------------------------	---------

Category

Brand:	Optics
--------	--------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	215 x 94 x 65 mm
Gross weight:	300 g
Net weight:	247,720 g
Shipping method:	Parcel service

Construction

Dimensions (WxDxH):	137 x 40 x 40 mm
---------------------	------------------

Piktogramme

STANDARD

