

# KERN ORA6HA

Area of application: Honey



## Measuring system

ATC:	yes
Division:	0,100 %
Measuring capacity:	12 % - 30 %
Scales:	Water content

## Category

Brand:	Optics
--------	--------

## Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	215 x 94 x 65 mm
Gross weight:	305 g
Net weight:	252,720 g
Shipping method:	Parcel service

## Construction

Dimensions (WxDxH):	168 x 40 x 40 mm
---------------------	------------------

## Piktogramme

### STANDARD



### OPTION

