

Product brochure

08.08.2026

KERN SET-311

EasyTouch Quick Key – select master data objects quickly and user-related


ENGLISH

Software EasyTouch SET-311 **Quick Key**

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- Note: Prerequisite for this set is the basic program SET-01 Base
- Quick Keys can be created in numerous sets, e. g.
 - SET-01 Base
 - SET-05 Dynamic
 - SET-06 Classification
 - SET-07 Variable
 - SET-08 Take Out
 - SET-31 Tolerance
 - SET-32 Count
 - SET-33 Target Count
- Quick selection keys for master data: Through ET Quick Keys, specific master data objects can be permanently displayed for quick selection. The Quick Keys are displayed on the screen immediately
- Saving: Whether selected in the conventional way or by using Quick Keys, every result of a weighing or piece count will be stored in relation to the selected master data object
- Selecting master data objects: The specific master data objects, which are displayed as Quick Keys, can either be selected individually or as a defined master data group, which contains the master data as Quick Keys
- User-related: The Quick Keys can be created as a scheme which relates to particular user roles. This means that the users of a particular role are shown the relevant master data objects for quick selection. This allows, for example, different departments in the same company to work with different master data without overlapping
- Saving and printing: For each scheme you can set whether, in addition to the being saved, it should also be automatically printed out or not
- Communication between balance/terminal device: The balance(s) can communicate with the PC, laptop or tablet by serial connection, USB, Bluetooth, Ethernet or WiFi

License information

- Licensing: One license can be operated on one terminal device (PC, laptop, Android-tablet)
- Balances: You can store and operate as many balances as you need in one license

PLEASE NOTE FOR YOUR ORDER

- Prerequisite for this set is the basic program SET-01 Base

PICTOGRAMS

Standard



Option



KERN Pictogram Overview



EasyTouch: Suitable for the connection, data transmission and control through PC or tablet



Central user administration for KERN EasyTouch functions. Primary master data, for example, first name, surname, image, is stored for a user. Each user is assigned a user role, which has individually configurable permissions in the KERN EasyTouch App. Through this structure you can create standard users, supervisor users, etc. You can log into the KERN EasyTouch App with a user name and a password.



All printed and stored weighing and measuring data is stored in this memory. The data is stored either on the local display device or centrally – for all connected weighing systems – in a server directory (Save-Data Server or Save-Data Cloud). Dynamic data can be recalled and printed out at any time and/or exported as a table.



All stored data is saved in a tamper-proof manner. Changes to master data are also saved in the dynamic data memory in a tamper-proof manner (Data Traceability)



Device qualification: standard-compliant qualification concept, which includes the following validation services: Installation qualification (IQ), function qualification (OQ)



Central data memory function for additional storage of all measuring data in a central, local server directory. This is where measuring data is stored from all weighing systems connected using KERN EasyTouch, as well as all installed KERN EasyTouch functions. The advantage of this, particularly for users with several weighing systems, is having all weighing data consolidated in just one database and only having to search for individual measuring data from different balances in one table. Save-Server data storage is also tamper-proof and cannot be changed.