

KERN 303-026

KERN

Weight set, E1, 1 mg - 50 g, Knob, Stainless steel polished, in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

Measuring System

Test weight - nominal value	1 mg - 50 g
OIML-compliant	✓
Weight OIML class	E1

Construction

Dimension (W×D×H)	165×120×75 mm
Shape weight	Knob
Material	Stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	963-201
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

Optional cases

Case - Aluminium article number	313-020-600
---------------------------------	-------------