

# KERN 303-036

**KERN**

Weight set, E1, 1 mg - 100 g, Knob, Stainless steel polished, in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Test weight - nominal value	1 mg - 100 g
OIML-compliant	✓
Weight OIML class	E1

## Construction

Dimension (W×D×H)	165×125×75 mm
Shape weight	Knob
Material	Stainless steel polished
Case of weight set	in aluminium case

## Services (optional)

Article number for DAkkS calibration	963-202
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×125×75 mm
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

## Optional cases

Case - Aluminium article number	313-030-600
---------------------------------	-------------