

# KERN 303-04

# KERN

Weight set, E1, 1 mg - 200 g, Knob, Stainless steel polished, in wooden case



## Optional cases

Case - Wood article number 313-040-100

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Test weight - nominal value	1 mg - 200 g
OIML-compliant	✓
Weight OIML class	E1

## Construction

Dimension (W×D×H)	152×152×66 mm
Shape weight	Knob
Material	Stainless steel polished
Number of weights	23
Case of weight set	in wooden case

## Services (optional)

Article number for DAkkS calibration	963-203
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,27 kg