

Product brochure

07.08.2026

KERN 303-08

Weight set, E1, 1 mg - 5 kg, Knob, stainless steel polished, in wooden case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-080-100	Wood weight case 313-080-100 (for nominal values 1 mg - 5 kg for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-716	Marking (Last digits of serial number) KERN 920-716
962-207R	DAkKS Recalibration 962-207R
963-207	DAkKS calibration certificate 963-207

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days