

KERN 303-086

KERN

Weight set, E1, 1 mg - 5 kg, Knob, Stainless steel polished, in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

Measuring System

Test weight - nominal value	1 mg - 5 kg
OIML-compliant	✓
Weight OIML class	E1

Construction

Dimension (W×D×H)	460×260×200 mm
Shape weight	Knob
Material	Stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	963-207
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	18 kg
Shipping weight	18 kg

Optional cases

Case - Aluminium article number	313-080-600
---------------------------------	-------------