

Product brochure

08.08.2026

KERN 304-026

Weight set, E1, 1 g - 50 g, Knob, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-020-600	Aluminium weight case 314-020-600 (for nominal values 1 g - 50 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-710	Marking (Last digits of serial number) KERN 920-710
921-215	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-215
962-215R	DAkKS Recalibration 962-215R
963-215	DAkKS calibration certificate 963-215

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days