

Product brochure

08.08.2026

KERN 304-046

Weight set, E1, 1 g - 200 g, Knob, stainless steel polished, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-040-600	Aluminium weight case 314-040-600 (for nominal values 1 g - 200 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-712	Marking (Last digits of serial number) KERN 920-712
921-217	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-217
962-217R	DAkKS Recalibration 962-217R
963-217	DAkKS calibration certificate 963-217

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days