

Product brochure

08.08.2026

KERN 304-056

Weight set, E1, 1 g - 500 g, Knob, stainless steel polished, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-050-600	Aluminium weight case 314-050-600 (for nominal values 1 g - 500 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-713	Marking (Last digits of serial number) KERN 920-713
921-218	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-218
962-218R	DAkkS Recalibration 962-218R
963-218	DAkkS calibration certificate 963-218

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days