

Product brochure

07.08.2026

KERN 304-06

Weight set, E1, 1 g - 1 kg, Knob, stainless steel polished, in wooden case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-060-100	Wood weight case 313-060-100 (for nominal values 1 mg - 1 kg for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-714	Marking (Last digits of serial number) KERN 920-714
921-219	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-219
962-219R	DAkKS Recalibration 962-219R
963-219	DAkKS calibration certificate 963-219

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days