

KERN 304-066

KERN

Weight set, E1, 1 g - 1 kg, Knob, Stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 314-060-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 1 kg
OIML-compliant	✓
Weight OIML class	E1

Construction

Dimension (W×D×H)	230×140×110 mm
Shape weight	Knob
Material	Stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	963-219
Article number for guaranteed plus-deviation	921-219

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	2,7 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg